


# eclipse & copernic

Compliant and hypercompliant balloon catheters

Designed to temporarily stop or control blood flow, to treat vasospasms, and for embolization of aneurysms with balloons.



balloon  
size from 4 to 6  
& 3 to 5mm



## Navigation

Hydrophilic coating  
on both the catheter & the balloon  
for a smooth progression

Compatibility with  
hybrid1214D guidewire  
which ensures support & trackability



## Compliant & hypercompliant balloons


Made from a very supple elastomer\*  
to easily conform to the anatomy

Rapid inflation/deflation  
with 1/2 contrast\*

Good visibility  
thanks to 2RX markers providing the  
balloon length



With courtesy of  
Dr Sukalyan Purkayastha, India



## Remodeling


Both Eclipse and Copernic  
can be used in parallel with  
a coiling microcatheter in a  
single 6F guiding catheter  
(FARGOMAX6F)

aneurysm treatment

# occlusion catheters with balloon


## copernic

Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
SCOPERNIC10	160	10	3 to 5	0.25	Compatible with .012" guidewire (packaging includes a HYBRID1214D)
SCOPERNIC15		15		0.35	
SCOPERNIC20		20		0.4	
SCOPERNIC30		30		0.5	
COPERNIC10	160	10	3 to 5	0.25	Compatible with .012" guidewire (not included in the package)
COPERNIC15		15		0.35	
COPERNIC20		20		0.4	
COPERNIC30		30		0.5	



## eclipse

Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
SECUPSE7	160	7	4 to 6	0.2	Compatible with .012" guidewire (packaging includes a HYBRID1214D)
SECUPSE9		9		0.25	
SECUPSE12		12		0.3	
SECUPSE15		15		0.4	
SECUPSE20		20		0.5	
ECLIPSE7	160	7	4 to 6	0.2	Compatible with .012" guidewire (not included in the package)
ECLIPSE9		9		0.25	
ECLIPSE12		12		0.3	
ECLIPSE15		15		0.4	
ECLIPSE20		20		0.5	



## baltacci

Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
BALTACCIB1	163	6 to 9	4 to 6	0.2	/



## magic

Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
MAGICB1	155	9	6	0.2	/
MAGICB2	155	18	6	0.5	/