

Phloxine B, C.I. 45410

Acid Red 92, Eosin 10B, Eosin Blue, BSC certified stain. For Methylene Blue-Phloxine staining and staining keratin and mucus according to Kreyberg.

25 g PB-P-25

Ponceau S, C.I. 27195

Acid Red 112, Fast Ponceau 2B. For use in electrophoresis.

25 g PS-P-25

Potassium permanganate

Synonym: Potassium (VII) manganate, Condy's crystals. Potassium permanganate for use in microscopy.

50 g KPH-P-50
100 g KPH-P-100

Pyronine B, C.I. 45010

Pyronine B, BSC certified stain. For bacterial, mold and RNA staining.

25 g PYB-P-25

Resorcinol, C.I. 76505

Resorcin, 1,3-benzendiol. Used for differentiation between aldoses and ketoses in Seliwanoff's test by staining ketose sugars dark red.

25 g RES-P-25

Rhodamine B, C.I. 45170

Basic Violet, Acid Red 52. For staining acid-fast bacteria.

25 g RHB-P-25

Rhodanine

DMAB-Rhodanine. For staining deposits of copper bound to proteins.

5 g RHO-P-5

Rose Bengal, C.I. 45440

Acid Red 94, BSC certified stain. For staining in microscopy and histology.

25 g RB-P-25

Safranin O, C.I. 50240

Basic Red 2, Safranin T, Phenazinium, BSC certified stain. Counterstain in the Gram method for detecting Gram-negative bacteria.

25 g SO-P-25

Sirius Red, C.I. 35780

Direct Red 80. Used for staining collagen and amyloid in histology.

25 g SR-P-25

Sudan III, C.I. 26100

Solvent Red 23, Sudan G, Fat Ponceau G, BSC certified stain. For lipid staining according to Lillie-Ashburn.

25 g S3-P-25

Sudan IV, C.I. 26105

Solvent Red 24, Scarlet Red, Biebrich Scarlet R, Fat Ponceau. For staining lipids according to Churukian's modification of Lillie-Ashburn's method.

25 g S4-P-25

Sudan Black B, C.I. 26150

Sudan Black, Solvent Black 3, Fat Black HB. For staining lipids and other hydrophobic materials.

5 g SB-P-5
10 g SB-P-10

Toluidine Blue O, C.I. 52040

Basic Blue 17, Tolonium Chloride, BSC certified stain. For staining mast cells, cartilage and Nissl bodies in the spinal cord.

25 g TB-P-25
100 g TB-P-100

Victoria Blue B, C.I. 44045

Basic Blue 26, Fat Blue B. For use in cytology and histology.

25 g VBB-P-25

Victoria Blue R, C.I. 42563

Basic Blue 8, Fat Blue 4R. For use in microbiology and histology.

25 g VBR-P-25

Water Blue, C.I. 42755

Acid Blue 22, BSC certified stain. For staining connective tissue according to Mallory and Gomori trichrome methods.

25 g WB-P-25

Wright's Stain

Colorants de Wright, BSC certified stain. Part of polychromatic Romanowsky stains commonly used in hematology.

25 g WS-P-25

Xylidine Ponceau, C.I. 16150

Ponceau 2R, Ponceau de Xylidine, BSC certified stain. For Masson trichrome staining.

25 g XP-P-25

