


UNDERPAD
REF.NO:800.003C-1.02
PRODUCT DATA SHEET



PRODUCT DATA SHEET
UNDERPAD-NON STERILE
800.003C-1.02

1) PRODUCT COMPONENTS



ITEM:1

Underpad
 90x180 cm
 Absorbent Area
 50x80 cm

2) DETAILS OF PRODUCT COMPONENTS

ITEM	QTY	PRODUCT	PRODUCTION TYPE	DIMENSION	COLOUR
1	1	Underpad	Total 60 gr (±3)	90x180 cm (±3)	White
			Backsheet – PE Film 18 gsm		
			Topsheet 8 gsm		
			Pulp-Absorbent Area-25.52 gr	50x80 cm (±3)	

3) TEST RESULT (ISO 11948-1)

ABSORPTION CAPACITY (ML)*	891
----------------------------------	-----

PRODUCT DATA SHEET
UNDERPAD-NON STERILE
800.003C-1.02

4)GENERAL INFORMATION

Product Description	UNDERPAD Non-Sterile
Certificates	EN ISO 13485:2016
Standards	EN ISO 14971:2019
	EN ISO 15223-1:2021
Intented Use	Underpad product is used as a mattress protector in health institutions such as hospitals, polyclinics, intensive care units, and home care patients.
Instructions For Usage	<p>Cut the products bag open with a scissor from the bottom side of the packaging. Doing so will present you a better place to hold onto the pad as you take it out of the packaging.</p> <p>Take the underpad out and place it in a folded-up state, onto the surface you'll be using it on.</p> <p>Unfold the pad. Take the edge of the product and unfurl it outward, away from you. Allow it to unfold completely.</p> <p>Place the pad down on the surface, with the non plastic side up. The soft side absorbs the urine, while the plastic-looking side will prevent any liquid from going through and onto the surface. The person will lie down on the softer side.</p>

5) PACKAGING

Packaging

Bag in 30 pcs