

常州市回春医疗器材有限公司

[Http://www.huichun.cn](http://www.huichun.cn) E-mail:huichun@huichun.cn

Changzhou Huichun Medical Equipment Co.,Ltd.


地址：江苏省常州市郑陆镇三河口 邮编：213115

Add: Sanhekou, Zhenglu Town, Changzhou, Jiangsu, China

TEL: 0519-88670653 88675093 88675214

FAX: 0519-88910533 88671528






®

连叶鱼

Huichun Medical



常州市回春医疗器材有限公司
Changzhou Huichun Medical Equipment Co.,Ltd.



H u i c h u n M e d i c a l

励精图治 务实创新

Make every effort to
pragmatism and innovation



连叶鱼

企业精神：励精图治、务实创新、与时俱进

发展理念：为客户提供满意的产品和服务是回春人永恒的追求

Entrepreneurial spirit: good governance, pragmatic innovation, keeping pace with the times.

Development concept: Your satisfaction is our eternal pursuit!

SINCE 1969

关于我们 HUICHUN MEDICAL

我们历史流长 不忘初心

常州市回春医疗器材有限公司位于风景秀丽、环境怡人、工业发达、物产富饶的长江三角洲，北临长江，南靠京沪铁路和沪宁高速公路。常州民航直达全国各大、中城市，通讯、水、陆、空交通便利。

公司现有占地面积 40000m²，10 万级净化车间达 10000m²，职工 800 余名，固定资产 10000 万元。公司从事医疗器械生产已有四十多年历史，从一九八五年就开始研制生产以医用高分子为材料的一次性使用医疗器械。现主要产品为一次性使用输液器、一次性使用无菌注射器、一次性使用输液袋（PVC）、一次性使用输血器、一次性使用采血针、一次性使用静脉血样采集容器、一次性使用静脉输液针、一次性使用活体取样钳、一次性使用负压引流装置、一次性使用输注泵、一次性使用膀胱冲洗器，并逐步形成从造粒、拉管、吹塑、注塑、装配、灭菌到成品检测一条龙生产规模。年生产一次性使用输液器 8000 万套，一次性使用输血器 300 万套，一次性使用无菌注射器 100000 万支。产品覆盖国内 20 多个省、直辖市，远销欧洲，南美洲，亚洲，非洲，中东等 110 多个国家。

公司坚持以质量为中心，重视企业生产经营和质量管理，注重产品质量和企业品牌，质量体系及服务网络健全，深受用户好评。1999 年“连叶鱼”牌一次性使用输液器、注射器被中国消费品市场公认著名品牌宣传活动组委会评为 99’ 消费品市场公认消费推荐品牌，99-2000 年连续两年被中国质量万里行评为无投诉单位，连续四年被评为武进市科技进步先进企业。2002 年 9 月通过中国 CMDCISO9001 产品质量及体系公证。2004 年通过了德国 TUV 公司的 CE 和 ISO 认证。公司将不断跟踪当今国内先进水平，应用先进的技术和科学的管理，完善公司的质量体系，提高企业整体素质，为客户提供优质的产品和服务。


Changzhou Huichun Medical Equipment Co., Ltd. is located in Yangtze River Delta, 200 kilometers away from the Shanghai city, 10 kilometers away from the Changzhou train station, and 28 kilometers away from the Changzhou Airport.

Our company has more than 40 years history to make the medical products. From 1986 we began to develop the disposable products. Our main products are disposable syringes, disposable infusion sets, disposable blood transfusion sets, hypodermic needles and blood collection systems etc. And we have the equipments to do all the work in forming raw material to the finished products.

Now our company covers an area of about 40000 square meters and 10000 square meters Class 100000 clean rooms. We have more than 800 employees and over 90 technicians. Our company is capable of producing annually 1000,000,000 pieces of syringes, 80,000,000 pieces of infusion sets. Most of the products are for exporting. We are now exporting our products to over 110 countries in Middle East, Africa, South America, Southeast Asia, East and West Europe.

We insist on the rule that quality is the core, pay high attention to produce management and quality management, and make sure the quality of the products and service. In 2002, the company passed ISO and CMD. In 2004, passed CE certification by TUV.

Huichun Medical considers the quality of products as a matter of life and death for the enterprise, and keeps abreast of new technology to improve the quality of products. Our company commits herself to the purpose "Quality first, customers first" and has set up a fast and high-efficient after-sales service system to meet demand of all customers.





公司将不断跟踪当今国内先进水平，应用先进的技术和科学的管理，完善公司的质量体系，提高企业整体素质，为客户提供优质的产品和服务。

Our company commits herself to the purpose "Quality first, customers first" and has set up a fast and high-efficient after-sales service system to meet demand of all customers.



我们重视高效设备的不断投入


Introduction of advanced equipment and technology, manufacturing high-quality products carefully return the old and new customers.

CLASS 100K GMP CLEANING WORKSHOPS, ADVANCED FACILITIES



Make every effort to
Pragmatism and innovation





SINCE 1969

我们重视规则与细节


无论是昨天，今天还是明天。。。。。

We pay attention to the rules and details Yesterday, today and tomorrow.....

Huichun Medical considers the quality of products as a matter of life and death for the enterprise, and keeps abreast of new technology to improve the quality of products. Our company commits herself to the purpose "Quality first, customers first" and has set up a fast and high-efficient after-sales service system to meet demand of all customers.

SINCE 1969







Http: www.huichun.cn E-mail: huichun@huichun.cn



Huichun

常州市回春医疗器材有限公司
Changzhou Huichun Medical Equipment Co.,Ltd.

一次性使用输液器系列 Infusion Sets

Disposable Infusion Sets: Vented Spike(Non-Vented Spike), 10/20/60 Drops= 1 ± 0.1 ml, Flexible Drip Chamber, Roller Clamp, Injection Port (Latex/Latex-free injection Site, Y-injection Port), Slip or Lock Adapter, For Gravity Use, Sterile





Bag-type infusion sets for single use

一次性使用袋式输液器



Transfusion Sets

一次性使用输血器



Transfusion Sets

一次性使用输血器

Huichun

Http: www.huichun.cn E-mail: huichun@huichun.cn

常州市回春医疗器材有限公司

一次性使用注射器系列 Syringes



3 件套注射器 (螺口) 3-Parts Syringe(luer lock)



3 件套注射器 (直口) 3-Parts Syringe(Luer Slip)



2 件套注射器 2-Parts Syringe



针管外径 Pipe Outside Diameter		色标 Colour Code	针管长度 Pipe Length
O.D.(mm)	Gauge	Colour	
0.25	31G	蓝色 Blue	1/2"
0.3	30G	淡黄 Light-yellow	1/2"
0.33	29G	橙红 Orange Red	1/2"
0.4	27G	中灰 Grey	5/8"
0.45	26G	褐 Brown	5/8"
0.5	25G	橙 Orange	3/4"
0.55	24G	中紫 Purple	1"
0.6	23G	深蓝 Blue	1" 11/4"
0.7	22G	黑 Black	11/8" 11/4"
0.8	21G	深绿 Green	13/8" 11/2"
0.9	20G	黄 Yellow	11/2"
1.1	19G	淡黄 Pea Yellow	11/2"
1.2	18G	粉红 Pink	11/8" 11/2"
1.6	16G	白 White	11/2"
1.8	15G	淡绿 Pea Green	11/2"



一次性使用胰岛素注射器
Insulin Syringe



1ML 注射器, 结核菌注射器
1ML Syringe, Tuberculin Syringe



自毁式注射器
Self destruction Syringe



冲洗器
Catheter Tip Syringe


Syringes

- Different Size, Meet Different Requirements.
- Central luer slip/Central Luer Lock/off-center Luer Slip.
- Special Design to Prevent Plunger Slipping Out.
- Latex/Later Free Piston
- PE/Blister Individual Package
- W/needle, W/O needle
- Custom Needle Size Per Specification

常州市回春医疗器材有限公司

Http: www.huichun.cn E-mail: huichun@huichun.cn

一次性使用静脉输液针 Scalp Vein Sets



针管外径 Pipe Outside Diameter		色标 Colour Code	针管长度 Pipe Length
O.D.(mm)	Gauge	Colour	
0.4	27G	中灰 Grey	5/8"
0.45	26G	褐 Brown	5/8"
0.5	25G	橙 Orange	3/4"
0.55	24G	中紫 Purple	3/4"
0.6	23G	深蓝 Blue	1"
0.7	22G	黑 Black	1"
0.8	21G	深绿 Green	11/8"
0.9	20G	黄 Yellow	11/8"
1.1	19G	淡黄 Pea Yellow	11/4"
1.2	18G	粉红 Pink	11/4"



一次性使用胰岛素针 Insulin Pen Needle

SIZE	Needle OD (mm)	Common needle effective length (mm)
28G	0.36	8,12
29G	0.33	10,12
30G	0.30	8,10
31G	0.25	4,5,6,8
32G	0.23	4
33G	0.21	4
34G	0.18	4

Remark: normal size in above table, we can also make special size as clients'request.



一次性使用采血针 Blood collection needle (Pen type / Scalp Vein Set type)

Huichun

常州市回春医疗器材有限公司

Http: www.huichun.cn

E-mail: huichun@huichun.cn



Make every effort to
Pragmatism and innovation

Huichun Medical [Http: www.huichun.cn](http://www.huichun.cn) E-mail: huichun@huichun.cn

常州市回春医疗器材有限公司

我们重视高效设备的不断投入

We attach great importance to the continuous input of efficient equipment


Introduction of advanced equipment and technology, manufacturing high-quality products carefully return the old and new customers.



尿袋 – 输液袋 – 扩张器 – 管类 – 氧气
面罩 – 活检钳 – 其他

Urine bag Infusion bag Disposable Vaginal Speculum Series
Tubing series
Others





一次性使用

尿袋系列 Urine bag



2000ml 豪华尿袋 | 2000ml Luxury Urine Bag |



可选配件

Different bottom drainage valves



2000ml 豪华尿袋 | 2000ml Luxury Urine Bag




2000ml 豪华尿袋 | 2000ml Luxury Urine Bag

Urine Collection Bags


- Medical Grade PVC
- 1000/1500/2000ml Collection Capacity
- Kink Resistant Tubing and Non-return Valve
- Conical Connector and Protective Cap
- Bottom Draining
- Sterile



腿袋 Urine Bag For Leg




婴儿尿袋 100ml Paediatric Urine Bag



Huichun

常州市回春医疗器材有限公司 Http: www.huichun.cn E-mail: huichun@huichun.cn



一次性使用


输液袋系列 Infusion bag



漏斗聚碳酸酯脂塞子
Crateriform polycarbonate cap



大号聚碳酸酯脂塞子
Large-size Polycarbonate cap




蝴蝶形易拉
Butterfly twist off

Infusion Bags

- Medical Grade PVC
- 50ml, 100ml, 250ml, 500ml, 1000ml Bag
- In Port and Out Port
- Injection Site


阴道扩张器系列 Disposable Vaginal Speculum Series




勾扣式阴道扩张器 Vaginal Speculum with Hook



中螺旋阴道扩张器 Vaginal Speculum with Middle Screw




侧螺旋阴道扩张器 Vaginal Speculum with Lateral Screw



法式阴道扩张器 Vaginal Speculum- French Type

Huichun

常州市回春医疗器材有限公司 Http: www.huichun.cn E-mail: huichun@huichun.cn



一次性使用

管类 Tubing series



肛门管 Rectal Tube



胃肠管 Stomach&Intestines Tube



喂食管 Feeding Tube



鼻氧管
Nasal Oxygen Cannula



连接管
Disposable Connection Tube

Huichun

Http: www.huichun.cn

E-mail: huichun@huichun.cn



氧气面罩 Oxygen Mask






雾化面罩 Nebulizer Oxygen Mask





识别带 I.D. Bracelet

活检钳 Biopsy Forceps



A	平口无针、不包塑 Uncoated Without spike		
	规格型号 Model	钳头直径 O.D.(mm)	长 度 Length(mm)
	QYQ-A2. 3X1600	2.3	1600
	QYQ-A2. 3X1800	2.3	1800
	QYQ-A2. 3X2300	2.3	2300
	QYQ-A1. 8X1050	1.8	1050
	QYQ-A1. 8X1200	1.8	1200
	QYQ-A1. 8X1600	1.8	1600

B	平口带针、不包塑 Uncoated With spike		
	规格型号 Model	钳头直径 O.D.(mm)	长 度 Length(mm)
	QYQ-B2. 3X1600	2.3	1600
	QYQ-B2. 3X1800	2.3	1800
	QYQ-B2. 3X2300	2.3	2300

AXC	平口无针、包塑 Coated Without spike		
	规格型号 Model	钳头直径 O.D.(mm)	长 度 Length(mm)
	QYQ-AXC2. 3X1600	2.3	1600
	QYQ-AXC2. 3X1800	2.3	1800
	QYQ-AXC2. 3X2300	2.3	2300
	QYQ-AXC1. 8X1050	1.8	1050
	QYQ-AXC1. 8X1200	1.8	1200
	QYQ-AXC1. 8X1600	1.8	1600

BXC	平口带针、包塑 Coated With spike		
	规格型号 Model	钳头直径 O.D.(mm)	长 度 Length(mm)
	QYQ-BXC2. 3X1600	2.3	1600
	QYQ-BXC2. 3X1800	2.3	1800
	QYQ-BXC2. 3X2300	2.3	2300

Huichun

常州市回春医疗器材有限公司 Http: www.huichun.cn E-mail: huichun@huichun.cn

一次性使用负压引流装置


The character of Silicone Suction Reservoir and Flat Drain:

- The product is made of medical silicone
- The capacity: 100ml, 200ml
- Size of draining tube: F10,F16
- Supplied sterile in individual sterile blister pack
- For single use only

本装置采用硅橡胶材料制成，由球囊、塞子、单向阀、球囊套、锁紧套、悬挂环、导管、引流管、管接头组成，按容量不同分两种。供临床术后负压引流和贮存引流液用。



其他产品 Others



一次性使用外伤冲流连接管 Disposable Irrigation Set

用于外部伤口冲洗时连接在外部冲洗防护罩与冲洗器之间



三通阀 Three way stopcock



一次性使用大便化验杯 Urine Cup



Huichun

常州市回春医疗器材有限公司 [Http: www.huichun.cn](http://www.huichun.cn) E-mail: huichun@huichun.cn