

8 Field Replaceable Unit

This chapter lists the details of the Field Replaceable Units (FRU) of the system. Please refer to the explosive view and the FRU table below.


Level 0 represents the main assemblies of the system, and the whole system is divided into 7 assemblies.

Level X represents the main parts of Level 0 assembly.

Level X-1 represents the sub-parts of Level X.


If order number is provided in the table, which means the part is an FRU and can be ordered; if the order number is not provided in the table, which means the part is not an FRU and cannot be ordered.

8.1 Explosive View



8.2 Assembly Explosive View

8.2.1 Monitor Assembly (A0)



No.	Order Number	Part Name	Qty.	Remark
A1	115-057736-00	21.5" Monitor Assembly(DC-33/FRU)	1	
A1	115-057739-00	21.5" Monitor Assembly(DC-35/FRU)	1	
A1	115-057740-00	21.5" Monitor Assembly(DC-36/FRU)	1	
A1	115-057741-00	21.5" Monitor Assembly(DC-38/FRU)	1	
A1	115-057742-00	21.5" Monitor Assembly(DC-39/FRU)	1	
A1	115-057743-00	21.5" Monitor Assembly(DC-40/FRU)	1	
A1	115-057744-00	21.5" Monitor Assembly(DC-41/FRU)	1	
A1	115-057745-00	21.5" Monitor Assembly(DC-42/FRU)	1	
A1	115-057746-00	21.5" Monitor Assembly(DC-43/FRU)	1	

No.	Order Number	Part Name	Qty.	Remark
A1	115-057747-00	21.5" Monitor Assembly(DC-44/FRU)	1	
A1	115-057748-00	21.5" Monitor Assembly(DC-45/FRU)	1	
A1	115-057749-00	21.5" Monitor Assembly(DC-40S/FRU)	1	
A1	115-057750-00	21.5" Monitor Assembly(DC-40 Pro/FRU)	1	
A1	115-057751-00	21.5" Monitor Assembly(DC-40 Exp/FRU)	1	
A1	115-057752-00	21.5" Monitor Assembly(DC-40T/FRU)	1	
A1	115-068382-00	21.5" Monitor Assy(DC-33/FHD/FRU)	1	
A1	115-068383-00	21.5" Monitor Assy(DC-35/FHD/FRU)	1	
A1	115-068384-00	21.5" Monitor Assy(DC-36/FHD/FRU)	1	
A1	115-068387-00	21.5" Monitor Assy(DC-38/FHD/FRU)	1	
A1	115-068388-00	21.5" Monitor Assy(DC-39/FHD/FRU)	1	
A1	115-068389-00	21.5" Monitor Assy(DC-40/FHD/FRU)	1	
A1	115-068390-00	21.5" Monitor Assy(DC-41/FHD/FRU)	1	
A1	115-068391-00	21.5" Monitor Assy(DC-42/FHD/FRU)	1	
A1	115-068393-00	21.5" Monitor Assy(DC-43/FHD/FRU)	1	
A1	115-068394-00	21.5" Monitor Assy(DC-44/FHD/FRU)	1	
A1	115-068395-00	21.5" Monitor Assy(DC-45/FHD/FRU)	1	
A1	115-068396-00	21.5" Monitor Assy(DC-40S/FHD/FRU)	1	
A1	115-068397-00	21.5" Monitor Assy(DC-40 Pro/FHD/FRU)	1	
A1	115-068398-00	21.5" Monitor Assy(DC-40 Exp/FHD/FRU)	1	
A1	115-068399-00	21.5" Monitor Assy(DC-40T/FHD/FRU)	1	
A1	115-068199-00	21.5" Monitor Assy(DC-Eva/FHD/FRU)	1	

8-4 Field Replaceable Unit


No.	Order Number	Part Name	Qty.	Remark
A1	115-068200-00	21.5" Monitor Assy(DC-N3 Elite/FHD/FRU)	1	
A1	115-068201-00	21.5" Monitor Assy(DC-N3 Exp/FHD/FRU)	1	
A1	115-068202-00	21.5" Monitor Assy(DC-35 Exp/FHD/FRU)	1	
A1	115-068203-00	21.5" Monitor Assy(DC-35 Pro/FHD/FRU)	1	
A1	115-068204-00	21.5" Monitor Assy(DC-39S/FHD/FRU)	1	
A1	115-068205-00	21.5" Monitor Assy(DC-38S/FHD/FRU)	1	
A1	115-068206-00	21.5" Monitor Assy(DC-36S/FHD/FRU)	1	
A1	115-068207-00	21.5" Monitor Assy(DC-35S/FHD/FRU)	1	
A1	115-068208-00	21.5" Monitor Assy(DC-33S/FHD/FRU)	1	
A1	115-068400-00	21.5" Monitor Assy(DC-33/Crystal/FRU)	1	
A1	115-068401-00	21.5" Monitor Assy(DC-35/Crystal/FRU)	1	
A1	115-068402-00	21.5" Monitor Assy(DC-36/Crystal/FRU)	1	
A1	115-068403-00	21.5" Monitor Assy(DC-38/Crystal/FRU)	1	
A1	115-068404-00	21.5" Monitor Assy(DC-39/Crystal/FRU)	1	
A1	115-068405-00	21.5" Monitor Assy(DC-40/Crystal/FRU)	1	
A1	115-068406-00	21.5" Monitor Assy(DC-41/Crystal/FRU)	1	
A1	115-068407-00	21.5" Monitor Assy(DC-42/Crystal/FRU)	1	
A1	115-068408-00	21.5" Monitor Assy(DC-43/Crystal/FRU)	1	
A1	115-068409-00	21.5" Monitor Assy(DC-44/Crystal/FRU)	1	
A1	115-068410-00	21.5" Monitor Assy(DC-45/Crystal/FRU)	1	
A1	115-068411-00	21.5" Monitor Assy(DC-40S/Crystal/FRU)		
A1	115-068412-00	21.5" Monitor Assy(DC-40 Pro/Crystal/FRU)	1	

No.	Order Number	Part Name	Qty.	Remark
A1	115-068413-00	21.5" Monitor Assy(DC-40 Exp/Crystal/FRU)	1	
A1	115-068414-00	21.5" Monitor Assy(DC-40T/Crystal/FRU)	1	
A1	115-068209-00	21.5" Monitor Assy(DC-Eva/Crystal/FRU)	1	
A1	115-068210-00	21.5" Monitor Assy(DC-N3 Elite/Crystal/FRU)	1	
A1	115-068211-00	21.5" Monitor Assy(DC-N3 Exp/Crystal/FRU)	1	
A1	115-068212-00	21.5" Monitor Assy(DC-35 Exp/Crystal/FRU)	1	
A1	115-068213-00	21.5" Monitor Assy(DC-35 Pro/Crystal/FRU)	1	
A1	115-068214-00	21.5" Monitor Assy(DC-39S/Crystal/FRU)	1	
A1	115-068215-00	21.5" Monitor Assy(DC-38S/Crystal/FRU)	1	
A1	115-068216-00	21.5" Monitor Assy(DC-36S/Crystal/FRU)	1	
A1	115-068217-00	21.5" Monitor Assy(DC-35S/Crystal/FRU)	1	
A1	115-068218-00	21.5" Monitor Assy(DC-33S/Crystal/FRU)	1	
A1	115-078457-00	23.8" Monitor Assy(FRU)	1	
A2	047-034592-00	Lable(full HD)	1	
A2	047-034593-00	Lable(Crystal)	1	
A3	047-025310-00	DC-39lable	1	
A3	047-025311-00	DC-41lable	1	
A3	047-025312-00	DC-42lable	1	
A3	047-025313-00	DC-43lable	1	
A3	047-025314-00	DC-44lable	1	
A3	047-025315-00	DC-45lable	1	
A3	047-025316-00	DC-40Slable	1	

8-6 Field Replaceable Unit

No.	Order Number	Part Name	Qty.	Remark
A3	047-025317-00	DC-40Prolable	1	
A3	047-025318-00	DC-40Explable	1	
A3	047-025319-00	DC-40Tlable	1	
A3	047-025320-00	DC-33lable	1	
A3	047-025321-00	DC-35lable	1	
A3	047-025322-00	DC-36lable	1	
A3	047-025323-00	DC-38lable	1	
A3	047-025324-00	DC-40lable	1	
A3	047-034727-00	lable (DC-33S)	1	
A3	047-034728-00	lable (DC-35S)	1	
A3	047-034729-00	lable (DC-36S)	1	
A3	047-034730-00	lable (DC-38S)	1	
A3	047-034731-00	lable (DC-39S)	1	
A3	047-034732-00	lable (DC-35 Pro)	1	
A3	047-034733-00	lable (DC-35 Exp)	1	
A3	047-034734-00	lable (DC-N3 Exp)	1	
A3	047-034735-00	lable (DC-N3 Elite)	1	
A3	047-034736-00	lable (DC-Eva)	1	


8.2.2 Monitor Support Arm Assembly (B0)



No.	Order Number	Part Name	Qty.	Remark
-----	--------------	-----------	------	--------


B1	115-057733-00	Monitor Arm Assembly(for fixed arm/FRU)	1	
B1-1	115-057836-00	arm cover(for fixed monitor arm)	1	
B2	115-057835-00	Monitor Arm Assembly(for adjustable arm/FRU)	1	
B2-1	115-057833-00	upper arm cover(for adjustable monitor arm)	1	
B2-2	115-057834-00	lower arm cover(for adjustable monitor arm)	1	

8.2.3 Main Control Panel Assembly (C0)




8-10 Field Replaceable Unit

No.	Order Number	Part Name	Qty.	Remark
C1	115-057726-00	Touch Screen Assembly(13.3"/FRU)	1	
C3	115-024667-01	Gel Warmer Module Kit	1	Including Gel Warmer machine and the bracket
C3-1	043-003498-00	Bracket of heater jack	1	Only the bracket, no warmer.
C4	020-000049-00	Single Speaker	1	
C5	043-008973-00	left speaker cover	1	Including cover and speaker net
C6	043-008974-00	right speaker cover	1	Including cover and speaker net
C7	115-038149-00	Small keyboard assembly(2136 FRU)	1	
C7	115-052165-00	Small keyboard assembly(Russia/2136 FRU)	1	
C8	115-037784-00	left probe holder assembly	1	
C9	115-037785-00	right probe holder assembly	1	
C10	115-035706-00	DVD Assembly(2136)	1	Only including DVD machine
C10	115-057734-00	DVD Module Kit(2148)	1	Including DVD machine, cover
C10	115-078624-00	DVD decorate cover Assy(Non-DVD)	1	Only cover
C11	043-007013-00	Small probe holder(2136)	1	
C11	043-009106-00	Small Probe Holder	1	Released in DC-40 V3.0
C12	043-007014-00	Middle probe holder(2136)	1	
C12	043-009105-00	Big Probe Holder	1	Released in DC-40 V3.0
C13	043-007015-00	Large probe holder(2136)	1	
C13	043-009107-00	4D Probe Holder	1	Released in DC-40 V3.0
C14	043-007036-00	Pencil probe holder(2136)	1	




No.	Order Number	Part Name	Qty.	Remark
C3	115-024667-01	Gel Warmer Module Kit	1	including Gel Warmer machine and the bracket
C3-1	043-003498-00	Bracket of heater jack	1	only the bracket, no warmer.



No.	Order Number	Part Name	Qty.	Remark
C15	043-003977-00	Small Encoder Cap	1	
C16	043-003978-00	Big Encoder Cap	1	


No.	Order Number	Part Name	Qty.	Remark
C17	043-004910-00	TGC Cap	1	
C18	115-057728-00	Top Cover of Keyboard Assembly(2148 FRU)	1	
C18	115-057728-02	Top Cover of Keyboard Assembly(2148 FRU)	1	EAF497F released. The Doppler software must be upgraded to V3.3.2 or above after repairing the old machine.
C18	115-068219-00	Top Cover of Keyboard Assy(2148/V3.0FRU)	1	Released in DC-40 V3.0
C18	115-068219-01	Top Cover of Keyboard Assy(2148/V3.0FRU)	1	EAF497F released. The Doppler software must be upgraded to V3.3.2 or above after repairing the old machine.
C19	023-001236-00	2" trackball(2136)	1	
C19	115-068417-00	2" trackball(2148/V3.0/FRU)	1	Released in DC-40 V3.0
C20	801-2119-00034-00	TGC Silica	1	
C21	115-057729-00	TGC Board	1	
C22	051-001471-00	Encoder Board PCBA	1	



No.	Order Number	Part Name	Qty.	Remark
C23	049-000590-00	control panel silica gel keys	1	
C24	801-2300-00008-00	Buzzer	1	
C26	047-024572-00	2136 overlay(RU)	1	
C26	047-024573-00	2136 overlay(DE)	1	
C26	047-024574-00	2136 overlay(FR)	1	
C26	047-024575-00	2136 overlay(IT)	1	


No.	Order Number	Part Name	Qty.	Remark
C26	047-024576-00	2136 overlay(PL)	1	
C26	047-024577-00	2136 overlay(ES)	1	
C26	047-024578-00	2136 overlay(PT)	1	
C26	047-024579-00	2136 overlay(CZ)	1	

8.2.4 Main Unit Assembly (E0)



No.	Order Number	Part Name	Qty.	Remark
E1	115-057890-00	ECG Module Assembly(2148/AHA)	1	
E1	115-068226-00	ECG Module Assembly(2148/V3.0/AHA)	1	Released in DC-40 V3.0
E1	115-057735-00	ECG Module Assembly(2148/IEC)	1	
E1	115-068225-00	ECG Module Assembly(2148/V3.0/IEC)	1	Released in DC-40 V3.0
E1-1	115-057819-00	Front-in Assembly(2148/FRU)	1	
E1-1	115-068419-00	Front-in Assembly(2148/v3.0//FRU)	1	Released in DC-40 V3.0
E1-2	115-045957-00	Pen probe upgrading kit(FRU)	1	extension cable+CW board
E2	043-010003-00	Top Cover (2148)	1	
E2	115-068231-00	Top Cover (2148/V3.0/FRU)	1	Released in DC-40 V3.0
E3	044-001270-00	Back Handle(2148)	1	
E3	115-068232-00	Back Handle(2148/V3.0/FRU)	1	Released in DC-40 V3.0
E4	043-002325-00	Tray-S	1	
E5	115-057349-00	Front Cover Assembly(2148/FRU)	1	
E5	115-068224-00	Front Cover Assembly(2148/V3.0/FRU)	1	Released in DC-40 V3.0
E6	115-057350-00	Back Cover Assembly(2148)	1	
E6	115-068223-00	Back Cover Assembly(2148/V3.0/FRU)		Released in DC-40 V3.0
E7	115-041710-00	Left Cover Assembly(2136,CE&FDA,FRU)	1	
E8	115-038151-00	Right cover assembly(2136 FRU)	1	
E9	043-002324-00	Tray-B	1	
E28	043-010250-00	2148 printer door	1	printer capsule plastic cover
E28	115-068422-00	printer cabinet cover(2148 V3.0/FRU)	1	Released in DC-40 V3.0

8-18 Field Replaceable Unit




No.	Order Number	Part Name	Qty.	Remark
E10	115-064358-00	HDD(2148/DC-35/FDA/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-064359-00	HDD(2148/DC-40/FDA/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-064360-00	HDD(2148/DC-45/FDA/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-064361-00	HDD(2148/DC-40S/FDA/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-064362-00	HDD(2148/DC-40 Pro/FDA/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057774-00	HDD(2148/DC-33/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057775-00	HDD(2148/DC-35/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057776-00	HDD(2148/DC-36/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057777-00	HDD(2148/DC-38/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057778-00	HDD(2148/DC-39/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then


No.	Order Number	Part Name	Qty.	Remark
				restore or reinstall the date after replacing the HDD
E10	115-057779-00	HDD(2148/DC-40/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057780-00	HDD(2148/DC-41/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057781-00	HDD(2148/DC-42/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057782-00	HDD(2148/DC-43/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057783-00	HDD(2148/DC-44/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057784-00	HDD(2148/DC-45/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057785-00	HDD(2148/DC-40S/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057786-00	HDD(2148/DC-40 Pro/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-057787-00	HDD(2148/DC-40 Exp/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD

No.	Order Number	Part Name	Qty.	Remark
E10	115-057788-00	HDD(2148/DC-40T/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068189-00	HDD(2148/DC-Eva/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068190-00	HDD(2148/DC-N3 Elite/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068191-00	HDD(2148/DC-N3 Exp/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068192-00	HDD(2148/DC-35 Exp/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068193-00	HDD(2148/DC-35 Pro/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068194-00	HDD(2148/DC-39S/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068195-00	HDD(2148/DC-38S/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068196-00	HDD(2148/DC-36S/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E10	115-068197-00	HDD(2148/DC-35S/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then

No.	Order Number	Part Name	Qty.	Remark
				restore or reinstall the date after replacing the HDD
E10	115-068198-00	HDD(2148/DC-33S/CE/FRU)	1	Please indicate the doppler software version when you apply. If possible, it is recommend to backup the old date such as preset date/patient date, and then restore or reinstall the date after replacing the HDD
E11	115-045956-00	CW Module(FRU)	1	Includes CW board and connecting cable for pen probe(009-002541-00)
E12	115-038200-00	Transmit & Receive Board(2136 FRU)	1	Please indicate the doppler software version when you apply
E13	115-035614-00	mother board assembly(2136)	1	
E14	115-057727-00	IO assembly(2148/FRU)	1	
E14	115-068227-00	IO assembly (2148/V3.0/FRU)	1	Released in DC-40 V3.0
E15	801-2120-00031-00	Small Fan Assembly(FRU)	1	outlet fan, one FRU contains 2 fans
E16	115-015308-00	Inlet Fan	1	Inlet Fan,one FRU contains 2 fans
E17	115-038202-00	DC-DC Power Board Assembly(2136 FRU)	1	
E18	115-038157-00	Power input assembly package (2136 FRU)	1	
E19	023-001528-00	WIFI Board(2148)	1	only including WIFI Board
E19	115-057730-00	WIFI Module Kit(2148)		including WIFI Board, WIFI antenna
E20	115-038143-00	4D Module Kit	1	installed on digital board
E21	115-018409-00	lithium-ion Battery Pack	1	14.8V 6600mAh(Except India)
E21	115-081622-00	lithium-ion Battery Pack(BIS)	1	14.8V 6600mAh(Only for India)
E21	115-018410-00	lithium-ion Battery Pack(FDA)	1	14.8V 6600mAh




No.	Order Number	Part Name	Qty.	Remark
E22	115-015306-00	AC-DC Power Board	1	
E23	115-015307-02	Power Auxiliary Output Link Board	1	




No.	Order Number	Part Name	Qty.	Remark
E24	115-057725-00	Digital Board Assembly(2148/FRU)	1	Not include CPU module. Please indicate the doppler software version when you apply
E24	115-057725-03	Digital Board Assembly(2148/FRU)	1	1. When using the new digital board FRU to repair the old digital board machine, you need to upgrade the Doppler software to V3.3.6 or above; 2. The old digital board FRU can directly repair the machine of the new and old digital board;
E25	115-015310-00	Li-ion Battery Assembly	1	The Li-ion battery and connect cable on digital board
E26	115-057724-00	PC Module Assembly(2148 FRU)	1	The MAC address will change with the replacement of CPU, so you may need apply the optional key the

No.	Order Number	Part Name	Qty.	Remark
				machine installed before
E26	115-057724-01	PC Module Assembly(2148 FRU)	1	The MAC address will change with the replacement of CPU, so you may need apply the optional key the machine installed before. Using this PN to replace 115-057724-00 needs recovery OS to enable the network function;
E26	115-057724-02	PC Module Assembly(2148 FRU)	1	The MAC address will change with the replacement of CPU, so you may need apply the optional key the machine installed before. Using this PN to replace 115-057724-00/115-057724-01 needs recovery OS V1.2.8 to enable the network function;
E20	115-038143-00	4D Module Kit	1	installed on digital board



No.	Order Number	Part Name	Qty.	Remark
E27	115-057731-00	Probe Board(2148/socket 3A+1E/FRU)	1	probe socket: 3 sockets of 156PIN & 1 socket of 260PIN
E27	115-057732-00	Probe Board(2148/socket 4A/FRU)	1	probe socket: 4 sockets of 156PIN

8.2.5 Base Assembly (F0)



No.	Order Number	Part Name	Qty.	Remark
F1	043-006967-00	Main unit pedal(2136)	1	
F1	115-068233-00	Main unit pedal(2148/V3.0/FRU)	1	Released in DC-40 V3.0
F2	115-035795-00	Main-base Assembly(2136)	1	
F2	115-068418-00	Main-base Assembly(2148/V3.0/FRU)	1	Released in DC-40 V3.0
F3	115-015314-00	Castor	1	

8.2.6 Cable (G0)

No.	Order Number	Part Name	Qty.	Remark
G1	009-002671-00	Signal cable of the Printer	1	
G2	009-009184-00	Touch Panel Cable Assembly(FRU)	1	
G3	115-057826-00	Cable for USB External Connection	1	
G4	009-002560-00	USB cable	1	
G5	009-004599-00	Trackball connecting cable	1	
G6	009-003635-00	encoder signal wires	1	
G7	2108-20-65853	STC connecting line	1	
G8	009-002540-00	pen probe cable	1	
G9	009-009186-00	USB cable for keyboard	1	
G10	009-009183-00	Main Cable Assembly(for fixed monitor arm)	1	for fixed monitor arm
G11	009-009182-00	Main Cable Assembly(for adjustable monitor arm)	1	for adjustable monitor arm

