


PAGINA 0

**CREATA PENTRU A
PUTEA SEMNA
ELECTRONIC
DOCUMENTUL**



ACTIVE
FLEXIBILITY

SUPPORTIVE
SECONDARY
CURVE

BALANCED PERFORMANCE THAT'S AHEAD OF THE CURVE

The Launcher™ guide catheter offers you a blend of flexibility, support and visualisation — the capability you need to respond to your challenging cases.

Advanced Platform Designed for Multiple Interventional Approaches


- Designed for radial and femoral approaches
- Soft distal segment enables you to actively engage the guide for backup support
- Full-Wall Technology construction provides kink resistance and stable torque control

Enhanced Visualisation

- Large 5 F I.D. (0.058 inches) allows for increased contrast flow for enhanced visualisation
- Vest-Tech Nylon outer jacket provides excellent hub to tip radiopacity

Full Range of Traditional and Specialty Curves

- Extensive collection of specialty transradial and brachial curves



MEDTRONIC FULL-WALL TECHNOLOGY:
Encapsulated flat-wire braid enables thinner, more robust walls without compromising support, retention, or torque response.

LAUNCHER™ 5 F GUIDE CATHETERS


LEFT CORONARY CURVES					Backup Support Left
Left Standard	Short Left	JCL	Amplatz Left	Short Amplatz Left	EBU
LA5JL20	LA5SL30	LA5JCL15	LA5AL75	LA5SAL75	LA5EBU20
LA5JL30	LA5SL35	LA5JCL30	LA5AL10	LA5SAL10	LA5EBU30
LA5JL35	LA5SL40	LA5JCL35	LA5AL15	LA5SAL15	LA5EBU35
LA5JL40	LA5SL45	LA5JCL40	LA5AL20	LA5SAL20	LA5EBU375
LA5JL45	LA5SL50	LA5JCL45	LA5AL25	LA5SAL25	LA5EBU40
LA5JL50	LA5SL50	LA5JCL50	LA5AL30	LA5SAL30	LA5EBU45
LA5JL60	LA5SL60		LA5AL40	LA5SAL40	LA5EBU50
RIGHT CORONARY CURVES					
Right Standard	Short Right	JCR	Amplatz Right	Shepherd's Crook Right	Backup Support Right
LA5JR30	LA5SR30	LA5JCR35	LA5AR10	LA5SCR35	LA5RBU35
LA5JR35	LA5SR35	LA5JCR40	LA5AR20	LA5SCR40	LA5RBU40
LA5JR40	LA5SR40	LA5JCR45	LA5ALR12	LA5SCR50	
LA5JR45	LA5SR50	LA5JCR50			ECR
LA5JR50	LA5SR60		Short Amplatz Right		LA5ECR35
LA5JR60			LA5SAR10		LA5ECR40
			LA5SAR20		LA5ECR45
					3-Dimensional
					LA53DRIGHT
OTHER CURVES					
Multipurpose	MAC (Multi-Aortic)	Bypass	Champ	Radial Specialty	NOTO
LA5MB1	LA5MAC30	LA5IMA	LA5CHAMP05	LA5LARA	LA5NOTO
LA5MB2	LA5MAC35	LA5LCB	LA5CHAMP10	LA5MRADIAL	
LA5MPST	LA5MAC375	LA5RCB	LA5CHAMP15	LA5MRESS	
LA5HSI	LA5MAC40	LA5RCBII	LA5CHAMP20	LA5RRAD	
LA5HSII	LA5MAC45		LA5CHAMP25		EL GAMAL
LA5HSIII	LA5MAC50		LA5CHAMP30		LA5ELGAMAL
LA5HSREL			LA5CHAMP35		
			LA5CHAMP40		

Shaft Diameter: 5 F Inner Diameter: 0.058 in. 5 packs available in select curves

Indications for Use: The Medtronic Guiding Catheter is designed to provide a pathway through which therapeutic devices are introduced. The guiding catheter is intended to be used in the coronary or peripheral vascular system.
Important Information: Refer to *Instructions for Use* supplied with the device for indications, warnings, precautions, and adverse effects.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
 The curve recommendations in this guide are not meant to replace patient-specific clinical judgement.



Medtronic
 3576 Unocal Place
 Santa Rosa, CA 95403
 USA
 Tel: 707.525.0111

LifeLine Customer Support
 Tel: 877.526.7890
 Tel: 763.526.7890

Product Services
 Tel: 888.283.7868
 Fax: 800.838.3103

Europe
Medtronic Intl. Trading SARL
 Route du Molliau 31
 Case Postale
 CH-1131 Tolochenaz
 Switzerland
 Tel: 41.21.802.7000
 Fax: 41.21.802.7900

Canada
Medtronic of Canada Ltd.
 99 Hereford Street
 Brampton, Ontario L6Y 0R3
 Canada
 Tel: 905.460.3800
 Fax: 905.460.3998
 Toll-free: 800.268.5346

Asia Pacific
Medtronic Intl. Ltd.
 49 Changi South Avenue 2
 Nasaco Tech Centre
 Singapore 486056
 Singapore
 Tel: 65.6436.5000
 Fax: 65.6776.6335

Latin America
Medtronic
 9850 Doral Blvd., 4th Floor
 Miami, FL 33178
 USA
 Tel: 786.709.4200
 Fax: 786.709.4244

medtronic.com