


Disposable Stone Extraction Basket

Indications

This product is operated under the endoscope during surgery to remove tissues, foreign bodies, and crush stones.



Features

- ※ The handle is equipped with an injection port for easy and smooth injection of contrast media
- ※ Good push performance, easier to reach the target location
- ※ Diamond-shaped mesh basket, easy to take and release, easy to repeat
- ※ With rotating function, the push-pull design of the control handle makes the stone grasping operation easier and more effective
- ※ A variety of product types, providing diamond type, spiral type, oval type, can meet the various clinical needs of doctors

Specifications

Code	Basket Type	Sheath outer diameter $\pm 0.2\text{mm}$	Working length $\pm 50\text{mm}$	Basket opening width $\pm 5\text{mm}$	Working Channel (mm)	Basket wire quantity
MB18-A-2420/07-4	Diamond Type A 	2.4	700	20	≥ 2.8	4
MB18-A-2420/07-6		2.4	700	20	≥ 2.8	6
MB18-A-2420/20-4		2.4	2000	20	≥ 2.8	4
MB18-A-2420/20-6		2.4	2000	20	≥ 2.8	6
MB18-A-2425/20-4		2.4	2000	25	≥ 2.8	4
MB18-A-2425/20-6		2.4	2000	25	≥ 2.8	6
MB18-A-2430/20-4		2.4	2000	30	≥ 2.8	4
MB18-A-2430/20-6		2.4	2000	30	≥ 2.8	6
MB18-B-2420/07-4	Spiral Type B 	2.4	700	20	≥ 2.8	4
MB18-B-2420/07-6		2.4	700	20	≥ 2.8	6
MB18-B-2420/20-4		2.4	2000	20	≥ 2.8	4
MB18-B-2420/20-6		2.4	2000	20	≥ 2.8	6
MB18-B-2425/20-4		2.4	2000	25	≥ 2.8	4
MB18-B-2425/20-6		2.4	2000	25	≥ 2.8	6
MB18-B-2430/20-4		2.4	2000	30	≥ 2.8	4
MB18-B-2430/20-6		2.4	2000	30	≥ 2.8	6
MB18-C-2420/07-4	Ellipse Type C 	2.4	700	20	≥ 2.8	4
MB18-C-2420/07-6		2.4	700	20	≥ 2.8	6
MB18-C-2420/20-4		2.4	2000	20	≥ 2.8	4
MB18-C-2420/20-6		2.4	2000	20	≥ 2.8	6
MB18-C-2425/20-4		2.4	2000	25	≥ 2.8	4
MB18-C-2425/20-6		2.4	2000	25	≥ 2.8	6
MB18-C-2430/20-4		2.4	2000	30	≥ 2.8	4
MB18-C-2430/20-6		2.4	2000	30	≥ 2.8	6