



## Cytology Brush

Used endoscopically to obtain mucosal tissue samples from gastrointestinal tract for microscopic-histopathological examination



CY  
not wire guided




CY-DL  
wire guided

Art. No.	Guidewire	Working Length (cm)	Working Channel Ø (mm)
<b>CY</b>	0.032"	230	2.8
<b>CY-DL</b>	0.032"	230	2.8

## Naso Jejunal Feeding Tube

To feed nutrition through mouth or nose with the tube leading into the jejunal

Tube Polyurethane



Art. No.	Size (Fr/CH)	Working Length (cm)
<b>NJFT-12120</b>	12	120
<b>NJFT-12</b>	12	155
<b>NJFT-14120</b>	14	120
<b>NJFT-14</b>	14	155
<b>NJFT-16120</b>	16	120
<b>NJFT-16</b>	16	155