

# TSUNAMI MEDICAL





## TVN

Targeting Needle 11G X 15 cm, Split Diamond and Bevel Tip

CE  
0123  
MADE IN ITALY





**New T-Shape Version**

**TVN**

Disposable Targeting Needle, 11 G X 15 cm composed by a rounded cannula and two stylets, with Split Diamond and Bevel Tip. The ergonomic handle has been designed to guarantee to ensure a safe grip and an optimal use of the hummer.

Available Codes	Description	Measures	Q.ty per box
TVN 11/15	1 Canula 1 Bevel Tip Mandrel 1 Split Diamond Mandrel	11 G X 15 cm	1

Represented by

NeuroKos SRL


mun. Chisinau, Republic of Moldova

tel. +373 7886 7886


e-mail: neurokos@yandex.ru

www. neurokos.md






**OSI POSTERIOR SPINE SYSTEM**  
**(5.5mm Rod System)**  
**Catalog Number:OPSS01**




**OSI MULTIAXIAL TRANSVERSE LINK**

Length	Reference Code
(33-40 MM)	ON0120-0035
(40-50 MM)	ON2120-0038
(50-60 MM)	ON0120-0050
(60-70 MM)	ON0120-0042




**OSI LATERAL CONNECTOR**

Length	Reference Code
10	ON0191-5510
15	OS0191-5515
20	OS0191-5520
25	OS0191-5525
30	OS0191-5530
40	OS0191-5540



**OSI DOMINO CONNECTOR**


Dimensions	Reference Code
21 x 14 x 11,5 mm	ON0195-6060



\* Optional Production (Delivery Time: 12 Weeks)  
 \*\*Visuals may vary

# X EXPANDABLE VERTEBRAL BODY REPLACEMENT

## Catalog Number: XXP01



X EXPANDABLE VERTEBRAL BODY REPLACEMENT		
Diameter	Lenght	Ref. Code
12	13-17MM – 8°	XP12H1317D8
	17-25MM – 8°	XP12H1725D8
	25-40MM – 8°	XP12H2540D8
	40-70MM – 8°	XP12H4070D8
14	13-17MM – 8°	X14H1317D8
	17-25MM – 8°	X14H1725D8
	25-40MM – 8°	X14H2540D8
	40-70MM – 8°	X14H4070D8
16	13-17MM – 8°	X16H1317D8
	17-25MM – 8°	X16H1725D8
	25-40MM – 8°	X16H2540D8
	40-70MM – 8°	X16H4070D8
18	13-17MM – 8°	X18H1317D8
	17-25MM – 8°	X18H1725D8
	25-40MM – 8°	X18H2540D8
	40-70MM – 8°	X18H4070D8


\* Optional Production (Delivery Time: 12 Weeks)

\*\*Visuals may vary



# BIA CERVICAL CAGE

## Catalog Number: BCC01



### BIA CERVICAL PEEK CAGE



Depth x Width	Height	Reference Code
12 x 14	H4	BOSCC412
	H5	BOSCC512
	H6	BOSCC612
	H7	BOSCC712
	H8	BOSCC812
13 x 15	H4	BOSCC413
	H5	BOSCC513
	H6	BOSCC613
	H7	BOSCC713
	H8	BOSCC813
14 x 16	H4	BOSCC414
	H5	BOSCC514
	H6	BOSCC614
	H7	BOSCC714
	H8	BOSCC814



\* Optional Production (Delivery Time: 12 Weeks)  
 \*\*Visuals may vary


# TLIF & PLIF CAGE

## Catalog Number:TPTFC01



### ZELOS PEEK PLIF CAGE (24 x 9 mm)

Degree	Height	Reference Code
0	H7	ZE240907
	H8	ZE240908
	H9	ZE240909
	H10	ZE240910
	H11	ZE240911
	H12	ZE240912
	H13	ZE240913
	H14	ZE240914
	H15	ZE240915
	H16	ZE240916
4	H7	ZE240907
	H8	ZE240908
	H9	ZE240909
	H10	ZE240910
	H11	ZE240911
	H12	ZE240912
	H13	ZE240913
	H14	ZE240914
	H15	ZE240915
	H16	ZE240916



\* Optional Production (Delivery Time: 12 Weeks)



\*\*Visuals may vary

# SENTINUS MESH CAGE

## Catalog Number:SCMC01

### SENTINUS CORPECTOMY MESH CAGE


Diameter	Height	Reference Code
10mm	60mm	SENTC1060
	80mm	SENTC1080
	100mm	SENTC10100
12mm	60mm	SENTC1260
	80mm	SENTC1280
	100mm	SENTC12100
14mm	60mm	SENTC1460
	80mm	SENTC1480
	100mm	SENTC14100
16mm	60mm	SENTL1660
	80mm	SENTL1680
	100mm	SENTL16100



\* Optional  
 \*\*Visuals may vary

# SENTINUS MESH CAGE

Catalog Number:SCMC01



---


	60mm	SENTL1860
18mm	80mm	SENTL1880
	100mm	SENTL18100

---

	60mm	SENTL2060
20mm	80mm	SENTL2080
	100mm	SENTL20100

---

\* Optional  
\*\*Visuals may vary






# PORTHOS POSTERIOR CERVICAL SPINE SYSTEM

## Catalog Number:PPCSS01




### PORTHOS POSTERIOR CERVICAL ROD (PRE-BEND)

Diameter (mm)	Length (mm)	Reference Code
3 - 4 - 3	240	PBR240




### PORTHOS DOUBLE DIAMETER ROD (HYBRID)

Diameter (mm)	Length (mm)	Reference Code
3 - 5,5	250	PB3250-55250
3 - 5,5	400	
3 - 6	400	



### PORTHOS TRANSVERSE CONNECTOR

Size	Reference Code
Small	PTC-S
Medium	PTC-M
Large	PTC-L




### PORTHOS DOMINO CONNECTOR

Size	Reference Code
Domino Connector S	PDOS
Domino Connector M/L	PDOML

\* Optional Production (Delivery Time: 12 Weeks)

\*\*Visuals may vary





8680571355116	OSI MONOAXIAL SCREW	Ø8,5X20	OS0155-8520	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355123	OSI MONOAXIAL SCREW	Ø8,5X25	OS0155-8525	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355130	OSI MONOAXIAL SCREW	Ø8,5X30	OS0155-8530	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355147	OSI MONOAXIAL SCREW	Ø8,5X35	OS0155-8535	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355154	OSI MONOAXIAL SCREW	Ø8,5X40	OS0155-8540	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355161	OSI MONOAXIAL SCREW	Ø8,5X45	OS0155-8545	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355178	OSI MONOAXIAL SCREW	Ø8,5X50	OS0155-8550	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355185	OSI MONOAXIAL SCREW	Ø8,5X55	OS0155-8555	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355192	OSI MONOAXIAL SCREW	Ø8,5X60	OS0155-8560	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE

MONOAXIAL REDUCTION SCREW



8680571355208	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X20	OS0161-4520	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355215	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X25	OS0161-4525	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355222	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X30	OS0161-4530	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355239	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X35	OS0161-4535	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355246	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X40	OS0161-4540	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355253	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X45	OS0161-4545	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355260	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X50	OS0161-4550	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355277	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X55	OS0161-4555	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355284	OSI MONOAXIAL REDUCTION SCREW	Ø4,5X60	OS0161-4560	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE




8680571355291	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X20	OS0161-5520	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355307	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X25	OS0161-5525	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355314	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X30	OS0161-5530	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355321	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X35	OS0161-5535	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355338	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X40	OS0161-5540	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355345	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X45	OS0161-5545	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355352	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X50	OS0161-5550	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355369	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X55	OS0161-5555	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355376	OSI MONOAXIAL REDUCTION SCREW	Ø5,5X60	OS0161-5560	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE



8680571355383	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X20	OS0161-6520	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355390	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X25	OS0161-6525	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355406	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X30	OS0161-6530	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355413	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X35	OS0161-6535	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355420	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X40	OS0161-6540	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355437	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X45	OS0161-6545	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355444	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X50	OS0161-6550	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355451	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X55	OS0161-6555	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE
8680571355468	OSI MONOAXIAL REDUCTION SCREW	Ø6,5X60	OS0161-6560	61325	Ti6Al4V Alloy	ISO 5832-3	NON-STERILE



**ARTHOS TITANIUM PLASMA COATED SCREW**  
**Catalog Number ART01**  
**(MONOAXIAL REDUCTABIL)**



ARTHOS TITANIUM PLASMA COATED SCREW		
Diameter	Lenght	Reference Code
4,5	25	ARTP4525*
	35	ARTP4535*
	40	ARTP4540*
	45	ARTP4545*
	50	ARTP4550*
5,5	30	ARTP5530
	35	ARTP5535
	40	ARTP5540
	45	ARTP5545
	50	ARTP5550
	55	ARTP5555
6,5	30	ARTP6530
	35	ARTP6535
	40	ARTP6540
	45	ARTP6545
	50	ARTP6550
	55	ARTP6555
7,5	30	ARTP7530
	35	ARTP7535
	40	ARTP7540
	45	ARTP7545
	50	ARTP7550
	55	ARTP7555

