


Technical Data Sheet

Model No.: N72A02

- **Product Name:** Surgical Drape With Adhesive Side
- **Classification:** Class Is, according to Annex VIII in MDR (EU) 2017/745
- **Description:** Surgical drape is made of nonwoven, without fenestration, with tape. It could be used for different operations according to its intended use. It is intended to be used to isolate the dust and bacteria, prevent the cross infection.
- **Intended use:** Used to isolate the dust and bacteria, prevent the cross infection.
- **Material:** PP+PE+PP 65gsm
- **Approvals:** CE, ISO 13485
- **Manufacturer:** Leboo Healthcare Products Limited
- **Country of origin:** China
- **Specification:** Detail see specification as below.



Material	Weight(gsm)	Size(cm)	Measurement Location(cm)				Color	Style
			A	B	C	D		
PP+PE+PP	65	100*75	100	75	100	5	Blue	Waterproof; Absorbent surface; Absorbtion 42g/m ² Celulose patient confort layer; Adhesive band on the edge of the drape; Latex free; Rosin free; EO Sterile

Design:


Check:

Approve:

Date:

●Packaging:

- (1) Carton Requirements: Five layers corrugated
- (2) Packing Quantity: 1pc/pouch, 100pcs/ctn
- (3) 5 Layers Carton(Paper)



●Storage:

- (1) Store in dry, clean conditions in original packaging.
- (2) Store away from direct sunlight, source of high temperature and solvent vapors.
- (3) Store with the temperature range -5°C to $+45^{\circ}\text{C}$ and with relative humidity below 65%.

●Shelf Life:

Shelf life is 36 months from the date of manufacture when stored as stated above.

●Cautions and disposal:

- (1) Do not contact with heavy oils, sparks or flames;
- (2) Do not use in the environment with conditions of excessive heat;
- (3) Replace if damaged, heavily contaminated or in accordance with local practice;
- (4) Hand and dispose of contaminated products with care and in accordance with national regulations.

CE 0197

CE Symbol



Medical Device



Keep away from sunlight




Do not reuse



Keep dry



Don't use when packing damaged
and consult instructions for use



Temperature limit



EO Sterile

LEBOO HEALTHCARE PRODUCTS LIMITED



Authorized Signature

Design:

Check:

Approve:

Date:



📍 H.O.: A-1506, Lar Valley International, No.168 Guang'anmen Wai St., Xicheng District, Beijing 100055, China.
☎ T: +86 10 63252699 📠 F: +86 10 63252899 ✉ E: leboo@leboo.com.cn 🌐 W: www.leboo.com.cn

Design:

Check:

Approve:

Date: