

MACH 1[™] Guide Catheter

Unique Polymer

• Enhances back-up support and curve retention

Large Lumen Design

- Offers the device compatibility needed for a multitude of treatment options and scenarios
- Enables catheter downsizing
- Provides improved dye flow for better visualization

Design and Performance Dynamics

- Softer shaft guide catheter designed for active users accustomed to deep seating
- Available in 6F (0.070"), 7F (0.081"), and 8F (0.091") sizes

PTFE liner

2 × 2 round braid pattern

Large lumen

Atraumatic, radiopaque tip

MACH 1[™] Guide Catheter

Left Femora	Annroa	ch Curves
Leitienioia	ΠΑΡΡΙΟά	

Amplatz Left (AL)
Amplatz Left Short Tip (AL ST)
CLS™ (Contralateral Left Support)
Femoral Curve Left (FCL)
Femoral Left (FL)
Femoral Left Short Tip (FL ST)
Judkins Left (JL)
Judkins Left Short Tip (JL ST)
KIESZ™ Left (KL)
Q-CURVE™ (Q)
VODA LEFT™ (VL)

Right Femoral Approach Curves
ALLRIGHT [™] (ART)
Amplatz Right (AR)
Femoral Curve Right (FCR)
Femoral Right (FR)
Femoral Right Short Tip (FR ST)
Hockey Stick (HS)
Judkins Right (JR)
Judkins Right Short Tip (JR ST)
KIESZ [™] Right (KR)
Multipurpose (MP)
VODA RIGHT [™] (VR)

Radial/Brachial Approach Curves
Brachial (BR)
KIMNY™
MANN™ Internal Mammary (MANN IM)
MUTA [™] Left (ML)
Radial
Radial Short Tip (Radial ST)
Radial Back-Up Left (RBL)
Radial Back-Up Right (RBR)
TIG

Additional Curves	
Internal Mammary (IM)	
Internal Mammary Cummings (IMC)	
Left Coronary Bypass (LCB)	
Right Coronary Bypass (RCB)	
Right Coronary Shepherd's Crook (RC SC)	

For a complete list of products, please contact your Boston Scientific representative.



www.bostonscientific-international.com

© 2016 Boston Scientific Corporation or its affiliates. All rights reserved. IC-374104-AB-EU JUL2016

All trademarks are the property of their respective owners.

Caution: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information contained herein is for use or distribution outside the US, France, and Japan only.