

Endoscopic Scissor

With monopolar dithermy
Used to cut and coagulate during
endoscopy mucosal resection

360° Rotatable



Art. No.	Sheath Ø (mm)	Working Channel Ø (mm)	Working Length (cm)
ESC-16	2.5	2.8	160
ESC-23	2.5	2.8	230

Class IIb device – for single use only

Spray Catheter

Used in conjunction with
an endoscope to spray
the gastrointestinal mucosa




Art. No.	Sheath Ø (mm)	Working Channel Ø (mm)	Working Length (cm)
WS-2418PB	2.4	2.8	180
WS-2423PB	2.4	2.8	230

Class IIa device – for single use only

Naso Jejunal Feeding Tube

To feed nutrition through mouth or nose
with the tube leading into the jejunal

Tube Polyurethane



Art. No.	Size (Fr/CH)	Working Length (cm)
NJFT-12120	12	120
NJFT-12	12	155
NJFT-14120	12	120
NJFT-14	12	155
NJFT-16120	12	120
NJFT-16	12	155

Class IIa devices – for single use only