



# ENDOASCULAR

Wide range of minimally invasive diagnostics devices dedicated to endovascular procedures enabling treatment of peripheral artery disease.



## STEEL GUIDEWIRE SUPER STIFF



- Extremely stiff wire
- Curved or straight tip configurations

Reference:	Diameter [inch]	Length [cm]	Type
P3290L	0.032"	90	STRAIGHT
P3580L			STRAIGHT
PJ3580L	0.035"	80	"J"
PJ3585L			"J"
P35100L		100	STRAIGHT
PJ35100L		150	"J"
PJ35150L			"J"
P3880L	0.038"	80	STRAIGHT
PJ3880L			"J"
P38100L		100	STRAIGHT
PJ38100L			"J"

When ordering use Ref. Number, e.g.: P35100L.

## STEEL GUIDEWIRE COATING PTFE SUPER STIFF



- PTFE coating
- Extremely stiff wire
- Curved or straight tip configurations

Reference:	Diameter [inch]	Length [cm]	Type
PL3590T	0.035"	90	STRAIGHT
PLJ3590T			"J"
PL35150T	0.035"	150	STRAIGHT
PLJ35150T			"J"
PL35180T		180	STRAIGHT
PLJ35180T			"J"
PL35260		260	STRAIGHT
PL35260T	"J"		
PLJ35260T	"J"		

When ordering use Ref. Number, e.g.: PL35150T.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

