

## EDTA Tube



Vacuum EDTA Tubes Vials (Lavender or Purple Top) contain EDTA as an anticoagulant and are widely used for determinations in EDTA whole blood clinical hematology tests, especially for a routine blood examination.

Vacuum Blood Collection Tube, EDTA Tube

Material: PET, Glass

Color: Lavender or Purple

| Cat. No. | Specification | Volume | Additive | Qty/case (Glass) | Qty/case (PET) |
|----------|---------------|--------|----------|------------------|----------------|
| 630801   | 13×75mm       | 1ml    | EDTAK2   | 1800             | 1800           |
| 630802   | 13×75mm       | 2ml    | EDTAK2   | 1800             | 1800           |
| 630803   | 13×75mm       | 3ml    | EDTAK2   | 1800             | 1800           |
| 630804   | 13×75mm       | 4ml    | EDTAK2   | 1800             | 1800           |
| 630805   | 13×75mm       | 5ml    | EDTAK2   | 1800             | 1800           |
| 630806   | 13×100mm      | 5ml    | EDTAK2   | 1200             | 1200           |
| 630807   | 13×100mm      | 6ml    | EDTAK2   | 1200             | 1200           |
| 630808   | 13×100mm      | 7ml    | EDTAK2   | 1200             | 1200           |

|        |          |      |         |      |      |
|--------|----------|------|---------|------|------|
| 630809 | 16×100mm | 10ml | EDTAK2  | 900  | 900  |
| 630901 | 13×75mm  | 1ml  | EDTAK3  | 1800 | 1800 |
| 630902 | 13×75mm  | 2ml  | EDTAK3  | 1800 | 1800 |
| 630903 | 13×75mm  | 3ml  | EDTAK3  | 1800 | 1800 |
| 630904 | 13×75mm  | 4ml  | EDTAK3  | 1800 | 1800 |
| 630905 | 13×75mm  | 5ml  | EDTAK3  | 1800 | 1800 |
| 630906 | 13×100mm | 5ml  | EDTAK3  | 1200 | 1200 |
| 630907 | 13×100mm | 6ml  | EDTAK3  | 1200 | 1200 |
| 630908 | 13×100mm | 7ml  | EDTAK3  | 1200 | 1200 |
| 630909 | 16×100mm | 10ml | EDTAK3  | 900  | 900  |
| 631001 | 13×75mm  | 2ml  | EDTANA2 | 1800 | 1800 |
| 631002 | 13×100mm | 10ml | EDTANA2 | 1200 | 1800 |