

## GROUNDING PADS



ADULT						
REF	DIMENSIONS	BAKING	ACTIVE SURFACE	REM/not REM	CABLE	CONNECTION
F7805	202x101 mm	FOAM	132 cm <sup>2</sup>	Not REM	No	Tab Connection
F7805W/V	202x101 mm	FOAM	132 cm <sup>2</sup>	Not REM	Yes 300 cm	Valleylab connector
F7805W/6.3	202x101 mm	FOAM	132 cm <sup>2</sup>	Not REM	Yes 300 cm	Ø 6.3 mm bipolar connector
F7820	202x101 mm	FOAM	118 cm <sup>2</sup>	REM	No	Tab Connection
F7820W/V	202x101 mm	FOAM	118 cm <sup>2</sup>	REM	Yes 300 cm	Valleylab connector
F7820W/6.3	202x101 mm	FOAM	118 cm <sup>2</sup>	REM	Yes 300 cm	Ø 6.3 mm bipolar connector



PAEDIATRIC						
REF	DIMENSIONS	BAKING	ACTIVE SURFACE	REM/not REM	CABLE	CONNECTION
F7805P	148x90 mm	FOAM	78 cm <sup>2</sup>	Not REM	No	Tab Connection
F7805PW/V	148x90 mm	FOAM	78 cm <sup>2</sup>	Not REM	Yes 300 cm	Valleylab connector
F7805PW/6.3	148x90 mm	FOAM	78 cm <sup>2</sup>	Not REM	Yes 300 cm	Ø 6.3 mm bipolar connector
F7820P	148x90 mm	FOAM	72 cm <sup>2</sup>	REM	No	Tab Connection
F7805PW/V	148x90 mm	FOAM	72 cm <sup>2</sup>	REM	Yes 300 cm	Valleylab connector

