

SOLUTIONS FOR ANEURYSM TREATMENT ECHELON™ MICROCATHETER TRUSTED IN ACTION

Large 0.017" ID¹
Low profile 0.028" OD¹

Fluoro-safe
marker band^{2,3}

Nitinol braid technology with progressive
picks-per-inch provides: ^{2,3}

- A balance of pushability and softness via graduated zones and transition points^{4,5}
- Optimal proximal pushability, stability, and support^{4,5}
- Kink resistance and minimal ovalization to maintain lumen durability for reliable device delivery^{4,5}

Steam-shapeable,
flexible tip with dual
marker bands ^{4,5}

PERFORMANCE AT EVERY TURN.


From navigation
through treatment.

Echelon™
Micro Catheter

Medtronic
Further, Together

Echelon™

Micro Catheter



ECHELON™ 10 MICRO CATHETER					
Product Catalog Number	Outer Diameter (F)	Inner Diameter (in)	Usable Length (cm)	Max. Guidewire (in)	Tip Configuration
105-5091-150	2.1>1.7	.017	150	0.014	Straight
145-5091-150	2.1>1.7	.017	150	0.014	45°
190-5091-150	2.1>1.7	.017	150	0.014	90°

ECHELON™ 14 MICRO CATHETER					
Product Catalog Number	Outer Diameter (F)	Inner Diameter (in)	Usable Length (cm)	Max. Guidewire (in)	Tip Configuration
105-5092-150	2.4>1.9	.017	150	0.014	Straight
145-5092-150	2.4>1.9	.017	150	0.014	45°
190-5092-150	2.4>1.9	.017	150	0.014	90°

ECHELON™ SYRINGE ADAPTER		
Product Catalog Number	Quantity	Units/Box
103-5095	1 Box	20

Medtronic

Europe
 Medtronic International Trading Sàrl.
 Route du Molliau 31
 Case postale
 CH-1131 Tolochenaz
 Tel: +41 (0)21 802 70 00
 Fax: +41 (0)21 802 79 00

REFERENCES:

1. DWGS105-5091-150. Echelon™ 10 catheter only
2. DWGS105-5092-150
3. DWGS105-5091-150
4. TR03-044
5. TR02-134A

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at medtronic.eu

UC201800632bEE © Medtronic 2020. All rights reserved. Printed in Europe.

medtronic.eu