

Indications, Safety, and Warnings

OVERVIEW

Ensure visibility before and after stent placement with the Visi-Pro™ balloon-expandable peripheral stent system for treatment of peripheral arterial disease (PAD). The stent's patented design balances strength and flexibility for proven patency as found in the VISIBILITY ILIAC study. Radiopaque markers on the stent ensure precise placement. Deliver over an 0.035" guidewire for use in the common or external iliac arteries or in the biliary system.

The Visi-Pro Balloon-Expandable Biliary Stent System is intended as a palliative treatment of malignant neoplasms in the biliary tree.

PRODUCT DETAILS

PREMIUM VISIBILITY

- Distinct 0.035" balloon-expandable stent with radiopaque markers ensures visibility during and after stent placement.

SOLID RADIAL STRENGTH AND FLEXIBILITY

- Patented design balances strength and flexibility.

PRECISE PLACEMENT

- Minimal dilation of healthy tissue is due to balloon to stent placement.
- Provides accurate catheter marker alignment to balloon taper.
- Progressive stent design delivers minimal shortening.

MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the product labeling supplied with each device or by calling technical support for cardiovascular services: +1 763-526-7890.

MODEL SPECIFICATIONS



Catheter Length		Stent Size		Balloon
Reference Number 80 cm	Reference Number 135 cm	Diameter (mm)	Length (mm)	Balloon Length (mm)
PXB35-05-12-080	-	5	12	15
PXB35-05-17-080	PXB35-05-17-135	5	17	20
PXB35-05-27-080	PXB35-05-27-135	5	27	30
PXB35-05-37-080	PXB35-05-37-135	5	37	40
PXB35-05-57-080	PXB35-05-57-135	5	57	60
PXB35-06-12-080	-	6	12	15
PXB35-06-17-080	PXB35-06-17-135	6	17	20
PXB35-06-27-080	PXB35-06-27-135	6	27	30
PXB35-06-37-080	PXB35-06-37-135	6	37	40
PXB35-06-57-080	PXB35-06-57-135	6	57	60
PXB35-07-12-080	-	7	12	15
PXB35-07-17-080	PXB35-07-17-135	7	17	20
PXB35-07-27-080	PXB35-07-27-135	7	27	30
PXB35-07-37-080	PXB35-07-37-135	7	37	40
PXB35-07-57-080	PXB35-07-57-135	7	57	60
PXB35-08-17-080	PXB35-08-17-135	8	17	20
PXB35-08-27-080	PXB35-08-27-135	8	27	30
PXB35-08-37-080	PXB35-08-37-135	8	37	40
PXB35-08-57-080	PXB35-08-57-135	8	57	60
PXB35-09-17-080	PXB35-09-17-135	9	17	20
PXB35-09-27-080	PXB35-09-27-135	9	27	30
PXB35-09-37-080	PXB35-09-37-135	9	37	40
PXB35-09-57-080	PXB35-09-57-135	9	57	60
PXB35-10-17-080	PXB35-10-17-135	10	17	20
PXB35-10-27-080	PXB35-10-27-135	10	27	30
PXB35-10-37-080	PXB35-10-37-135	10	37	40
PXB35-10-57-080	PXB35-10-57-135	10	57	60

* 6 F = 0.085" I.D.

Compliance Chart					
Visi-Pro Diameter (mm)	Inflation Pressure (atm)				
	8	9	10	11	12
5.0	5.00	5.06	5.16	5.22	5.28
6.0	6.00	6.11	6.22	6.31	6.39
7.0			7.00	7.09	7.17
8.0			8.00	8.15	8.26
9.0			9.00	9.15	9.28
10.0			10.00	10.11	10.21

 Diameter at nominal pressure
 Diameter at RBP

See product catalog for complete, detailed product information.

RELATED PAGES

- Clinical Outcomes

RELATED PRODUCTS

- Peripheral and Biliary Stent Systems
- Cardiovascular Products



PERIPHERAL CATALOG

Choose from a full portfolio of peripheral products for endovascular interventions – products that help you reach your PAD patient treatment goals.

DOWNLOAD CATALOG

CONTACT US

CardioVascular Lifeline Customer Support

+1 763-526-7890
Toll free +1 877-526-7890
rs.cstechsupport@medtronic.com

Orders

+1 763-514-8510
Toll free +1 800-716-6700
rs.cusvasorders@medtronic.com

24-hour Technical Support

Toll free +1 800-328-2518
+1 763-514-4000 (Worldwide)



Our Company


- Careers
- Contact & Support
- Locations
- ID Cards

Order Products

- Diabetes Supplies
- MRI Tools
- Recycling
- Product Manuals

Product Security

- Electromagnetic Compatibility for Cardiac Devices
- California Transparency in Supply Chains Act
- Information About Proposition 65 for California Customers
- Digital Millennium Copyright Act (DMCA) Notification



Indications, Safety, and Warnings

OVERVIEW

Ensure visibility before and after stent placement with the Visi-Pro™ balloon-expandable peripheral stent system for treatment of peripheral arterial disease (PAD). The stent's patented design balances strength and flexibility for proven patency as found in the VISIBILITY ILIAC study. Radiopaque markers on the stent ensure precise placement. Deliver over a 0.035" guidewire for use in the common or external iliac arteries or in the biliary system.

The Visi-Pro balloon-expandable biliary stent system is intended as a palliative treatment of malignant neoplasms in the biliary tree.

PRODUCT DETAILS

PREMIUM VISIBILITY

- Distinct 0.035" balloon-expandable stent with radiopaque markers ensures visibility during and after stent placement.

SOLID RADIAL STRENGTH AND FLEXIBILITY

- Patented design balances strength and flexibility.

PRECISE PLACEMENT

- Minimal dilation of healthy tissue is due to balloon to stent placement.
- Provides accurate catheter marker alignment to balloon taper.
- Progressive stent design delivers minimal shortening.

MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the product labeling supplied with each device or by calling technical support for cardiovascular services: +1 763-526-7890.

MODEL SPECIFICATIONS

Stent Size		Balloon	Compatibility	
Diameter (mm)	Length (mm)	Balloon Length (mm)	Recommended Introducer Sheath (F)	Guidewire (in)
5	12	15	6*	0.035"
5	17	20	6*	0.035"
5	27	30	6*	0.035"
5	37	40	6*	0.035"
5	57	60	6*	0.035"
6	12	15	6*	0.035"
6	17	20	6*	0.035"
6	27	30	6*	0.035"
6	37	40	6*	0.035"
6	57	60	6*	0.035"
7	12	15	6*	0.035"
7	17	20	6*	0.035"
7	27	30	6*	0.035"
7	37	40	6*	0.035"
7	57	60	6*	0.035"
8	17	20	6*	0.035"
8	27	30	6*	0.035"
8	37	40	6*	0.035"
8	57	60	6*	0.035"
9	17	20	7	0.035"
9	27	30	7	0.035"
9	37	40	7	0.035"
9	57	60	7	0.035"
10	17	20	7	0.035"
10	27	30	7	0.035"
10	37	40	7	0.035"
10	57	60	7	0.035"

* 6 F = 0.085" I.D.

Compliance Chart					
Visi-Pro Diameter (mm)	Inflation Pressure (atm)				
	8	9	10	11	12
5.0	5.00	5.06	5.16	5.22	5.28
6.0	6.00	6.11	6.22	6.31	6.39
7.0			7.00	7.09	7.17
8.0			8.00	8.15	8.26
9.0			9.00	9.15	9.28
10.0			10.00	10.11	10.21

	Diameter at nominal pressure
	Diameter at RBP

See product catalog for complete, detailed product information.

RELATED PAGES

- Clinical Outcomes

RELATED PRODUCTS

- Peripheral and Biliary Stent Systems
- Cardiovascular Products



PERIPHERAL CATALOG

Choose from a full portfolio of peripheral products for endovascular interventions – products that help you reach your PAD patient treatment goals.

[DOWNLOAD CATALOG](#)

CONTACT US

CardioVascular Lifeline Customer Support

+1 763-526-7890
Toll free +1 877-526-7890
rs.cstechsupport@medtronic.com

Orders

+1 763-514-8510
Toll free +1 800-716-6700
rs.cusvasorders@medtronic.com

24-hour Technical Support

Toll free +1 800-328-2518
+1 763-514-4000 (Worldwide)



Our Company

- Careers
- Contact & Support
- Locations
- ID Cards

Order Products

- Diabetes Supplies
- MRI Tools
- Recycling
- Product Manuals

Product Security

- Electromagnetic Compatibility for Cardiac Devices
- California Transparency in Supply Chains Act
- Information About Proposition 65 for California Customers
- Digital Millennium Copyright Act (DMCA) Notification