




Biliary Stone Extraction Balloon Catheter



Diameter of inflated Balloon 16mm Balloon max. Air Inflation 4.0cc Balloon max. Fluid Injection 3.0ml Recommended Guidewire 0.035"

Used for endoscopic


during procedure



Biliary Nitinol Stent Self Expandable

Used for stenosis or obstructed biliary tract caused by benign or malignant changes

With Introducer System
Radiotherapy compatible
Uncovered stent
Expanded Stent Diameter – 10mm
Compatible with instrument channel size 3.2mm



Art. No.	Stent Length (mm)	Stent Ø (mm)	Guidewire
S-MBS-6	60	10	0.038"
S-MBS-8	80	10	0.038"
S-MBS-10	100	10	0.038"
S-MBS-12	120	10	0.038"