








Eyepiece




**CFI 10X
MAK10110**




**CFI UW 10X
MAK30105**




**CFI 12.5X
MAK10120**




**CFI 15X
MAK10150**




**C-W 10XB
MAK30102**



**C-W 15X
MAK30155**



**C-W 20X
MMK30205**




**C-W 30X
MMK30300**


Code.No	Eyepieces	Field number	Weight (g)*
MAK10110	CFI 10X	22	75
MAK30105	CFI UW 10X	25	100
MAK10120	CFI 12.5X	16	63
MAK10150	CFI 15X	14.5	48
MMK30102	C-W 10XB	22	75
MMK30155	C-W 15X	16	50
MMK30205	C-W 20X	12.5	60
MMK30300	C-W 30X	7	70

*Approximately


Eyepiece Tube





**C-TB Binocular Tube
MBB92106**



**C-TF Trinocular Tube F
MBB93106**



**LV-TI3 Trinocular Tube
MBB63425**



**C-TT Trinocular Tube T
MBB93115**

Code.No	Eyepiece Tubes	Field number	Observation Inclination	Observation Image	Interpupillary adjustment	Intensity split ratio (Eyepiece:Port)	Weight (g)*
MBB92106	C-TB Binocular Tube	22	25°	Inverted image	50 to 75mm		1400
MBB63425	LV-TI3 Trinocular Tube	25	20°	Upright image	50 to 75mm	100:0 or 0:100 two-way switching	1800
MBB93106	C-TF Trinocular Tube F	25	25°	Inverted image	50 to 75mm	100:0 or 0:100 two-way switching	2000
MBB93115	C-TT Trinocular Tube T	25	25°	Inverted image	50 to 75mm	100:0, 20:80 or 0:100 three-way switching	2500
MBB93800	C-TE2 Ergonomic Binocular Tube	22	10°-30°	Inverted image	50 to 75mm	100:0 or 50:50 two-way switching	2000

*Approximately

Enter the “Microscopy University” on the web and discover a whole new world.

MicroscopyU

<https://www.microscopyu.com>

Nikon’s International Small World Photomicrography Competition



<https://www.nikonsmallworld.com>

N.B. Export of the products* in this brochure is controlled under the Japanese Foreign Exchange and Foreign Trade Law. Appropriate export procedure shall be required in case of export from Japan.

*Products: Hardware and its technical information (including software)

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. February 2026 ©2026 NIKON CORPORATION



WARNING

TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.



NIKON CORPORATION

Head office/Manufacturer
1-5-20, Nishioji, Shinagawa-ku, Tokyo 140-8601, Japan
<https://www.healthcare.nikon.com/en/>

Nikon Instruments Inc.

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A.
phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only)
fax: +1-631-547-0299
<https://www.microscope.healthcare.nikon.com/>

Nikon Europe B.V.

Stroombaan 14, 1181 VX Amstelveen, The Netherlands
phone: +31-20-7099-000
https://www.microscope.healthcare.nikon.com/en_EU/

Nikon Precision (Shanghai) Co., Ltd.

CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060
(Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026
(Guangzhou branch) phone: +86-20-3882-0550 fax: +86-20-3882-0580
<https://www.nikon-precision.com.cn/>

Nikon Canada Inc.

CANADA phone: +1-905-625-9910 fax: +1-905-602-9953

Nikon France, Succursale de Nikon Europe B.V.

FRANCE phone: +33-1-4516-4516

Nikon Deutschland, Zweigniederlassung der Nikon Europe B.V.

GERMANY phone: +49-211-9414-888

Nikon Italy, Branch of Nikon Europe B.V.

ITALY phone: +39-055-300-9601

Nikon Europe B.V., Amstelveen, Zweigniederlassung Schweiz (Egg/ZH)

SWITZERLAND phone: +41-43-277-2867

Nikon UK, Branch of Nikon Europe B.V.

UNITED KINGDOM phone: +44-208-247-1717

Nikon Österreich, Zweigniederlassung der Nikon Europe B.V.

AUSTRIA phone: +43-1-972-6111

Nikon Singapore Pte. Ltd.

SINGAPORE phone: +65-6559-3651 fax: +65-6559-3668

Nikon Australia Pty Ltd

AUSTRALIA phone: +61-2-8767-6900

Nikon Instruments Korea Co., Ltd.

KOREA phone: +82-2-6288-1900 fax: +82-2-555-4415

Nikon India Private Limited

INDIA phone: +91-124-4688-500