






27005 FA



HOPKINS® Telescope 12°,
enlarged view, diameter 4 mm, length 30 cm, autoclavable, fiber optic light transmission incorporated,
color code: black



495 GO




Adaptor,
to connect KARL STORZ HOPKINS® telescopes to Olympus/Winter & Ibe light cable, with attached lens


Image View – 495 GO



27068 D




SACHSE Urethrotome Sheath,
21 Fr., with channel for FILIFORM bougies




27068 DO



Obturator,
for Urethrotome Sheath 21 Fr., 27068 D



27050 E




Working Element,



(also for use with optical urethrotome), motion by means of a spring. The thumb support is movable. In rest position the electrode is inside the sheath.




27069 M





LUDVIK Knife,
straight, serrated, not to be used with HF current



27069 K



SACHSE Knife,
straight, not to be used with HF current



27069 L



Cold Knife,
round, not to be used with HF current

