



DİLAMED
Tıbbi Ürünler A.Ş

BROCHURE TECNICAL DATA SHEET



info@dilamedtibbi.com
sales@dilamedtibbi.com



+90 232 281 11 25



www.dilamedtibbi.com



DİLAMED
Tıbbi Ürünler A.Ş.

DİLEX SURGICAL SCRUB BRUSH IMPREGNATED WITH 7.5% POVIDONE-IODINE

Features

- It is a broad spectrum skin antisepsis containing povidone iodine.
- The iodine contained in the polyvinyl pyrrolidone iodine complex is slowly released and provides an antiseptic effect.
- Foams easily.
- Can be cleaned with water without leaving stains.
- It has a dark reddish brown appearance with a characteristic odour.
- Easy to use with ergonomic brush design.
- Nail interiors are thoroughly cleaned with soft bristles that do not irritate the skin.
- For professional users (healthcare personnel).
- For single use only. Dispose as medical waste after use.
- Do not use if the packaging is opened or damaged.
- Does not contain latex or phthalates.
- Options with or without nail pick are available.
- For external use only.
- Biocompatible.



PRODUCT NAME	REF NO	BARCODE	PACKAGE
SURGICAL SCRUB BRUSH IMPREGNATED WITH 7.5% POVIDONE-IODINE	DLM009-1	8682835111926	INNER BOX: 40 PCS
			OUTER BOX: 400 PCS

Quality Standard

Manufactured according to ISO 13485:2016
Quality Management Standards.
Biocidal license from the Turkish Ministry of Health.



info@dilamedtibbi.com
sales@dilamedtibbi.com



+90 232 281 11 25



www.dilamedtibbi.com



DİLAMED
Tıbbi Ürünler A.Ş

DİLEX SURGICAL SCRUB BRUSH IMPREGNATED WITH 4% CHLORHEXIDINE GLUCONATE

Features

- It is a broad spectrum skin antiseptis containing 4% Chlorhexidine Gluconate.
- Foams easily.
- Can be cleaned with water without leaving stains.
- It has a dark pink appearance with a characteristic odour.
- Ph is compatible with skin.
- Easy to use with ergonomic brush design.
- Nail interiors are thoroughly cleaned with soft bristles that do not irritate the skin.
- For professional users (healthcare personnel).
- For single use only. Dispose as medical waste after use.
- Do not use if the packaging is opened or damaged.
- Does not contain latex or phthalates.
- Options with or without nail pick are available.
- For external use only.
- Biocompatible.



PRODUCT NAME	REF NO	BARCODE	PACKAGE
SURGICAL SCRUB BRUSH IMPREGNATED WITH 4% CHLORHEXIDINE GLUCONATE	DLM0010-1	8682835112572	INNER BOX: 40 PCS
			OUTER BOX: 400 PCS

Quality Standard

Manufactured according to ISO 13485:2016
Quality Management Standards.
Biocidal license from the Turkish Ministry of Health.



info@dilamedtibbi.com
sales@dilamedtibbi.com



+90 232 281 11 25



www.dilamedtibbi.com



DILAMED
Tıbbi Ürünler A.Ş

Surgical Scrub Brush Description

Scrub Brushes consist of a polyethylene brush and sponge. In line with customer requests, a nail pick made of polystyrene can be added to the packaging, which is used to clean in and around the nail.

Indications

Scrub Brushes are used for cleaning and disinfecting forearms, hands and nails before and after surgical intervention, dressing and examination.

Physical Characteristics

Scrub Brushes consist of a plastic brush on one side and a sponge on the other. It is produced according to the following dimensions: 8 cm (Length) x 5 cm (Width) x 4 cm (Height). The brush is made of low density polyethylene. The sponge is 22±2 kg/ m³ dense and is polyurethane based. The sponge and brush are glued together via a hotmelt adhesive.

Indications

Scrub Brushes are used for cleaning and disinfecting forearms, hands and nails before and after surgical intervention, dressing and examination.

Use of the Device

Surgical Handwashing: Washing in 3 Steps

Step 1 – Preliminary Wash

- Wet hands, wrists and forearms.
- Apply a sufficient amount of soap/ antiseptic agent and wash from the fingertips to the elbow for a period of 1 minute.
- Hands should be held in an upright position, higher than the elbows at all times.
- Rinse hands, wrists and forearms with plenty of running water.

Step 2 – First Wash/Scrub

- Repeat the same wash from Step 1.
- Scrub each hand and nails individually while counting for 30 seconds per hand; 1 minute in total.
- If the brush is impregnated, the sponge can be used to distribute the antiseptic agent over the hands, wrists and forearms in a downward motion.
- Using the brush bristles, scrub the nails of one hand, using a downward motion to avoid spray contamination (30 seconds).
- Rinse the brush under running water and repeat for the second hand (30 seconds).
- Using the provided nail pick, clean underneath and around the edges of the nails.
- Properly discard of the used brush.
- Rinse hands, wrists and forearms with plenty of running water.

Step 3 – Second Wash

- Repeat the same wash from Step 1 for 2 minutes.
- Dry off from the fingertips to the elbow with a single-use sterile towel.
- Hold hands in an upright position.
- Maintain this position while donning a gown.



info@dilamedtibbi.com
sales@dilamedtibbi.com



+90 232 281 11 25



www.dilamedtibbi.com



DİLAMED
Tıbbi Ürünler A.Ş

Scrub Brush components are specified in Table 1.

Name of the product	Components	Raw materials
Surgical Scrub Brush	Brush	Polyethylene 122 19T Petkim
	Sponge	Polyurethane
	File	Polystyrene
	Packaging bottom foil	PE/PA Film
	Packaging top foil	Medical Paper
	Box	Cardboard
	Parcel	Cardboard (C+B Wave)

Table 1. Nail Brush components

SDS (Safety Data Sheets) and Technical Information Documents of raw materials are in Annex 1.

Surgical Scrub Brush production processes;

- Plastic Injection
- Sponge Cutting
- Gluing
- Packaging (packaging closing, boxing)

Specifications

- Production at DİLAMED is carried out in accordance with TS EN ISO 13485 Quality Management System standard. Certification was carried out by QRS Quality Registration.
- Scrub Brushes are suitable for use by professional users (medical personnel) only. It is not suitable for use by untrained people.
- The integrity of the packaging should be checked before use.
- The product is disposable and should not be reused.
- It should not be used for purposes other than the declared purpose.
- It is not suitable and not recommended for use on third parties other than for medical purposes.
- Latex and Phthalate free.
- Symbols regarding usage warnings of the product are available on the label.

User Information Storage Conditions:

Store in a dry and cool place at room temperature below 25 °C. Do not expose to direct sunlight for long periods of time.

Packaging

Box Qty : 400 pcs
Inner Bag Qty : 40 pcs



info@dilamedtibbi.com
sales@dilamedtibbi.com



+90 232 281 11 25



www.dilamedtibbi.com



	COMPANY LOGO
Surgical Scrub Brush	Product Name
	Lot Number
	Catalog Reference Number
	Instruct
	Use-by Date
	Single Use
	Keep Out of Direct Sunlight
	Do not use if the package is open or damaged
	Manufacturer's Information DİLAMED Tıbbi Ürünler San. ve Tic. A.Ş. Yazıbaşı Mah. Fatma Tonkuş Cad. No: 24 Torbalı - İZMİR Tel: +90 232 281 11 25 sales@dilamedtibbi.com info@dilamedtibbi.com www.dilamedtibbi.com
	Max Temperature Storage



info@dilamedtibbi.com
sales@dilamedtibbi.com



+90 232 281 11 25



www.dilamedtibbi.com