

Pharmapore® Sterile

Post-operative Dressings



Description:

Pharmapore® Family is a self-adhesive dressing having an absorbent pad with a non-adherent wound contact layer, placed centrally on different backing materials coated with hypoallergenic adhesive, to allow fixation without secondary dressing.

Post-operative dressings:

a- **Pharmapore®**: The dressing consists of a pad placed centrally on a non-woven backing material for patients with sensitive skin.

Pharmapore® Ultra: The pad is placed centrally on a nonwoven backing material laminated to breathable, virus-proof, bacteria-proof and waterproof polyurethane film.

b- **Pharmapore® PU**: The pad is placed centrally on transparent, breathable, virus-proof, bacteria-proof and waterproof polyurethane backing material.

Pharmapore® PU frame style: same as **Pharmapore® PU**, but with frame style delivery system for easier application.

Features:

- Non-adherent wound contact layer: allows minimal trauma upon removal.
- Breathable: prevents maceration and helps to maintain a moist wound environment.
- Soft and Conformable: as it conforms to skin surface even the difficult areas like underarm and shoulders.

In case of Pharmapore® PU variants:

- Decreases risk of infection: as PU film is waterproof, bacteria-proof and virus-proof.
- Transparent: allows visual inspection of the area around the wound.

In case of Pharmapore® Ultra:

- Additional flexibility and strength: due to the unique combination of non-woven with PU.
- Colored PU film: tan color offers a cosmetically pleasant appearance.

Intended Purpose:

Primary dressings indicated for low to moderate exuding wounds such as minor cuts, abrasions, lacerations and puncture sites. Besides, they can be used post-operatively to protect surgical wounds and incisions.

- Amria - Alexandria - Egypt

Tel. : 002 03 4500264 Fax : 002 03 4500263

- Kafr-El-Zayat P.O.Box 38 Egypt

Tel. : 002 040 2545643 Fax : 002 040 2545645

E-mail : sales@pharmaplast-online.com

 **Pharmaplast**

www.pharmaplast-online.com

Pharmapore® Sterile

Post-operative Dressings

Pharmapore®	Size	Pieces/Box	Boxes/Case	Case Dimensions
PORE5070	5x7cm	200	8	44x30x30 cm
PORE6070	6x7cm	200	8	44x30x30 cm
PORE6080	6x8cm	200	8	44x30x30 cm
PORE5090	5x9cm	200	8	44x30x30 cm
PORE6083	6x8.3cm	200	8	44x30x30 cm
PORE50100	5x10cm	200	8	44x30x30 cm
PORE6082	6x8.25cm	200	8	44x30x30 cm
PORE6085	6x8.5cm	200	8	44x30x30 cm
PORE7080	7x8cm	200	8	44x30x30 cm
PORE6580	6.5x8cm	200	8	44x30x30 cm
PORE7575	7.5x7.5cm	200	8	44x30x30 cm
PORE6090	6x9cm	200	8	44x30x30 cm
PORE60100	6x10cm	200	8	44x30x30 cm
PORE5082	5x8.25cm	200	8	44x30x30 cm
PORE75100	7.5x10cm	100	8	44x30x30 cm
PORE80100	8x10cm	100	8	44x30x30 cm
PORE82100	8.25x10cm	100	8	44x30x30 cm
PORE100100	10x10cm	100	8	44x30x30 cm
PORE70120	7x12cm	60	8	44x30x30 cm
PORE80120	8x12cm	60	8	44x30x30 cm
PORE100120	10x12cm	60	8	44x30x30 cm
PORE83120	8.3x12cm	60	8	44x30x30 cm
PORE90150	9x15cm	60 10	8 24	43X38x25 cm 44x30x30 cm
PORE80150	8x15cm	60	8	44x30x30 cm
PORE82150	8.25x15cm	60	8	44x30x30 cm
PORE100140	10x14cm	60	8	44x30x30 cm
PORE100150	10x15cm	60	8	44x30x30 cm
PORE60150	6x15cm	60	8	
PORE83180	8.3x18cm	60	8	44x30x30 cm

Pharmapore® Sterile

Post-operative Dressings

PORE82200	8.25x20cm	40	8	44x30x30 cm
PORE90200	9x20cm	40	8	53x27x33 cm
		10	24	44x30x30 cm
PORE100200	10x20cm	40	8	44x30x30 cm
PORE83240	8.3x24cm	40	8	44x30x30 cm
PORE82250	8.25x25cm	40	8	44x30x30 cm
PORE90250	9x25cm	40	8	44x30x30 cm
PORE100250	10x25cm	40	8	44x30x30 cm
PORE82300	8.25x30cm	40	8	44x30x30 cm
PORE90300	9x30cm	40	8	44x30x30 cm
PORE100300	10x30cm	40	8	44x30x30 cm
PORE90350	9x35cm	40	8	44x30x30 cm
PORE82350	8.25x35cm	40	8	44x30x30 cm
PORE100350	10x35cm	40	8	44x30x30 cm
PORE83415	8.3x41.5cm	40	8	44x30x30 cm
PORE100160	10x16cm	60	8	44x30x30 cm
PORE100210	10x21cm	40	8	44x30x30 cm
PORE100260	10x26cm	40	8	44x30x30 cm
PORE100400	10x40cm	40	8	44x30x30 cm
PORE120160	12x16cm	60	8	44x30x30 cm
PORE150150	15x15cm	10	24	43x38x25 cm
PORE200200	20x20cm	10	24	53x27x33 cm
PORE2572	2.5x7.2cm	400	8	44x30x30 cm
PORE83100	8.3x10cm	100	8	44x30x30 cm
PORE5072	5x7.2 cm	200	8	56x30x30 cm
		50	60	44x30x30 cm
PORE5075	5x7.5cm	200	8	44x30x30 cm
PORE70100	7x10cm	100	8	44x30x30 cm
PORE6075	6x7.5cm	200	8	44x30x30 cm
PORE85155	8.5x15.5cm	60	8	44x30x30 cm
PORE85150	8.5x15cm	60	8	44x30x30 cm
PORE82120	8.2x12cm	60	8	44x30x30 cm
PORE8595	8.5x9.5cm	100	8	44x30x30 cm
PORE90100	9x10cm	100	8	44x30x30 cm
PORE125125	12.5x12.5cm	10	24	44x30x30 cm

Pharmapore® Sterile

Post-operative Dressings

Pharmapore® Ultra	Size	Pieces/Box	Boxes/Case	Case Dimensions
POREU5070	5x7cm	200	8	44x30x30 cm
POREU6070	6x7cm	200	8	44x30x30 cm
POREU6080	6x8cm	200	8	44x30x30 cm
		150	8	44x30x30 cm
POREU5090	5x9cm	200	8	44x30x30 cm
POREU6083	6x8.3cm	200	8	44x30x30 cm
POREU50100	5x10cm	200	8	44x30x30 cm
POREU6082	6x8.25cm	200	8	44x30x30 cm
POREU6085	6x8.5cm	200	8	44x30x30 cm
POREU7080	7x8cm	200	8	44x30x30 cm
POREU6580	6.5x8cm	200	8	44x30x30 cm
POREU7575	7.5x7.5cm	200	8	44x30x30 cm
POREU6090	6x9cm	200	8	44x30x30 cm
POREU60100	6x10cm	200	8	44x30x30 cm
POREU5082	5x8.25cm	200	8	44x30x30 cm
POREU75100	7.5x10cm	100	8	44x30x30 cm
POREU80100	8x10cm	100	8	44x30x30 cm
POREU82100	8.25x10cm	100	8	44x30x30 cm
POREU100100	10x10cm	100	8	44x30x30 cm
POREU70120	7x12cm	60	8	44x30x30 cm
POREU80120	8x12cm	60	8	44x30x30 cm
POREU100120	10x12cm	60	8	44x30x30 cm
POREU83120	8.3x12cm	60	8	44x30x30 cm
POREU90150	9x15cm	60	8	44x30x30 cm
		10	24	43X38x25 cm
POREU80150	8x15cm	60	8	44x30x30 cm
POREU82150	8.25x15cm	60	8	44x30x30 cm
POREU100140	10x14cm	60	8	44x30x30 cm
POREU100150	10x15cm	60	8	44x30x30 cm
POREU83180	8.3x18cm	60	8	44x30x30 cm
POREU82200	8.25x20cm	40	8	44x30x30 cm
POREU90200	9x20cm	40	8	53x27x33 cm
		10	24	44x30x30 cm
POREU100200	10x20cm	40	8	44x30x30 cm
POREU83240	8.3x24cm	40	8	44x30x30 cm
POREU82250	8.25x25cm	40	8	44x30x30 cm
POREU90250	9x25cm	40	8	44x30x30 cm
POREU100250	10x25cm	40	8	44x30x30 cm

Pharmapore® Sterile

Post-operative Dressings

POREU82300	8.25x30cm	40	8	44x30x30 cm
POREU90300	9x30cm	40	8	44x30x30 cm
POREU100300	10x30cm	40	8	44x30x30 cm
POREU90350	9x35cm	40	8	44x30x30 cm
POREU82350	8.25x35cm	40	8	44x30x30 cm
POREU100350	10x35cm	40	8	44x30x30 cm
POREU83415	8.3x41.5cm	40	8	44x30x30 cm
POREU100160	10x16cm	60	8	44x30x30 cm
POREU100210	10x21cm	40	8	44x30x30 cm
POREU100260	10x26cm	40	8	44x30x30 cm
POREU100400	10x40cm	40	8	44x30x30 cm
POREU120160	12x16cm	60	8	44x30x30 cm
POREU150150	15x15cm	10	24	43X38x25 cm
POREU200200	20x20cm	10	24	53x27x33 cm
POREU2572	2.5x7.2cm	400	8	44x30x30 cm
POREU83100	8.3x10cm	100	8	44x30x30 cm
POREU5072	5x7.2 cm	50	60	56x30x30 cm
		200	8	44x30x30 cm
POREU5075	5x7.5cm	200	8	44x30x30 cm
POREU70100	7x10cm	100	8	44x30x30 cm
POREU6075	6x7.5cm	200	8	44x30x30 cm
POREU85155	8.5x15.5cm	60	8	44x30x30 cm
POREU85150	8.5x15cm	60	8	44x30x30 cm
POREU7090	7x9cm	100	8	44x30x30 cm
POREU8595	8.5x9.5cm	100	8	44x30x30 cm
POREU90100	9x10cm	100	8	44x30x30 cm
POREU125125	12.5x12.5cm	10	24	44x30x30 cm

Pharmapore® Sterile

Post-operative Dressings

Pharmapore® PU	Size	Pieces/Box	Boxes/Case	Case Dimensions
POREPU5070	5x7cm	200	8	44x30x30 cm
POREPU6070	6x7cm	200	8	44x30x30 cm
POREPU6080	6x8cm	200	8	44x30x30 cm
POREPU5090	5x9cm	200	8	44x30x30 cm
POREPU6083	6x8.3cm	200	8	44x30x30 cm
POREPU50100	5x10cm	200	8	44x30x30 cm
POREPU6082	6x8.25cm	200	8	44x30x30 cm
POREPU6085	6x8.5cm	200	8	44x30x30 cm
POREPU7080	7x8cm	200	8	44x30x30 cm
POREPU6580	6.5x8cm	200	8	44x30x30 cm
POREPU7575	7.5x7.5cm	200	8	44x30x30 cm
POREPU6090	6x9cm	200	8	44x30x30 cm
POREPU60100	6x10cm	200	8	44x30x30 cm
POREPU5082	5x8.25cm	200	8	44x30x30 cm
POREPU75100	7.5x10cm	100	8	44x30x30 cm
POREPU80100	8x10cm	100	8	44x30x30 cm
POREPU82100	8.25x10cm	100	8	44x30x30 cm
POREPU100100	10x10cm	100	8	44x30x30 cm
POREPU70120	7x12cm	60	8	44x30x30 cm
POREPU80120	8x12cm	60	8	44x30x30 cm
POREPU100120	10x12cm	60	8	44x30x30 cm
POREPU83120	8.3x12cm	60	8	44x30x30 cm
POREPU90150	9x15cm	60	8	44x30x30 cm
		10	24	43X38x25 cm
POREPU80150	8x15cm	60	8	44x30x30 cm
POREPU82150	8.25x15cm	60	8	44x30x30 cm
POREPU100140	10x14cm	60	8	44x30x30 cm
POREPU100150	10x15cm	60	8	44x30x30 cm
POREPU83180	8.3x18cm	60	8	44x30x30 cm
POREPU82200	8.25x20cm	40	8	44x30x30 cm
POREPU90200	9x20cm	40	8	53x27x33 cm
		10	24	53x27x33 cm
POREPU100200	10x20cm	40	8	44x30x30 cm
POREPU83240	8.3x24cm	40	8	44x30x30 cm

Pharmapore® Sterile

Post-operative Dressings

POREPU82250	8.25x25cm	40	8	44x30x30 cm
POREPU90250	9x25cm	40	8	44x30x30 cm
POREPU100250	10x25cm	40	8	44x30x30 cm
POREPU82300	8.25x30cm	40	8	44x30x30 cm
POREPU90300	9x30cm	40	8	44x30x30 cm
POREPU100300	10x30cm	40	8	44x30x30 cm
POREPU90350	9x35cm	40	8	44x30x30 cm
POREPU82350	8.25x35cm	40	8	44x30x30 cm
POREPU100350	10x35cm	40	8	44x30x30 cm
POREPU83415	8.3x41.5cm	40	8	44x30x30 cm
POREPU100160	10x16cm	60	8	44x30x30 cm
POREPU100210	10x21cm	40	8	44x30x30 cm
POREPU100260	10x26cm	40	8	44x30x30 cm
POREPU100400	10x40cm	40	8	44x30x30 cm
POREPU120160	12x16cm	60	8	44x30x30 cm
POREPU150150	15x15cm	10	24	44x30x30 cm
POREPU200200	20x20cm	10	24	43X38x25 cm
POREPU2572	2.5x7.2cm	400	8	53x27x33 cm
POREPU83100	8.3x10cm	100	8	44x30x30 cm
POREPU5072	5x7.2cm	50	60	44x30x30 cm
		200	8	44x30x30 cm
POREPU5075	5x7.5cm	200	8	44x30x30 cm
POREPU70100	7x10cm	100	8	44x30x30 cm
POREPU6075	6x7.5cm	200	8	44x30x30 cm
POREPU85155	8.5x15.5cm	60	8	44x30x30 cm
POREPU85150	8.5x15cm	60	8	44x30x30 cm
POREPU82120	8.2x12cm	60	8	44x30x30 cm
POREPU8595	8.5x9.5cm	100	8	44x30x30 cm
POREPU90100	9x10cm	100	8	44x30x30 cm
POREPU125125	12.5x12.5cm	10	24	38x23x26cm

Pharmapore® Sterile

Post-operative Dressings

Pharmapore® PU Frame Style	Size	Pieces/Box	Boxes/Case	Case Dimensions
POREPUFS5070	5x7cm	200	8	44x30x30 cm
POREPUFS6070	6x7cm	200	8	44x30x30 cm
POREPUFS6080	6x8cm	200	8	44x30x30 cm
		150	8	44x30x30 cm
POREPUFS5090	5x9cm	200	8	44x30x30 cm
POREPUFS6083	6x8.3cm	200	8	44x30x30 cm
POREPUFS50100	5x10cm	200	8	44x30x30 cm
POREPUFS6082	6x8.25cm	200	8	44x30x30 cm
POREPUFS6085	6x8.5cm	200	8	44x30x30 cm
POREPUFS7080	7x8cm	200	8	44x30x30 cm
POREPUFS6580	6.5x8cm	200	8	44x30x30 cm
POREPUFS7575	7.5x7.5cm	200	8	44x30x30 cm
POREPUFS6090	6x9cm	200	8	44x30x30 cm
POREPUFS60100	6x10cm	200	8	44x30x30 cm
POREPUFS5082	5x8.25cm	200	8	44x30x30 cm
POREPUFS75100	7.5x10cm	100	8	44x30x30 cm
POREPUFS80100	8x10cm	100	8	44x30x30 cm
POREPUFS82100	8.25x10cm	100	8	44x30x30 cm
POREPUFS100100	10x10cm	100	8	44x30x30 cm
POREPUFS70120	7x12cm	60	8	44x30x30 cm
POREPUFS80120	8x12cm	60	8	44x30x30 cm
POREPUFS100120	10x12cm	60	8	44x30x30 cm
POREPUFS83120	8.3x12cm	60	8	44x30x30 cm
POREPUFS90150	9x15cm	60	8	44x30x30 cm
		10	24	43X38x25 cm
POREPUFS80150	8x15cm	60	8	44x30x30 cm
POREPUFS82150	8.25x15cm	60	8	44x30x30 cm
POREPUFS100140	10x14cm	60	8	44x30x30 cm
POREPUFS100150	10x15cm	60	8	44x30x30 cm
POREPUFS83180	8.3x18cm	60	8	44x30x30 cm
POREPUFS82200	8.25x20cm	40	8	44x30x30 cm

Pharmapore® Sterile Post-operative Dressings

POREPUFS90200	9x20cm	40	8	44x30x30 cm
		10	24	53x27x33 cm
POREPUFS100200	10x20cm	40	8	44x30x30 cm
POREPUFS83240	8.3x24cm	40	8	44x30x30 cm
POREPUFS82250	8.25x25cm	40	8	44x30x30 cm
POREPUFS90250	9x25cm	40	8	44x30x30 cm
POREPUFS100250	10x25cm	40	8	44x30x30 cm
POREPUFS82300	8.25x30cm	40	8	44x30x30 cm
POREPUFS90300	9x30cm	40	8	44x30x30 cm
POREPUFS100300	10x30cm	40	8	44x30x30 cm
POREPUFS90350	9x35cm	40	8	44x30x30 cm
POREPUFS82350	8.25x35cm	40	8	44x30x30 cm
POREPUFS100350	10x35cm	40	8	44x30x30 cm
POREPUFS83415	8.3x41.5cm	40	8	44x30x30 cm
POREPUFS100160	10x16cm	60	8	44x30x30 cm
POREPUFS100210	10x21cm	40	8	44x30x30 cm
POREPUFS100260	10x26cm	40	8	44x30x30 cm
POREPUFS100400	10x40cm	40	8	44x30x30 cm
POREPUFS120160	12x16cm	60	8	44x30x30 cm
POREPUFS150150	15x15cm	10	24	43X38x25 cm
POREPUFS200200	20x20cm	10	24	53x27x33 cm
POREPUFS2572	2.5x7.2cm	400	8	44x30x30 cm
POREPUFS83100	8.3x10cm	100	8	44x30x30 cm
POREPUFS5072	5x7.2cm	200	8	44x30x30 cm
POREPUFS5075	5x7.5cm	200	8	44x30x30 cm
POREPUFS70100	7x10cm	100	8	44x30x30 cm
POREPUFS6075	6x7.5cm	200	8	44x30x30 cm
POREPUFS85155	8.5x15.5cm	60	8	44x30x30 cm
POREPUFS85150	8.5x15cm	60	8	44x30x30 cm
POREPUFS8595	8.5x9.5cm	100	8	44x30x30 cm
POREPUFS90100	9x10cm	100	8	44x30x30 cm
POREPUFS125125	12.5x12.5cm	10	24	44x30x30 cm

- Amria - Alexandria - Egypt
Tel. : 002 03 4500264 Fax : 002 03 4500263
- Kafr-El-Zayat P.O.Box 38 Egypt
Tel. : 002 040 2545643 Fax : 002 040 2545645
E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : pharmapore family (post operative dressings)						
Issue date : 01/04/2020		Prepared by : Marina Medhat						
Colors : <input type="checkbox"/>		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	13/02/2020	Amen table of size	CRS request					
3	19/02/2020	Amen table of size and add box size	customers needs					
4	06-05-2020	- adding case dimension to table of sizes - modifying sizes - modifying codes - modifying number of units - Inteded Purpose instead of Inteded use	- file maintenance - file update					
5	18-05-2020	adding new packaging configuration according to CRS	Customer request					
6	13 - 12 - 2020	adding new packaging configuration according to CRS	Customer request					
7	02-03-2021	amendments in table of sizes	- CRS Request					
8	25 - 08 - 2021	updating codes	To be aligned with the new coding system					
9	16 - 10 - 2021	Amending size and its related code	Amending typo to maintain the file					
10	17 - 09 - 2022	Adding new code	CRS request					
11	28-05-2023	Replacing "low-adherent" with "non-adherent" in all languages	File maintenance					