

## Thermo Scientific™ Mark-It™ Tissue Marking Dyes

Provide an efficient means of tissue margin identification. User-friendly system uses permanent dyes that maintain their original brilliance throughout the entire routine tissue processing cycle; can be used on all fresh or formalin-fixed tissues.

- Various uses for this system include:
- Marking surgical margins
- Identifying serial segments of large tissues to ensure proper orientation and embedding
- Maintaining proper tissue orientation for Moh's surgery specimens
- Identifying multiple pieces of tissues processed in the same cassette

### Includes

Kit includes one base, one of each of the six colour dyes and applicator sticks.



Description	Qty.	Order Number
Kit	ea.	5000
Black dye, 59ml	ea.	5000BK
Blue dye, 59ml	ea.	5000BL
Green dye, 59ml	ea.	5000G
Orange dye, 59ml	ea.	5000O
Red dye, 59ml	ea.	5000R
Yellow dye, 59ml	ea.	5000Y
Applicator Sticks	864	5001

## Frozen Section Staining Kit

Ideal for laboratories that process a limited number of slides. Works well for staining tissue mounted on coverglass.

- Sturdy cover with piano hinge seals beakers to minimise evaporation
- Rubber feet on rack protect table surface
- Glass beakers have a capacity of 200ml and are 10cm tall
- Slide holder will hold ten 76 x 26mm slides
- Stainless steel support rack measures 21.6 x 48.3 x 22.6mm

### Includes

Kit includes 12 beakers, a stainless steel slide holder and a heavy-gauge stainless steel support rack with six round openings to hold beakers.



Description	Qty.	Order Number
Frozen section staining kit	ea.	145
Glass beakers only	12/pk.	146
Slide holder only	ea.	147
Stainless steel support rack only	ea.	23145