







Trendelenburg Position



Height Adjustment (Low)



Reverse Trendelenburg Position



Height Adjustment (High)

OPTIUM MEDIKAL LTD. STI.

♥ Airport Industrial Zone Hidayet Street No:23 Oguzeli / Gaziantep / TURKEY ↓ +90 850 807 6284 ∞ info@optium.com.tr





- Hydraulic height adjustment by dual sided pedals
- Manual backrest adjustment by gas spring
- Trendelenburg & Reverse Trendelenburg by gas spring
- Epoxy powder panged steel main construction
- Synthetic leather coated mattress
- Steel sheet mattress platform
- Full length lockable safety side rails
- Undercarriage storage for equipment
- IV pole
- Protective bumpers
- Swivel castors, 2 with brake
- Patient safety belt (optional)
- X-ray translucent HPL mattress platform (optional)
- Central locking mechanism (optional)
- Oxygen Holder (optional)

C		•	
	*		

Patient Stret

IV Pole

TECHNICAL DETAILS				
Width:	78 cm			
Length:	207 cm			
Trendelenburg Angle:	0°-12°			
Reverse Trendelenburg Angle:	0°-12°			
Sleep Surface Dimensions:	60x190 cm			
Height Range:	62- 87 cm			
Castor diameter:	15 cm			

