

## Hyper Medozon comfort

### For normobaric and hyperbaric ozone therapy

The **Hyper Medozon comfort** provides the whole spectrum of the ozone therapy.

It is especially developed for the hyperbaric big autohaemotherapy. During this therapy procedure a repeating function is available.

### Highlights of Hyper Medozon comfort

- equipped with an air trap for more safety at the reinfusion
- excellent operating comfort through an intuitive Touchscreen with icons
- multiple setting options for reacting to individual patients requirements
- multilingual operator guidance
- Luer-Lock connections allow a maximum compatibility to the disposables and the accessories, e.g. the Aquazon
- integrated UV-Analyzer for concentration measurement
- suitable for the hyperbaric and normobaric big autohaemotherapy, syringe withdrawal, permanent gassing, vacuum function, water-ozonization (via Aquazon) and intestinal insufflation



\*cart 18140 optional available.

## Technical data:

|                          |                                   |                  |
|--------------------------|-----------------------------------|------------------|
| Name of appliance:       | <b>Hyper Medozon comfort</b>      |                  |
| Power supply:            | 100-240 V / 50-60 Hz              |                  |
| Power input:             | 190 VA                            |                  |
| Ozone concentration:     | 2-80 µg/ml (control via analyser) |                  |
| Flow rate:               | 1,0 l/min                         |                  |
| Dimensions (B x H x T):  | 424 x 185 x 388 mm                |                  |
| Weight:                  | 13,5 kg (with airtrap)            | ■ Ord. no. 16675 |
| Product category acc. to |                                   | ■ Ord. no. 16676 |
| MP-directives 93/42 EWG: | II b                              |                  |

HERRMANN Apparatetechnik GmbH • Im Höning 3 • 63820 Elsenfeld • Germany

☎ +49 (0) 6022-6581-3 • ☎ +49 (0) 6022-6581-59

🌐 www.h-a-b.de • @ info@hab.gmbh

[Home](#)   [Profile](#)
[News](#)   [Production](#)   [Development](#)
[Products](#)
[Accessories](#)
[Ozone](#)   [Colon](#)   [Veterinary](#)   [Cosmetic](#)
[Service](#)
[Shop](#)   [Katalog](#)
[Distribution](#)   [Support](#)   [Dates & Seminars](#)   [Videos](#)   [Media partners](#)
[Contact](#)   [Home](#)
[Profile](#) • [Products](#) • [Accessories](#) • [Service](#) • [Shop](#) • [Katalog](#) • [Contact](#)


## Our accessory for the ozone therapy

Here you will find all the accessories for the normobaric and hyperbaric ozone therapy.


Colon accessory you will find [here](#).

### Medozon i-Set, 300 ml

100% sterile and Ozone resistant  
without latex and DEHP-free  
certified for the normobaric and  
hyperbaric autohaemotherapy,

consisting of:

- connection line with disposable filter and roll clamp
- 300 ml blood chamber with blood filter and injection port
- infusion line with roll clamp
- with hydrophic filter included



### Vacuum bottle 250ml, Glas

Suitable for the Ozone Oxygen Therapy and all common Ozone Therapy devices of company HERRMANN.

### Vacuum bottle 500ml, Glas

Suitable for the Ozone Oxygen Therapy and all common Ozone Therapy devices of company HERRMANN, as well as for use in the veterinary medicine.



### Sodium Citrate 3.13% 10ml Ampoules

Solution to prevent the blood  
coagulation in the Ozone, HOT and  
Own blood Therapy.

### Natrium Citrat 3.13% 100 x 10 ml Ampullen

Lösung zur Verhinderung der  
Blutgerinnung bei der Ozon-, HOT-  
und Eigenbluttherapie



zur Katheterblockung u.a. in der Dialyse.  
VPE 100 x 10 ml Ampullen

**>> direkt bestellen**

### **Natrium Citrat 3.13% 100 ml Stechflasche**


Lösung zur Verhinderung der Blutgerinnung bei der Ozon-, HOT- und Eigenbluttherapie,  
zur Katheterblockung u.a. in der Dialyse.  
VPE 100 ml Stechflaschen

**>> direkt bestellen**

### **PPS Transfusion Set (Sangodrop) with Luer-Lock connection acc. to ISO 594 with steel spike and blood filter 200 my**

Suitable for the Ozone Oxygen Therapy and all common Ozone Therapy devices of company HERRMANN.

Transfusion Set with steel spike serves e. g. for reinfusing the treated native blood from the vacuum bottle in the course of an ozone-oxygen-therapy. Tested for sterility and non pyrogenicity.



### **PPS Heidelberger Extension line 140cm Luer-Lock neutral for the extension of Infusion and Transfusion Sets acc. to DIN58362VLP**

Suitable for the Ozone Oxygen Therapy and all common Ozone Therapy devices of company HERRMANN.  
Connection line for attachment to infusion and transfusion systems. Tested for sterility and non-pyrogenicity.




### **Disposable filter „hyperbaric“ Lock connection**

Suitable for the Ozone Oxygen Therapy and all common Ozone Therapy devices of company HERRMANN.  
Ventilation and Syringe filter also suitable for the hyperbaric Ozone Therapy.



### **Disposable filter "taking by syringe"**

Filter for taking by syringe for the Ozone Oxygen Therapy.



### **Accessory for the colon therapy**

Here you will find our accessories for the colon therapy.

#### **„Standard“ Specula Set for the Colon Hydrotherapy and Colon irrigation**

Suitable for all Colon Hydrotherapy devices existing on the market.  
Set contains: 1 Speculum, corrugated tube and fresh water tube  
Available in two lengths: 1,50 m and 1,20 m  
Packing unit 60 Pcs.

#### **„Oliven“ Specula Set for the Colon Hydrotherapy and Colon irrigation**

Suitable for all Colon Hydrotherapy devices existing on the market.  
Set beinhaltet: 1 Spekulum, corrugated tube and fresh water tube  
Available in two lengths: 1,50 m and 1,20 m  
Packing unit 60 Pcs.

