

# TECHNICAL DATA SHEET Disposable Sterile Standard Surgical Gown

- 1. Product name
  Disposable Sterile Standard Surgical Gown.
- 2. Product code Model 101-502.


3. Product description

Sterile surgical gown made of SMS (Spunbond+Meltblown+Spunbond) three-layer non-woven material according to EN 20811:1981 the resistance to liquid penetration is greater than 50 cm H2O. With long sleeves and elasticated knitted cuffs, neck fastening system with "velcro" tape and back overlap fastening with ties.

Product complies with: ISO 13485, EN 13795, EN 20811:1981, TS EN ISO 811(ISO 811:2018), TS 393 EN ISO 13938-1:2002.

Length of gown according to table. The gown is sterile, disposable and cannot be resterilized. It retains its sterility as long as the package is unopened and undamaged.

- Material
   Fibre composition 100% polypropylene
   Material density 35 gr/m2
   Sterilized Ethylene Oxide.
- 5. Construction features



Tek +90 342 137 24 Fax +90 142 327 31 Sandar 10 001 001 542 54

| Size | Model number | Gown length (cm) | Breast perimeter (cm) |
|------|--------------|------------------|-----------------------|
| S    | 101-502-01   | 117              | 116                   |
| M    | 101-502-02   | 125              | 124                   |
| L    | 101-502-03   | 132              | 132                   |
| XL   | 101-502-04   | 140              | 140                   |
| 2XL  | 101-502-05   | 150              | 148                   |

#### 6. Instructions for use

The equipment is worn over the work coat. It is fastened at the back with straps and at the neck with "velcro" fasteners (hedgehog).

### 7. Packing method

The product is individually packed in a polythene bag with the size and design of the dressing gown

### 8. Storage instructions

The storage is in the original packaging. Store in a dry and clean place away from moisture, from 0°C to +20°C, away from direct sunlight, water droplets and solvent vapors.

9. Warranty term 3 years from date of production.

## FACTORY CENTER

2. Organized industrial Zone No:1 Başpınar/Gaziantep Tel.: +90 (342) 337 24 16; +90 (342) 337 21 11

E-mail: info@3teks.com.tr

