

Adjustable Pipette



Model	Volume range(μ l)	Increment(μ l)	Measure volume (μ l)	Inacc.%	Impr.%	Tip
APT1-10	0.5-10 μ l	0.1 μ l	10 μ l	± 3.0	≤ 12.01	10 μ l
			5 μ l	± 6.0	≤ 13.01	
			1 μ l	± 10.0	≤ 13.01	
APT1-20	2-20 μ l	0.5 μ l	20 μ l	± 3.0	≤ 11.51	200 μ l
			10 μ l	± 3.0	≤ 12.01	
			2 μ l	± 10.0	≤ 13.01	
APT1-50	5-50 μ l	0.5 μ l	50 μ l	± 2.4	≤ 11.51	200 μ l
			25 μ l	± 2.8	≤ 12.01	
			5 μ l	± 6.0	≤ 13.01	
APT1-100	10-100 μ l	1 μ l	100 μ l	± 1.5	≤ 11.01	200 μ l
			50 μ l	± 2.2	≤ 11.51	
			10 μ l	± 3.0	≤ 12.01	
APT1-200	20-200 μ l	1 μ l	200 μ l	± 1.5	≤ 11.01	200 μ l
			100 μ l	± 1.5	≤ 11.01	
			20 μ l	± 3.0	≤ 11.51	
APT1-1000	100-1000 μ l	5 μ l	1000 μ l	± 1.2	≤ 10.51	1000 μ l
			500 μ l	± 1.5	≤ 10.51	
			100 μ l	± 1.5	≤ 11.01	
APT1-5000	1000-5000 μ l	50 μ l	5000 μ l	± 1.0	≤ 10.51	5000 μ l
			2500 μ l	± 1.0	≤ 10.51	
			1000 μ	± 1.2	≤ 10.51	

Multichannel-volume pipette 8 channel/12 channel



8 channel / Model	Volume range (μl)	increment (μl)	Tip
APT8-10	0.5-10 μl	0.1 μl	10 μl
APT8-50	5-50 μl	0.5 μl	200 μl
APT8-300	50-300 μl	5 μl	300 μl /350 μl
12 channel / Model			
APT12-10	0.5-10 μl	0.1 μl	10 μl
APT12-50	5-50 μl	0.5 μl	200 μl
APT12-300	50-300 μl	5 μl	300 μl /350 μl