


Product-Information-Sheet

Polypectomy, EMR and ESD


Cold and Hot
Snares



- For polyp ligation prior to polypectomy
- Prevents postpolypectomy bleeding
- Preloaded loop, ready to use with handle
- Nylon-loop ligation keeps tissue damage to a minimum
- Opening diameter of 30 mm enables ligation of large polyps




Article number	Article Name	Quantity value	Min. working channel ø value	Working length value	Loop ø	Area Compatibility
N5382130	HX-400U-30	5	● 2.8	2300 mm	30 mm	



- L-shaped hook for incision and dissection in longitudinal and lateral directions
- To minimise invasion, the use of a distal attachment is recommended (D-201 series) to help maintain the field of view
- Rotatable



Article number	Article Name	Quantity value	Min. working channel ø value	Working length value	Knife length value	Features	Area Compatibility
N3046130	KD-620UR	1	● 2.8	2300 mm	4.5 mm	Hook length: 1.3 mm; rotatable	



- Triangle-tipped knife for marking, incision and dissection
- Conductive triangle tip for cutting the mucosa
- No need to rotate the knife; cutting can be performed in any direction
- Tip design suitable for controlling moderate bleeding
- To ensure a clear field of view, a distal attachment (D-201 series) may be used



Article number	Article Name	Quantity value	Min. working channel \varnothing value	Working length value	Knife length value	Features	Area Compatibility
N2119630	KD-640L	1	● 2.8	1650 mm	4.5 mm	Triangle-tip plate length: 0.4 mm	