




GE- Marquette	E9004GE; 2006644- 001	ADP-10010A	i9500 with Tram module 1/451N/851/851N CardioLab®/Mac- Lab® with Tram Module 451/451N Corometrics® 120/2120 series, Eagle 1000/3000/4000 Responder 3000 with Nellcor option; use with Nellcor spo2 sensor, such as RSP- 10041A		GE11P-DB9F
GE- Marquette Oximax	2021406- 001	ADP-10010B	Dash 3000/4000/5000 with Nellcor OxiMax option Solar 8000/8000M, i9500 with Tram module 451N5/851N5		GE11P-DB9F
GE Masimo	2027263- 002	ADP-10010C	Dash 3000/4000/5000 with Masimo option Solar 8000/8000M, i9500 with Tram Module 451M/851M CardioLab®/Mac-Lab® with Tram Module 451M Corometrics® 259		GE11P- DB9F(Masimo)
	2017002- 001 (3.6m) 2017002- 003 (2.4m); 2002592- 001(3.6m); 2002592- 002 (2.4m)	ADP-10010D			Interconnect cable, Masimo LNOP, SpO2
GE- Ohmed a	OXY-ES3	ADP-10010E	S/5 FM port, Anesthesia monitor; All CARESCAPE™ modular monitors and S/5, AS/3 and CS/3 modular monitors with E-PRESTN and E-PSM modules. CARESCAPE V100 with TruSignal® option.		OxyTip+® Interconnect cable with GE to UN connector
	TS-G3	ADP-10010F	CARESCAPE™ modular monitors and S/5, AS/3 and CS/3 modular monitor with E-(P)RESTN or E- PSM(P) modules. CARESCAPE V100 with TruSignal® option, Responder 2000.		TruSignal Interconnect cable with GE to TruSignal connector
Goldwa y		ADP-10011A	UT4000A,UT4000C,UT4000E,UT4 000F,UT4000Apro,UT4000FPro,UT 6000A		Lemo5P-DB9F
		ADP-10011B	Goldway UT4000A,Pace Tech		Redel5P-DB9F
		ADP-10011C	Goldway UT4000B		Redel7P-DB9F
		ADP-10011D	Goldway OXIMAX type, 7 pins		Redel7P-DB9F