

**PAGINA 0**

**CREATA PENTRU A  
PUTEA SEMNA  
ELECTRONIC  
DOCUMENTUL**



HEALTHCARE PROFESSIONALS  
**Cragg-McNamara Valved Infusion Catheters**  
 Infusion Therapy

RELATED PRODUCTS

- Infusion Therapy Products
- Cardiovascular Products

REQUEST MORE INFORMATION

[CONTACT US](#)

PERIPHERAL CATALOGUE

Choose from a full portfolio of peripheral products for endovascular interventions – products that help you reach your PAD patient treatment goals.

[DOWNLOAD CATALOGUE](#)

LIFELINE CARDIOVASCULAR TECHNICAL SUPPORT

**Worldwide Technical Support - 24/7**  
 +1 763-514-4000

[rcstechsupport@medtronic.com](mailto:rcstechsupport@medtronic.com)

[Indications, Safety, and Warnings](#)

OVERVIEW

Infuse with ease, with or without a guidewire in place, thanks to the Cragg-McNamara™ valved infusion catheter. Offered in 4 and 5 F sizes, the valved tip single lumen catheter gives you the option to deliver anti-thrombus agents without a tip-occluding guidewire, for streamlined patient care.

INDICATIONS

The Cragg-McNamara™ Micro Therapeutics Infusion Catheter is intended to be used for the controlled selective infusion of physician-specified pharmacologic agents or radiopaque contrast media into the general vasculature. All pharmacologic agents utilised with the Micro Therapeutics Infusion Catheter should be fully prepared and used according to the instructions for use of the specified pharmacologic agent. The Micro Therapeutics Infusion Catheter is not intended for coronary, pediatric or neonatal use.

area available for infusion.

Infusion Lumen Comparison

Compare the infusion lumen areas of the Cragg-McNamara™ (below) and other 5 F infusion catheters.



STREAMLINE PATIENT CARE

- Infuse without a guidewire in place, which eliminates the risk of guidewire movement, for simplified patient care.
- Treat without a tip-occluding guidewire and sidearm adapter, which may contribute to cost savings.

MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the product labeling supplied with each device.

MODEL SPECIFICATIONS

CRAGG-McNAMARA™ VALVED INFUSION CATHETER

| Reference Number | Diameter (F) | Usable Length (cm) | Infusion Length (cm) |
|------------------|--------------|--------------------|----------------------|
| 41032-01         | 4            | 40                 | 10                   |
| 41033-01         | 4            | 40                 | 20                   |
| 41034-01         | 4            | 65                 | 5                    |
| 41035-01         | 4            | 65                 | 10                   |
| 41036-01         | 4            | 65                 | 20                   |
| 41037-01         | 4            | 100                | 5                    |
| 41038-01         | 4            | 100                | 10                   |
| 41039-01         | 4            | 100                | 20                   |
| 41040-01         | 4            | 135                | 5                    |
| 41041-01         | 4            | 135                | 10                   |
| 41042-01         | 4            | 135                | 20                   |
| 41043-01         | 5            | 40                 | 5                    |
| 41044-01         | 5            | 40                 | 10                   |
| 41045-01         | 5            | 40                 | 20                   |
| 41046-01         | 5            | 65                 | 5                    |
| 41047-01         | 5            | 65                 | 10                   |
| 41048-01         | 5            | 65                 | 20                   |
| 41049-01         | 5            | 100                | 5                    |
| 41050-01         | 5            | 100                | 10                   |
| 41051-01         | 5            | 100                | 20                   |
| 41052-01         | 5            | 100                | 30                   |
| 41053-01         | 5            | 100                | 40                   |
| 41054-01         | 5            | 100                | 50                   |
| 41055-01         | 5            | 135                | 5                    |
| 41056-01         | 5            | 135                | 10                   |
| 41057-01         | 5            | 135                | 20                   |
| 41058-01         | 5            | 135                | 30                   |
| 41059-01         | 5            | 135                | 40                   |
| 41060-01         | 5            | 135                | 50                   |

See product catalog for complete, detailed product information.

- HEALTHCARE PROFESSIONALS
- Overview
  - Products
  - Therapies & Procedures
  - Integrated Health Solutions

- PATIENTS
- Overview

- OUR COMPANY
- Who We Are
  - Key Facts
  - Leadership
  - History
  - Mission
  - Locations
  - Governance
  - Careers

- OUR IMPACT
- Our Impact
  - Healthcare access
  - Inclusion, Diversity & Equity
  - Protecting our planet
  - ESG Report



HEALTHCARE PROFESSIONALS

## Cragg-McNamara Valved Infusion Catheters

Infusion Therapy

RELATED PRODUCTS

- Infusion Therapy Products
- Cardiovascular Products

REQUEST MORE INFORMATION

CONTACT US

PERIPHERAL CATALOGUE

Choose from a full portfolio of peripheral products for endovascular interventions – products that help you reach your PAD patient treatment goals.

DOWNLOAD CATALOGUE

LIFELINE CARDIOVASCULAR TECHNICAL SUPPORT

Worldwide Technical Support - 24/7  
+1 763-514-4000

rs-technicalsupport@medtronic.com

Indications, Safety, and Warnings

OVERVIEW

Infuse with ease, with or without a guidewire in place, thanks to the Cragg-McNamara™ valved infusion catheter. Offered in 4 and 5 F sizes, the valved tip single lumen catheter gives you the option to deliver anti-thrombus agents without a tip-occluding guidewire, for streamlined patient care.

INDICATIONS

The Cragg-McNamara™ Micro Therapeutics Infusion Catheter is intended to be used for the controlled selective infusion of physician-specified pharmacologic agents or radiopaque contrast media into the general vasculature. All pharmacologic agents utilised with the Micro Therapeutics Infusion Catheter should be fully prepared and used according to the instructions for use of the specified pharmacologic agent. The Micro Therapeutics Infusion Catheter is not intended for coronary, pediatric or neonatal use.

area available for infusion.

Infusion Lumen Comparison

Compare the infusion lumen areas of the Cragg-McNamara™ (below) and other 5 F infusion catheters.



STREAMLINE PATIENT CARE

- Infuse without a guidewire in place, which eliminates the risk of guidewire movement, for simplified patient care.
- Treat without a tip-occluding guidewire and sidearm adapter, which may contribute to cost savings.

MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the product labeling supplied with each device.

MODEL SPECIFICATIONS

CRAGG-McNAMARA™ VALVED INFUSION CATHETER

| Diameter (F) | Usable Length (cm) | Infusion Length (cm) | Maximum Guidewire (in) |
|--------------|--------------------|----------------------|------------------------|
| 4            | 40                 | 10                   | 0.035                  |
| 4            | 40                 | 20                   | 0.035                  |
| 4            | 65                 | 5                    | 0.035                  |
| 4            | 65                 | 10                   | 0.035                  |
| 4            | 65                 | 20                   | 0.035                  |
| 4            | 100                | 5                    | 0.035                  |
| 4            | 100                | 10                   | 0.035                  |
| 4            | 100                | 20                   | 0.035                  |
| 4            | 135                | 5                    | 0.035                  |
| 4            | 135                | 10                   | 0.035                  |
| 4            | 135                | 20                   | 0.035                  |
| 5            | 40                 | 5                    | 0.038                  |
| 5            | 40                 | 10                   | 0.038                  |
| 5            | 40                 | 20                   | 0.038                  |
| 5            | 65                 | 5                    | 0.038                  |
| 5            | 65                 | 10                   | 0.038                  |
| 5            | 65                 | 20                   | 0.038                  |
| 5            | 100                | 5                    | 0.038                  |
| 5            | 100                | 10                   | 0.038                  |
| 5            | 100                | 20                   | 0.038                  |
| 5            | 100                | 30                   | 0.038                  |
| 5            | 100                | 40                   | 0.038                  |
| 5            | 100                | 50                   | 0.038                  |
| 5            | 135                | 5                    | 0.038                  |
| 5            | 135                | 10                   | 0.038                  |
| 5            | 135                | 20                   | 0.038                  |
| 5            | 135                | 30                   | 0.038                  |
| 5            | 135                | 40                   | 0.038                  |
| 5            | 135                | 50                   | 0.038                  |

See product catalog for complete, detailed product information.

HEALTHCARE PROFESSIONALS

- Overview
- Products
- Therapies & Procedures
- Integrated Health Solutions

PATIENTS

- Overview

OUR COMPANY

- Who We Are
- Key Facts
- Leadership
- History
- Mission
- Locations
- Governance
- Careers

OUR IMPACT

- Our Impact
- Healthcare access
- Inclusion, Diversity & Equity
- Protecting our planet
- ESG Report