

TECHNICAL SHEET

Page 2/2

DISPOSABLE SYRINGE

Scope of application

The suction liquid medicine dissolving after intradermal, subcutaneous, for human body, muscle, intravenous injection.

REGULATORY INFORMATION


Name	Sterile Syringes for Single Use
size	1ml,2ml,3ml,5ml,10ml,20ml,30ml,50ml,60ml
Luer type	Luer lock, luer slip
Material	Medical grade PP
Needle	With or without needle
Piston	Latex or latex free
Sterile	Sterilized by EO gas, non-toxic, non-pyrogenic
Packing	PE package, carton packing
Syringe type	3-part
Brand Name	Changzhou Jinlong Medical Plastic Appliance Co.,Ltd
Certificated	CE ISO13485



TECHNICAL SHEET

Page 2/2

DEVICES DESCRIPTION:



MAIN COMPONENTS:

1. BARREL:PP
2. PLUNGER:PP
3. PISTON:NATURE RUBBER
4. PROTECTIVE CAP OF NEEDLE:PP
5. NEEDLE HUB:PP
6. NEEDLE TUBE:STAINLESS STEEL

TECHNICAL FEATURES:

Device description with variant.

Model: 1mL、2mL、3mL、5mL、10mL、20mL、30mL、50mL、60mL.

needle: 0.3、0.33、0.36、0.4、0.45、0.5、0.55、0.6、0.7、0.8、0.9、1.1、1.2.

Component: Plunger, Barrel, Piston, needle.

life time of the product: After being sterilized by EO, the effective period is five years.

Disposable syringe with 3 parts, luer lock or luer slip

CE&ISO certificate

transparent barrel is easy to master the volumn

transparent barrel is also easy to wipe bubble out

the plunger is moved smoothly inside the barrel

latex or latex free

the material is medical grade PP

sterilized by EO gas

DISPOSABLE SYRINGE

PRECAUTIONS FOR USE:

use for once, don't reuse

if the PE bag is broken, don't use it

throw the used syringes properly

store in clean and dry space

PRODUCTION :1 MILLION ONE DAY

STERILIZATION:STERILIZED BY EO GAS, NON-TOXIC, NON-PYROGENIC

PRESERVATION:

Shelf life in months:60 Months

Storage:store in clean and dry space

UNIT OF USE:PE PACKAGE:1 PIECE