



Home » Products » Interventional Product » Balloon Inflation Device BID1 Type



Balloon Inflation Device BID1 Type

· It is used in PTCA operation to pressurize the balloon dilatation catheter, inflate and contract the balloon accurately, so as to achieve the purpose of dilating blood vessels or indwelling stents in blood vessels.

INQUIRY

FEATURES

- The maximum output pressure is 30atm, and the maximum capacity is 30ml, which is convenient for medical staff to use
- A safety lockout controls pressure and releases it quickly
- The pressure gauge reads accurately and is easy to use



- Locking design enables one handed operation in deflation

PRODUCT QUICK LOOK

Product at a glance -



PARAMETER


Product No.	Volume (ml)	Contents
42.04.30000	20	Balloon inflation device, Stopcock("RIGHT ON")
42.04.30001	20	Balloon inflation device, Stopcock("RIGHT OFF")
42.04.30002	30	Balloon inflation device, Stopcock("RIGHT ON")
42.04.30003	30	Balloon inflation device, Stopcock("RIGHT OFF")
42.04.40000	20	Balloon inflation device, Stopcock("RIGHT ON")+ Hemostasis valve set (Twist)
42.04.40001	20	Balloon inflation device, Stopcock("RIGHT OFF")+ Hemostasis valve set (Twist)
42.04.40002	30	Balloon inflation device, Stopcock("RIGHT ON")+ Hemostasis valve set (Twist)
42.04.40003	30	Balloon inflation device, Stopcock("RIGHT OFF")+ Hemostasis valve set (Twist)



2.04.40014	30	Balloon inflation device, Stopcock("RIGHT ON")+ Hemostasis valve set (Push Pull)
42.04.40015	30	Balloon inflation device, Stopcock("RIGHT OFF")+ Hemostasis valve set (Push Pull)
42.04.40004	20	Balloon inflation device, Stopcock("RIGHT ON")+ Hemostasis valve set (Push Pull with 25cm extension line & stopcock)
42.04.40005	20	Balloon inflation device, Stopcock("RIGHT OFF")+ Hemostasis valve set (Push Pull with 25cm extension line & stopcock)
42.04.40006	30	Balloon inflation device, Stopcock("RIGHT ON")+ Hemostasis valve set (Push Pull with 25cm extension line & stopcock)
42.04.40007	30	Balloon inflation device, Stopcock("RIGHT OFF")+ Hemostasis valve set (Push Pull with 25cm extension line & stopcock)
42.04.40026	20	Balloon inflation device, Stopcock("RIGHT OFF")+ Hemostasis valve set (Twist type with 25cm extension line & stopcock)
42.04.40027	20	Balloon inflation device, Stopcock("RIGHT ON")+ Hemostasis valve set (Twist type with 25cm extension line & stopcock)
42.04.40017	30	Balloon inflation device, Stopcock("RIGHT ON")+ Hemostasis valve set (Twist type with 25cm extension line & stopcock)
42.04.40016	30	Balloon inflation device, Stopcock("RIGHT OFF")+ Hemostasis valve set (Twist type with 25cm extension line & stopcock)



What else do you want to know? Contact us right now!

*Your Name	Country
*Your Email	Your Phone
*Message	
*Verify Code	

*We respect your privacy. When you submit your contact information, we agree to only contact you in accordance with our [Privacy Policy](#).

Submit

We'd Love to Hear From You



Contact Us

Contact

Email: info@scwmed.com

Mobile Phone: +8613691980679

Telephone: 86-755-89926266 ext 4

Fax: 86-755-89668199 ext 827

NO.4 Baolong 6th Road, Baolong industrial Town,
Longgang District, Shenzhen, P.R. China



Copyright © SCW Medicath Ltd All rights reserved. [Privacy](#) [Site Map](#) [Terms of use](#)