

Innovative range of masks N.I.V. Bluestar Plus



The air at service of breathing

Lightness

Seal

Comfort

KOO MEDICAL EQUIPMENT www.koomedical.com

The lightness of the air
The seal of the air cushion
The comfort of the anatomic cushion

Facial Mask - Intensive

bi level CPAP Non-vented

Its design makes it widely versatile, the presence in the Kit of 3 connectors allows it to operate with different configurations of circuits. The connector with hole for exhalation allows to transform the mask from **Non-vented** to **vented**. A port on the connector allows the O2 enrichment or monitoring.



BLUESTAR PLUS Intensive

- Replaceable anatomic air cushion
- Wide range of fixing points
- Headgear mask fixer in highly breathable material and with wide bearing surface for better forces distribution
- Operative openings for O2 enrichment and for monitoring
- Port for the passage of the gastric tube, or the introduction of the bronchoscope, or for simple needs of the patient
- Front support fully adjustable and removable



The lightness of the air
The seal of the air cushion
The comfort of the anatomic cushion

Facial Mask - Secure

bi level CPAP Non-vented

The presence of a safety valve Anti-Asphyxia, permits a reasonable airflow even in case of the flow lack coming from the fan.

There is a rotating connector a 360° blue color with a diameter of 22 mm F.



BLUESTAR PLUS Secure with Security valve A.A.V.

- Replaceable anatomic air cushion
- Wide range of fixing points
- Headgear mask fixer in highly breathable material and with wide bearing surface for better forces distribution
- Operative openings for O2 enrichment and monitoring.
- Port for the passage of the gastric tube, or the introduction of the bronchoscope, or for simple needs of the patient
- Safety valve with immediate opening in the event of flow lack
- Front support fully adjustable and removable



The lightness of the air
The seal of the air cushion
The comfort of the anatomic cushion

Facial Mask - Home CARE

Bi-level CPAP Vented

The presence of exhalation holes, the safety valve Anti Asphyxia and the replaceable air cushion, make this mask the most versatile for use in the ventilation of the patient at home-care.





- Replaceable anatomic air cushion
- Wide range of fixing points
- Headgear mask fixer in highly breathable material and with wide bearing surface for better forces distribution
- Operative openings for O2 enrichment and for monitoring.
- Port for the passage of the gastric tube, or the introduction of the bronchoscope, or for simple needs of the patient
- Safety valve with immediate opening in the event of flow lack
- Exhalation holes for a correct evacuation of the exhaled breath
- Front support fully adjustable and removable



Our range





KM -	310 B	Extra Large	- 11::-
KM -	311 B	Large	Each kit includes: Mask with air cushion, transparent connector 22 mm M,
KM -	312 B	Medium	blue connector 22 mm F, exhalation valve, headgear
KM -	313 B	Small	mask fixer, bulb pump.
KM -	314 B	Pediatric	



Bluestar Plus Secure AVV

KM – 310 V	Extra Large	
KM – 311 V	Large	Each kit includes: Mask with air cushion, blue connector 22 mm F with
KM – 312 V	Medium	
KM – 313 V	Small	security valve, headgear mask fixer, bulb pump.
KM – 314 V	Pediatric	



Bluestar Plus Home CARE

KM – 310 E	Extra Large	
KM – 311 E	Large	Each kit includes:
KM – 312 E	Medium	Mask with air cushion, transparent connector 22 mm M with security valve and exhalation holes, headgear mask
KM – 313 E	Small	fixer, bulb pump.
KM – 314 E	Pediatric	



Accessories

KM – 3101	Extra Large	
KM – 3111	Large	
KM – 3121	Medium	Cushion replacement
KM – 3131	Small	
KM - 3141	Pediatric	
KM - 3114	Transparent connector 22 mm M	
KM – 3113	Blue connector 22 mm F	
KM - 3104	Forehead Space	
KM – P009	Headgear Large	
KM – Po10	Headgear Small	
KM - 903	Exhalation valve	



The use of the air cushion allows a reduction of weight.

A mask lighter requires a lower pulling force
from the headgear, to keep it stable

The seal of the air cushion

The air cushion, anatomically designed, ensures an excellent control of the seal, avoiding losses peripherals unpleasant

The comfort of the anatomical cushion

The anatomical air cushion allows a distribution of the pressure sealing eliminating dangerous pressure's drops.

The complete adjustment of the mask and the use of a headgear made of breathable material, makes the mask highly accepted by most patients, preventing the interruption of ventilation, during use of the mask.

The best handy and low-cost solution for non-invasive ventilation with a device that allows the replacement of the cushion and the shell



Responsible for **EMEA (Europe Middle East Africa) Market Koo Europe s.r.l.**

viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market Koo Asia Ltd**

Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Koo Medical Equipment (Shanghai) Co.,Ltd. www.koomedical.com

1 Songjiang Industrial Park - Songjiang, Shanghai, 201613 P.R.China