

Abbreviation

AEC	: Automatic Exposure Control
APR	: Anatomical Program
HU	: Heat Unit
IBS	: Image Brightness Stabilizer
I.I.	: Image Intensifier
kV	: kilovolt; unit of tube voltage
LED	: Light Emitting Diode
mA	: milliamperes; unit of tube current
mAs	: milliamperes second; unit of tube current time product
min	: minute; unit of fluoroscopy time
ms (msec)	: millisecond; unit of radiography time
sec	: second; unit of radiography time and fluoroscopy time
RIS	: Radiology Information System
CR	: Computed Radiography

Options

The "Options" described in this manual may be provided as standard, or may be unavailable, depending on the sales format.

System Activation

This equipment needs activation during installation.

The activation work needs authorization by Shimadzu.

*This is available for X-ray High Voltage Generator shipped after January, 2020.

Related Documents

M516-E009 High Speed Rotation Starter SA-60 Operation Manual