

CHIRASORB Plus

Antibacterial

CHIRASORB Plus is a mid-term absorbable, antibacterial coated, synthetic, multifilament suture made of 90% glycolide and 10% L-lactide. The braided suture is covered with antibacterial coating, which consists of copolymer glycolide and lactide, calcium stearate and chlorhexidine diacetate.

Chlorhexidine is a cationic antiseptic with a wide range of antimicrobial effects. It is effective on both Gram-positive and Gramnegative bacteria; the mechanism of action is cytoplasmic membrane disruption.

Sterile Chirasorb Plus is available in violet or natural. The product elicits a minimal acute inflammatory reaction in tissues.

Progressive loss of tensile strength and absorption of Chirasorb Plus occurs by hydrolysis, when the polymer is degraded to glycolic and lactic acids, which are subsequently absorbed and metabolised by the body.

Indication: General surgery, Gynaecology, Skin Closure (e.g. intracutaneous, subcutaneous),

Muscular layer and Ligation.

Contraindication: Chirasorb Plus should not be used in tissues subject to expansion, stretching or disten-

sion, or which might require long-term mechanical support. Chirasorb Plus should not be

used in patients with known allergic reactions to Chlorhexidine.

Sizes: 6/0 to 4 USP

Structure: Multifilament

Absorption method: Hydrolysis

Strength Retention

Profile: MID-TERM HEALING (35 days)

1 Day 100%

14 Days 75%

21 Days 50%

Mass Absorption

profile: 56 - 70 days

UK03

85



Absorbable - Chirasorb Plus



CHIRASORB Plus

eedle	EP	USP	Length	Colour	Cat. No	Other information
RT: ½ circle, roun	d body needle	- tapercut				
HRT 17	1.5	4/0	1 x 45 cm	Violet	CV 0566	II THE
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	CV 0561 CV 0562	
HRT 25		3/0	1 x 75 cm	Violet	CV 0581	
, ⊗ п	2 3	2/0	1 x 75 cm	Violet	CV 0582	
	3.5	0	1 x 75 cm	Violet	_CV 0583	
HRT 30	3.5	0	1 x 75 cm	Violet	CV 0586	
	4	1	1 x 75 cm	Violet	CV 0587	
	4	1	1 x 90 cm	Violet	CV 0589	
HRT 31		Book.				
(0)	3	2/0	1 x 75 cm	Violet	CV 0606	
			10			
LIDT 27	3	2/0	1 x 75 cm	Violet	CV 0592	
HRT 37	3.5	0	1 x 75 cm	Violet	CV 0593	
(®)) (; -1 -A	CV 0507	
	3	2/0	1 x 90 cm 1 x 90 cm	Violet Violet	CV 0597 CV 0598	
	3.5 4	0 1	1 x 90 cm	Violet	CV 0599	
HRT 40	3.5	0	1 x 75 cm	Violet	CV 0653	
	4	1 -	1 x 90 cm	Violet	CV 0651	D III
HRT 45	\		NO CLINE			
⊗ <i>∏</i>	\ 2	2/0	1 x 75 cm	Violet	CV 0663 CV 0665	
	3 5	2/0	1 x 75 cm 1 x 75 cm	Violet	CV 0665	
HRT 48	3	2/0	1 x 75 cm	Violet	CV 1671	
$\ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ $	5	2	1 x 90 cm	Violet	CV 1668	1

UK03 **91**



