

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

## Bacterial/Viral Filter, with HME

Bacterial/Viral Filter with HME and CO<sub>2</sub>-port. The purpose of the filter is to protect the patient and to humidify the inhaled air. The filter is placed between the endotracheal tube/tracheal cannula and the Y-piece of the breathing system.

**Shape (Ventilation Filter):** Round

**ISO style (patient end):** Male/Female

**Inner diameter (patient end):** 20/15 mm

**Sterile:** No

**Country availability:** Please note: not all products are available in all countries. Feel free to contact our international sales team about any questions relating to our sales and services.

Jutta Dorenberg

Director International Sales

Tel.: +49 2236 3913 - 405

eMail: [jutta.dorenberg@dahlhausen.de](mailto:jutta.dorenberg@dahlhausen.de)

shorten text

<b>Article no.</b>	<b>Application area</b>	<b>Inner diameter (machine end)</b>	<b>ISO style (machine end)</b>	<b>Max. tidal volume (Breathing Filter)</b>	<b>minimal quantity</b>
4300906120	Child	13	Male	900	50
4300906020	Adult	22/13	Male/Female	1.500	25
4300905708	Adult	22/13	Male/Female	1.500	50