

Chemila
Chemila



Chemila, spol. s r.o., Za Dráhou 4386/3, Hodonín 69501, Phone +420518340919, chemila@chemila.cz
Chemical and Microbiological Laboratory, Testing Laboratory No. 1273 certified by Czech Accreditation Institute according to ČSN EN ISO/IEC 17025:2005.

CONFIDENTIAL

Copy No.: 1
Issue No.: 1

Test report No. S186/2019

DETERMINATION OF FUNGICIDAL (EN 13624:2013, EN 14562:2006)
ACTIVITY OF THE PRODUCT **GLOBACID SF 2%**

Sample ID: S186/2019

Sample name: **GLOBACID SF 2%**

Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404

Producer: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404

Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404

Page: 1
From pages: 9

Incoming date:
29.5.2019

Delivery date:
22.7.2019

Hodonín, 22.7.2019



.....
Ing. Jana Šlitrová, Head of Laboratory

The report may be reproduced only as a whole, in parts only upon written permission of the laboratory. The test results relate only to the samples stated in the Test Report. The Lab does not take any guarantee for the identity of samples not taken by the lab personnel.

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019
Rep No: 93
Sample name: GLOBACID SF 2%
Sampled: by client
Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Batch No: 16-322A

Sampling date: 24.5.2019
Sample delivered: 29.5.2019
Testing date: 28.6. – 17.7.2019
Delivered amount: 500 ml
Page: 2

Subject of testing:

Determination of fungicidal activity of the product.

CONFIDENTIAL

Identification of the sample:

Name of the product: GLOBACID SF 2%
Batch number: 16-322A
Date of manufacture: 09.2018
Expiry date: 09.2021
Manufacturer: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Incoming date: 29.5.2019
Storage conditions: stated by the manufacturer
Active ingredients:
CAS 2372-82-9 N-(3aminopropyl)-N-dodecyl propane-1,3-diamine <1%
CAS 94667-33-1 quaternary ammonium compounds <1%

Experimental conditions:

Testing of disinfecting efficiency of chemical disinfecting and antiseptic agents by suspension method

SOP-M-19-00 (EN 13624:2013)

Period of analysis: 28.6. – 1.7.2019
Test temperature: 20 °C ± 1 °C
Test method: dilution neutralization method
Neutralization medium: Dey-Engley Neutralizing Broth M 1062
Appearance of the product: colourless liquid
Test concentration: 100% (concentrated)*
Contact time: 5 min, 15 min, 30 min
Interfering substances: 3 g/l BSA and 3 ml/l sheep erythrocytes (dirty conditions)
Test organisms: *Candida albicans* ATCC 10231
Aspergillus brasiliensis (niger) ATCC 16404
Incubation conditions: 30 °C ± 1 °C, 48 hours and additional period of 24 or 48 hours

Test procedure:

1. Preparation of test suspension
2. Preparation of product test solutions
3. Quantitative suspension test
4. Incubation and calculation
5. Expression and interpretation of results

Note:

Presence of a high concentration (at least 75%) of *Aspergillus brasiliensis* spiny spores in the test suspension – yes.

Fungicidal activity – the capability of a product to produce a reduction in the number of viable fungi belonging to reference strains under defined conditions by at least a 4 lg reduction (10^4).

Yeasticidal activity – the capability of a product to produce a reduction in the number of viable yeast cells of relevant test organisms under defined conditions by at least a 4 lg reduction (10^4).

$R = N_0 / N_a$ = the reduction in viability, or $\lg R = \lg N_0 - \lg N_a$

* Product can only be tested at a concentration of 97% (RTU product) or less, as some dilution is always produced by adding the test organisms and interfering substance.

The standard:

EN 13624:2013 Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of fungicidal or yeasticidal activity in the medical area - Test method and requirements (phase 2, step 1) September 2013

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019
 Rep No: 93
 Sample name: GLOBACID SF 2%
 Sampled: by client
 Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
 Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
 Batch No: 16-322A

Sampling date: 24.5.2019
 Sample delivered: 29.5.2019
 Testing date: 28.6. – 17.7.2019
 Delivered amount: 500 ml

Page: 3

CONFIDENTIAL

The Number of CFU in the tested product: 0 CFU/ml

1. Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Candida albicans* ATCC 10231
 Tab No. 1.1 Verification of methodology, dirty conditions

Validation of suspension (N _{vo})			Validation of selected experimental conditions (A)			Neutralizer toxicity control (B)			Method validation (C) Product conc. 100%*		
V _{e1}	36	Φ _{Nvo} = 33	V _{e1}	28	Φ _A = 30.5	V _{e1}	25	Φ _B = 28	V _{e1}	24	Φ _C = 28
V _{e2}	30		V _{e2}	33		V _{e2}	31		V _{e2}	32	
30 ≤ Φ _{Nvo} ≤ 160			Φ _A > 0.5 Φ _{Nvo}			Φ _B ≥ 0.5 Φ _{Nvb}			Φ _C ≥ 0.5 Φ _{Nvc}		
x	yes	no	x	yes	no	x	yes	no	x	yes	no
Validation of suspension (N _{vb})			V _{e1}	29	V _{e2}	35	Φ _{Nvb}	32	30 ≤ Φ _{Nvb} (N _{vb} /1000) ≤ 160		
									x	yes	no

Tab No. 1.2 Test suspension

Test suspension N	N	V _{e1}	V _{e2}	Test suspension N ₀ (time = 0)
Φ = 249 × 10 ⁶ = lg 8.40	10 ⁻⁶	235	255	lg N ₀ = lg N/100 = lg 6.40
8.17 ≤ lg N ≤ 8.70	10 ⁻⁷	24	34	6.17 ≤ lg N ₀ ≤ 6.70
				x
				yes
				no

Tab No. 1.3 Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Candida albicans* ATCC 10231

Test concentration (%)*/contact time (min)/conditions	Dilution after test procedure	V _{e1}	V _{e2}	lg N _a = lg (Φ _a × 10)	lg R (lg N ₀ = lg 6.40)
100 / 5 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.25
100 / 15 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.25
100 / 30 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.25

2. Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Aspergillus brasiliensis (niger)* ATCC 16404

Tab No. 2.1 Verification of methodology, dirty conditions

Validation of suspension (N _{vo})			Validation of selected experimental conditions (A)			Neutralizer toxicity control (B)			Method validation (C) Product conc. 100%*		
V _{e1}	43	Φ _{Nvo} = 49	V _{e1}	37	Φ _A = 44	V _{e1}	54	Φ _B = 41.5	V _{e1}	46	Φ _C = 38.5
V _{e2}	55		V _{e2}	51		V _{e2}	29		V _{e2}	31	
30 ≤ Φ _{Nvo} ≤ 160			Φ _A > 0.5 Φ _{Nvo}			Φ _B ≥ 0.5 Φ _{Nvb}			Φ _C ≥ 0.5 Φ _{Nvc}		
x	yes	no	x	yes	no	x	yes	no	x	yes	no
Validation of suspension (N _{vb})			V _{e1}	50	V _{e2}	49	Φ _{Nvb}	49.5	30 ≤ Φ _{Nvb} (N _{vb} /1000) ≤ 160		
									x	yes	no

Tab No. 2.2 Test suspension

Test suspension N	N	V _{e1}	V _{e2}	Test suspension N ₀ (time = 0)
Φ = 49 × 10 ⁷ = lg 8.69	10 ⁻⁶	> 165	> 165	lg N ₀ = lg N/100 = lg 6.69
8.17 ≤ lg N ≤ 8.70	10 ⁻⁷	47	51	6.17 ≤ lg N ₀ ≤ 6.70
				x
				yes
				no

Tab No. 2.3 Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Aspergillus brasiliensis (niger)* ATCC 16404

Test concentration (%)*/contact time (min)/conditions	Dilution after test procedure	V _{e1}	V _{e2}	lg N _a = lg (Φ _a × 10)	lg R (lg N ₀ = lg 6.69)
100 / 5 / dirty	10 ⁰	17	18	2.24	4.45
100 / 15 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.54
100 / 30 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.54

* Product can only be tested at a concentration of 97% (RTU product) or less, as some dilution is always produced by adding the test organisms and interfering substance.

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019
Rep No: 93
Sample name: GLOBACID SF 2%
Sampled: by client
Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Batch No: 16-322A

Sampling date: 24.5.2019
Sample delivered: 29.5.2019
Testing date: 28.6. – 17.7.2019
Delivered amount: 500 ml
Page: 4

CONFIDENTIAL

3. Evaluation of fungicidal activity of the product GLOBACID SF 2%

Tab No. 3.1 The efficacy of chemical disinfectant GLOBACID SF 2% on test strains – fungicidal activity

Strain	Fungicidal activity of the product (EN 13624:2013)					
	Test temperature [°C]	Contact time [min]	Product test concentrations [%]*	Interfering substances - conditions	lg R EN 13624:2013	lg R
<i>Candida albicans</i> ATCC 10231	20	5	100	dirty	≥ 4	> 4
<i>Aspergillus brasiliensis (niger)</i> ATCC 16404	20	5	100	dirty	≥ 4	> 4
<i>Candida albicans</i> ATCC 10231	20	15	100	dirty	≥ 4	> 4
<i>Aspergillus brasiliensis (niger)</i> ATCC 16404	20	15	100	dirty	≥ 4	> 4
<i>Candida albicans</i> ATCC 10231	20	30	100	dirty	≥ 4	> 4
<i>Aspergillus brasiliensis (niger)</i> ATCC 16404	20	30	100	dirty	≥ 4	> 4

Note: V_c = value is the number of cfu per ml, Φ = average V_{c1} a V_{c2} (1. + 2. duplicate V_c values), N = the number of cfu/ml of the test suspension, N_{V0} = the number of cfu/ml of the test suspension for validation, N_a = the number of surviving fungi per ml in the test mixture, A, B, C = the number of surviving fungi per ml in control tests (A – experimental conditions validation, B – neutralizer toxicity validation, C – method validation), $R = N_0 / N_a$ = the reduction in viability, or $lg R = lg N_0 - lg N_a$

* Product can only be tested at a concentration of 97% (RTU product) or less, as some dilution is always produced by adding the test organisms and interfering substance.

Prepared by: Mgr. Karolína Světlíková, Lab Technician

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019
Rep No: 93
Sample name: GLOBACID SF 2%
Sampled: by client
Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Batch No: 16-322A

Sampling date: 24.5.2019
Sample delivered: 29.5.2019
Testing date: 28.6. – 17.7.2019
Delivered amount: 500 ml
Page: 5

Experimental conditions:

Period of analysis:
Test temperature:
Test method:
Neutralization medium:
Appearance of the product:
Test concentration:
Contact time:
Interfering substances:
Test organisms:
Incubation conditions:

Testing of disinfecting efficiency of chemical disinfecting and antiseptic agents on carriers SOP-M-22-12 (EN 14562:2006)
15.7. – 17.7.2019
20 °C ± 1 °C
dilution neutralization method
Dey-Engley Neutralizing Broth M 1062
colourless liquid
100% (concentrated)
5 min, 15 min, 30 min
3 g/l BSA and 3 ml/l sheep erythrocytes (dirty conditions)
Candida albicans ATCC 10231
Aspergillus brasiliensis (niger) ATCC 16404
30 °C ± 1 °C, 48 hours and additional period of 24 or 48 hours

Test procedure:

1. Preparation of the test suspension
2. Preparation of the product test solutions
3. Quantitative carrier test
4. Incubation and calculation
5. Expression and interpretation of the results

CONFIDENTIAL

Note:

Presence of a high concentration (at least 75%) of *Aspergillus brasiliensis* spiny spores in the test suspension – yes.

Fungicidal activity – the capability of a product to produce a reduction in the number of relevant organisms on carriers under defined conditions by at least a 4 lg reduction (10^4).

Yeasticidal activity – the capability of a product to produce a reduction in the number of viable fungi belonging to reference strain *Candida albicans* on carriers under defined conditions by at least a 4 lg reduction (10^4).

The drying time: 20-38 min

$R = N_w / N_a$ or $\lg R = \lg N_w - \lg N_a$ = the reduction in viability

The standard:

EN 14562:2006 Chemical disinfectants and antiseptics – Quantitative carrier test for the evaluation of fungicidal or yeasticidal activity for instruments used in the medical area - Test method and requirements (phase 2, step 2)
May 2006

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019

Rep No: 93

Sample name: GLOBACID SF 2%

Sampled: by client

Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404

Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404

Batch No: 16-322A

Sampling date: 24.5.2019

Sample delivered: 29.5.2019

Testing date: 28.6. – 17.7.2019

Delivered amount: 500 ml

Page: 6

4. Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Candida albicans* ATCC 10231 on carriers

Tab No. 4.1 Verification of methodology, dirty conditions

Validation of suspension (N ₀)			Validation of selected experimental conditions (A)			Neutralizer toxicity control (B)			Method validation (C) Product conc.: 100%		
V _{e1}	26	Φ _{N₀} = 31	V _{e1}	28	Φ _A = 28	V _{e1}	31	Φ _B = 28	V _{e1}	25	Φ _C = 25.5
V _{e2}	36		V _{e2}	28		V _{e2}	25		V _{e2}	26	
30 < Φ _{N₀} < 160			Φ _A > 0.5 Φ _{N₀}			Φ _B > 0.5 Φ _{N₀}			Φ _C > 0.5 Φ _{N₀}		
x	yes	no	x	yes	no	x	yes	no	x	yes	no

Tab No. 4.2 Test suspension

Test suspension (N)	N	V _{e1}	V _{e2}	Φ = 34.5 x 10 ⁷ = lg 8.54 8.17 ≤ lg N ≤ 8.70		
	10 ⁻⁶	> 330	> 330			
	10 ⁻⁷	34	35			
				x	yes	no

Tab No. 4.2.1 The control test suspension, dirty conditions

Test suspension (N _w)	N _w	V _{e1}	V _{e2}	Φ x 10 = 1250 x 10 ⁴ = lg 7.10 lg N _w = lg 7.10 6.15 ≤ lg N _w ≤ (lg N - 1.3) 7.24		
	10 ⁻⁴	131	119			
				x	yes	no

Tab No. 4.3 Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Candida albicans* ATCC 10231 on carriers, dirty conditions

Test concentration (%) / contact time (min) / conditions	Dilution after test procedure	V _{e1}	V _{e2}	lg N _s = lg (Φ _s x 10)	lg R (lg N _w = lg 7.10)
100 / 5 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.95
100 / 15 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.95
100 / 30 / dirty	10 ⁰	<14	<14	< 2.15	≥ 4.95

Note: V_e = value is the number of cfu per ml, Φ = average V_{e1} a V_{e2} (1. + 2. duplicate V_e values), N = the number of cfu/ml of the fungal test suspension, N_w = the number of cfu/ml of the control fungal test suspension, N_s = the number of surviving fungi per ml in the test mixture at the end of the contact time, N₀ = the number of cfu/ml of the fungal test suspension in the mixture A,B,C at the beginning of the contact time (time „0“), A,B,C = the number of surviving fungi per ml in control tests (A – experimental conditions validation, B – neutralizer toxicity validation, C – method validation), R = N_w / N_s or lg R = lg N_w - lg N_s = the reduction in viability

CONFIDENTIAL

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019
 Rep No: 93
 Sample name: GLOBACID SF 2%
 Sampled: by client
 Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
 Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
 Batch No: 16-322A

Sampling date: 24.5.2019
 Sample delivered: 29.5.2019
 Testing date: 28.6. – 17.7.2019
 Delivered amount: 500 ml
 Page: 7

5. Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Aspergillus brasiliensis (niger)* ATCC 16404 on carriers

Tab No. 5.1 Verification of methodology, dirty conditions

Validation of suspension (N_{t0})		Validation of selected experimental conditions (A)		Neutralizer toxicity control (B)		Method validation (C) Product conc.: 100%	
V_{c1}	47	V_{c1}	19	V_{c1}	38	V_{c1}	30
V_{c2}	22	V_{c2}	34	V_{c2}	27	V_{c2}	33
$\Phi_{N_{t0}} = 34.5$		$\Phi_A = 26.5$		$\Phi_B = 32.5$		$\Phi_C = 31.5$	
$30 \leq \Phi_{N_{t0}} \leq 160$		$\Phi_A > 0.5 \Phi_{N_{t0}}$		$\Phi_B > 0.5 \Phi_{N_{t0}}$		$\Phi_C > 0.5 \Phi_{N_{t0}}$	
x	yes	x	yes	x	yes	x	yes
	no		no		no		no

Tab No. 5.2 Test suspension

Test suspension (N)	N	V_{c1}	V_{c2}	
	10^{-6}	> 165	> 165	$\Phi = 46 \times 10^7 = \lg 8.66$ $8.17 \leq \lg N \leq 8.70$
	10^{-7}	51	41	
				x yes no

Tab No. 5.2.1 The control test suspension, dirty conditions

Test suspension (N_w)	N_w	V_{c1}	V_{c2}	
	10^{-4}	64	52	$\Phi \times 10 = 580 \times 10^4 = \lg 6.76$ $\lg N_w = \lg 6.76$ $6.15 \leq \lg N_w \leq (\lg N - 1.3) 7.36$
				x yes no

Tab No. 5.3 Testing the efficacy of chemical disinfectant GLOBACID SF 2% on *Aspergillus brasiliensis (niger)* ATCC 16404 on carriers, dirty conditions

Test concentration (%) / contact time (min) / conditions	Dilution after test procedure	V_{c1}	V_{c2}	$\lg N_s = \lg (\Phi_s \times 10)$	$\lg R$ ($\lg N_w - \lg 6.76$)
100 / 5 / dirty	10^{-1}	86	65	3.88	2.88
100 / 15 / dirty	10^{-1}	14	27	3.31	3.45
100 / 30 / dirty	10^0	27	25	2.41	4.35

Note: V_c = value is the number of cfu per ml, Φ = average V_{c1} a V_{c2} (1. + 2. duplicate V_c values), N = the number of cfu/ml of the fungal test suspension, N_w = the number of cfu/ml of the control fungal test suspension, N_s = the number of surviving fungi per ml in the test mixture at the end of the contact time, N_{t0} = the number of cfu/ml of the fungal test suspension in the mixture A,B,C at the beginning of the contact time (time „0“), A,B,C = the number of surviving fungi per ml in control tests (A – experimental conditions validation, B – neutralizer toxicity validation, C – method validation), $R = N_w / N_s$ or $\lg R = \lg N_w - \lg N_s$ = the reduction in viability

CONFIDENTIAL

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019
Rep No: 93
Sample name: **GLOBACID SF 2%**
Sampled: by client
Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404
Batch No: 16-322A

Sampling date: 24.5.2019
Sample delivered: 29.5.2019
Testing date: 28.6. – 17.7.2019
Delivered amount: 500 ml
Page: 8

6. Evaluation of fungicidal activity of the product **GLOBACID SF 2%**

Tab No. 6.1 The efficacy of chemical disinfectant **GLOBACID SF 2%** on test strains – fungicidal activity

activity	Strain	Fungicidal activity of the product (EN 14562:2006)					
		Test temperature [°C]	Contact time [min]	Product test concentrations [%]	Interfering substances - conditions	lg R EN 14562:2006	lg R
	<i>Candida albicans</i> ATCC 10231	20	5	100	dirty	≥ 4	> 4
	<i>Aspergillus brasiliensis (niger)</i> ATCC 16404	20	5	100	dirty	≥ 4	< 4
	<i>Candida albicans</i> ATCC 10231	20	15	100	dirty	≥ 4	> 4
	<i>Aspergillus brasiliensis (niger)</i> ATCC 16404	20	15	100	dirty	≥ 4	< 4
	<i>Candida albicans</i> ATCC 10231	20	30	100	dirty	≥ 4	> 4
	<i>Aspergillus brasiliensis (niger)</i> ATCC 16404	20	30	100	dirty	≥ 4	> 4

Note: V_c = value is the number of cfu per ml, Φ = average V_{c1} a V_{c2} (1. + 2. duplicate V_c values), N = the number of cfu/ml of the fungal test suspension, N_w = the number of cfu/ml of the control fungal test suspension, N_a = the number of surviving fungi per ml in the test mixture at the end of the contact time, N_{t0} = the number of cfu/ml of the fungal test suspension in the mixture A,B,C at the beginning of the contact time (time „0“), A,B,C = the number of surviving fungi per ml in control tests (A – experimental conditions validation, B – neutralizer toxicity validation, C – method validation), $R = N_w / N_a$ or $\lg R = \lg N_w - \lg N_a$ = the reduction in viability

Prepared by: Ing. Eva Kremlová, Lab Technician

CONFIDENTIAL

Description: Testing the efficacy of chemical disinfectants and antiseptics

Sample ID: S186/2019

Rep No: 93

Sample name: GLOBACID SF 2%

Sampled: by client

Sampling point: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404

Client: Goodpoint Chemicals Ltd, Urda tee 3, Jälgimäe, Saku vald, Harjumaa, Estonia 76404

Batch No: 16-322A

Sampling date: 24.5.2019

Sample delivered: 29.5.2019

Testing date: 28.6. – 17.7.2019

Delivered amount: 500 ml

Page: 9

CONFIDENTIAL

Interpretation:

Results of tests are in Tabs.

According to EN 13624:2013 the tested concentrated* product GLOBACID SF 2%, batch No. 16-322A, in the contact times 5 min, 15 min and 30 min under dirty conditions at temperature $20\text{ }^{\circ}\text{C} \pm 1\text{ }^{\circ}\text{C}$ by the dilution neutralization method decreased the number of viable yeast cells of *Candida albicans* ATCC 10231 and the number of mould spores of *Aspergillus brasiliensis* (*niger*) ATCC 16404 by at least a 4 lg reduction.

* Product can only be tested at a concentration of 97% (RTU product) or less, as some dilution is always produced by adding the test organisms and interfering substance.

The tested concentrated product GLOBACID SF 2%, batch No. 16-322A, in the contact times 5 min, 15 min and 30 min under dirty conditions at temperature $20\text{ }^{\circ}\text{C} \pm 1\text{ }^{\circ}\text{C}$ by the dilution neutralization method decreased on carriers the number of viable yeast cells of *Candida albicans* ATCC 10231 by at least a 4 lg reduction (EN 14562:2006).

The tested concentrated product GLOBACID SF 2%, batch No. 16-322A, in the contact time 30 min under dirty conditions at temperature $20\text{ }^{\circ}\text{C} \pm 1\text{ }^{\circ}\text{C}$ by the dilution neutralization method decreased on carriers the number of mould spores of *Aspergillus brasiliensis* (*niger*) ATCC 16404 by at least a 4 lg reduction (EN 14562:2006).

Conclusion:

The product GLOBACID SF 2% is capable of reducing the number of viable vegetative yeast cells and the number of mould spores of the relevant organisms under defined conditions (EN 13624:2013 – min. 97%, 5 min, 15 min and 30 min, dirty, $20\text{ }^{\circ}\text{C}$) to the declared values, and consequently, can be called fungicidal.

The product GLOBACID SF 2% is capable of reducing the number of viable vegetative yeast cells of the relevant organism on carriers under defined conditions (EN 14562:2006 – C.a., 100%, 5 min, 15 min and 30 min, dirty, $20\text{ }^{\circ}\text{C}$) to the declared values, and consequently, can be called yeasticidal on carriers.

The product GLOBACID SF 2% is capable of reducing the number of mould spores of the relevant organism on carriers under defined conditions (EN 14562:2006 – A.b., 100%, 30 min, dirty, $20\text{ }^{\circ}\text{C}$) to the declared values, and consequently, can be called fungicidal on carriers.

22.7.2019, Hodonín

Ing. Barbora Stoklásková, Leader of Study

