

MICROGAS

CF38 PURE
GAS FOR INTRAOCULAR USE



MICROGAS C3F8 PURE

MicroC3F8 (pure or mixed) is injected into the vitreous chamber for the purpose of replacing vitreous and promoting retinal adhesion. During the average residence time in the eye (about 7-8 weeks), the gas is gradually replaced by aqueous humor and exhaled.

PRODUCT ASSORTMENT

The **MicroC3F8 Puro** device is available in one variant in which C3F8 is 100% pure. Besides the gas canister, the kit consists of a Connector, a 0.2 μ Sterilizing Filter, a 60 ml PP Syringe, a 3-way Tap, a 27 G Needle and a 30 G Needle.

Pure gases should be mixed with nitrogen before use.

CODE	DESCRIPTION	NAME
MMD-788	Pre-filled canister with • C3F8 Pure Gas and Connector • 0.2 μ sterilizing filter • 60 ml PP syringe • 3-way tap • Needle 27 G and Needle 30 G	MICROC2F6 Pure

