












Email: sales@panswmed.com
 Tel: 86-20-28961025
 Skype: panmedical

PANSW NIBP Connectors

| Part No | Description | Picture | Compatibility | Original part number |
|--------------------------|-------------|---|--|----------------------|
| CNT-01 | plastic |  | all kinds of disposable nibp cuffs; Male Slip Luer to 5/32 in. ID Tube, Plastic | 300666 |
| CNT-02 | plastic |  | Nihon Kohden dual tube system interconnect hoses; other dual tube systems | GE Part # 300829 |
| CNT-03 (F) CNT-03 (M) | |  | Datex Ohmeda submin connectors Mated Pair Submin to 1/8" ID Barbed, Plastic Compatibility: Datex-Ohmeda Part # 330058 Artema NIBP tubing system Female Submin to 1 /8" ID Barbed, Plastic Compatibility: GE, Marquette Protocol Part # 330064 Male Submin to 1 /8" ID Barbed, Plastic Compatibility: GE, Marquette, Datex-Ohmeda Part # 330091 | |
| CNT-04 | |  | Fukuda Kontron, Schiller | Use with CNT-06 |

| | | | | |
|------------------------|--|--|--|-----------------|
| CNT-05 | neonate NIBP connector for all the marks neonate system |  | <p>5/32" ID Barbed to Female Slip Luer , Metal Compatibility: Part # 300669 Use with CNT-18</p> | |
| CNT-06 |  | | | Use with CNT-04 |
| CNT-07(F) CNT-07(M) |  <p>CNT-07 (F)</p>  <p>CNT-07 (M)</p> |  | <p>Compatibility: Philips/HP, Siemens, Datascope, Colin Part # 330060 Compatibility: Philips/HP, Siemens, Datascope, Colin; Mindray Part # 330059 Spacelabs 714-0045-00;</p> <p>5/32" ID Barbed to Female Bayonet, Metal Male Bayonet to 5/32" ID Barbed, Metal Spacelabs: 714-0038-00</p> | |
| CNT07FF |  | | <p>outer diameter: 5.3mm</p> | |
| CNT-08 |  | | <p>Nihon Kohden adult & pediatric</p> | |