
Standard Instruments for Rhinology and Rhinoplasty


RHINOFORCE® II Nasal Scissors and Nasal Cutting Forceps





449202	Same, curved to right
449203	Same, curved to left
449204	RHINOFORCE® II Nasal Scissors, cutting blades curved upwards, to cut the anterior part of the middle turbinate in case of narrow nose and to prevent postoperative synechia, with cleaning connector, working length 13 cm
449205	Same, cutting blades curved downwards

STRUYCKEN RHINOFORCE® II Nasal Cutting Forceps, working length 13 cm



Recommended Cleaning Rack 39220 CA for Instruments with Cleaning Connectors see page 28

11-983