## CLEANASCOPET ADVANTAGE

Transport & Short-term Storage System



## CLEAR IDENTIFICATION OF CLEAN VERSUS CONTAMINATED ENDOSCOPES

Know it's Clean at a Glance. The CLEANASCOPE™ ADVANTAGE Transport and Short-term Storage System helps to protect your endoscopes by minimizing direct handling and reducing contact with airborne contaminants.

The long established design creates a protective environment for each disinfected endoscope. Green (clean) and red (contaminated) color-coded liners immediately identify clean versus contaminated scopes.



## OPTIMIZED ENDOSCOPE WORKFLOW



<sup>&</sup>quot;Transport the soiled endoscope to the reprocessing area in a closed, puncture-resistent container that prevents exposing staff, patients, or the environment to potentially infectious organisms..."

Trays fit into the contoured cart top, creating an ideal work surface

Provides enclosed short-term storage of reprocessed scopes

Hard lids available for added protection when transporting outside department

Available in 5, 6 and 10 tray options









MANUAL CLEANING



CASSETTE FITS SECURELY IN ADVANTAGE PLUS™ REPROCESSOR

## ${\sf CLEANASCOPE^{\tiny{\sf TM}}}~{\sf ADVANTAGE}~{\sf System}~{\sf Ordering}~{\sf Information}$

CLEANASCOPE™ CART		
PART NUMBER	DESCRIPTION	DIMENSIONS
103470	Five tray cart (trays included)	35.8"H x 25.1"W x 22.4"D (in)
103466	Six tray cart (trays included)	41.7"H x 25.1"W x 22.4"D (in)
103467	Ten tray cart (trays included)	64.6"H x 25.1"W x 22.4"D (in)







64.6"

CLEANASCOPE™ LINERS & ACCESSORIES			
PART NUMBER	DESCRIPTION	UNITS PER BOX	
103585	Sterile Liner Set	200 Sterile liners/case	
103604	Dispenser, (wall mounted)	1	
103456	Tray	1	
103457	Tray Lid, clear	1	
CAS-1000	Advantage Cassette	1	

