



THUASNE

Breastcare overview

SILIMA & Serena

Breast forms and lingerie following mamma ablatio

2024 - THUASNE DEUTSCHLAND GmbH



Confidentiality

This document is the property of Thuasne®, is for the exclusive attention of its intended recipients and the information it contains is strictly confidential.

No part of this document may be disclosed, reproduced or transmitted by any means whether, electronic, mechanical, photocopying or otherwise, without the prior written authorization of the company's management. If you have received the present document in error by e-mail, please notify us immediately by responding to this e-mail and then delete the present document as well as the concerned e-mail from your system.

THUASNE® SOLUTION
SILIMA breast forms
and lingerie



SILIMA® soft & light super soft

Natural and particularly smooth-fitting



66 379 – symmetric
Standard A/B Size 1 – 11

- **Super soft cup nestles closely in the bra cup** and thus resembles even more to particularly soft breast tissue
- Adapts to fit the thoracic wall with ease, ideally suited for sensitive scar tissue
- For an optimum weight **compensation 15 % lighter** than the remaining breast
- Adapts naturally to all body movements, whether standing, sitting or lying down
- Lightweight silicone core encased in particularly soft silicone

SILIMA® soft & light –Especially soft lightweight breast forms



Without lateral extension:
66 378 – symmetric
Standard B Gr. 1 – 11
Full C Gr. 3 – 14

- Anatomically shaped, 3-layer breast form
- The **extra soft material** ensures a perfect fit and natural drape, ideal for softer breast tissue and larger breasts
- It naturally follows the wearer's rhythm and movement, even when lying down
- Soft and conformable layer of silicone next to the body
- Perception of coolness and comfort, **30% weight reduced**

SILIMA® direct – Lightweight adhesive form




Rear view with and without
adhesive pad

66 377 - symmetrical
Standard B size 1 – 11
Full C size 3 – 12

- Lightweight breast form - **up to 30 % lighter** vs. standard
- Concept: **2 in 1 - Can be worn with or without adhesive pad**
- Adhesion properties are maintained through a simple cleaning procedure
- Cuddly, flexible adhesive pad, adapts optimally to the scar area for a snug fit between breastform and body – feels like a ‚second skin‘
- Specially designed inside with an edge recess for **easy and secure placement of the adhesive pad**

SILIMA® ultra light – Extra-lightweight form



66 374 – symmetric
Standard B Size 1 – 11
Full C Size 3 – 14

- Extra light breast form: **60% lighter** than standard silicone gel breast forms
- Ideal for lymphedema patients (arm or thorax) and
- Deeply **hollowed underside to protect the sensitive scar area**

- **Unique design with concentric ring feature recreates natural body movement** ①
- Bottom edge finish ensures good ventilation. When worn in water, drainage is effortlessly achieved ②
- Perfect fit ③, Large open shell structure relieves sensitive scar tissue ④

SILIMA® soft & light – Extra soft lightweight breast form – Heartshape




66 328 – symmetric
Standard B Size 1 – 8
Full C Size 5 - 12

- Symmetric shape – **for left and right or bilateral fittings**
- Versatile in use and easy to fit
- Up **to 30% lighter** than conventional gel-forms
- The curved outline towards the center ensures that it does not touch the natural breast
- Gently tapers to the neckline and underarm area

SILIMA® elite – Soft lightweight form

- Natural compensation for **larger breasted women**
- Available in D/E cup in four double sizes: 8, 10, 12, 14
- Symmetrically designed for left and right fitting
- Anatomical shape



66 398 – symmetric
Extra Full D/E Size 8/9 – 14/15

SILIMA® classic – standard gel form (Basic care)

- Consist of standard weighted silicone gel and thus nearly match the weight of the natural breast. particularly natural and have a silky soft surface, resistant to sea water and chlorine.
- They adjust anatomically to the natural rounding of the ribcage without wrinkling



66 370 – symmetric
Standard B size 1 – 11
Full C size 3 – 14

SILIMA® soft & light – Especially soft lightweight breast forms



With lateral extension:
66 368 – asymmetric
A Cup Size 3 – 7
Standard B Size 3 – 10
Full C Size 5 – 12

- Anatomically shaped, 3-layer breast form
- Extra soft material **ensures a perfect fit and natural drape, ideal for softer breast tissue and larger breasts**
- It naturally follows the wearer's rhythm and movement, even when lying down
- Soft and conformable layer of silicone next to the body
- Perception of coolness and comfort, **30% weight reduced**

SILIMA® – FOR EVEN MORE NATURALNESS

Product	light tinted	medium tinted	dark tinted
SILIMA Fiberfill triangular	x		
SILIMA classic asymmetric	x	x	x
SILIMA classic symmetric, Art. 66370	x	x	x
SILIMA soft&light, Super soft	x	x	x
SILIMA elite	x	x	x
SILIMA direct	x	x	x
SILIMA soft&light asymmetric	x	x	x
SILIMA soft&light symmetric	x	x	x
SILIMA soft&light Heartshape	x	x	x
SILIMA ultra light	x	x	x
SILIMA shell	x	x	x
SILIMA conform	x	x	x
SILIMA Xtra	x		
Serena	x		



**Available in three colors:
light, medium, dark tinted**

INSTRUCTIONS FOR USE (IFU)



- ✓ Instructions for use with
- ✓ indication, precautions,
- ✓ Dressing instruction,
- ✓ Care instructions and material composition
- ✓ Weight information per bf and size

Gewichte der SILIMA® Epithesen (in Gramm - Stand 2021-11) ¹

Art.-Nr. ²	Bezeichnung/Name ³	Cups ⁴	Größen ⁵				
			0	1	2	3	4
T66379	SILIMA® soft & light	A/B		85	102	155	190
T66328	SILIMA® soft & light	B		60	90	130	170
T66328	SILIMA® soft & light	C					
T66378	SILIMA® soft & light	B		70	100	150	185
T66378	SILIMA® soft & light	C				185	225
T66368	SILIMA® soft & light	A				100	150
T66368	SILIMA® soft & light	B				155	220
T66368	SILIMA® soft & light	C					
T66398	SILIMA® elite	D/E					
T66377	SILIMA® direct (ohne Pad) ⁶	B		75	90	125	160
T66377	SILIMA® direct (ohne Pad) ⁶	C				160	195
T66374	SILIMA® ultra light	B		50	65	100	105
T66374	SILIMA® ultra light	C				110	125
T66375*	SILIMA® light	B				140	160
T66375*	SILIMA® light	C				170	205
T66370	SILIMA® classic	B		100	125	180	225
T66370	SILIMA® classic	C				230	290
T66360	SILIMA® classic	A		90	120	145	195
T66360	SILIMA® classic	B		120	160	200	280
T66360	SILIMA® classic	C					
T66307	SILIMA® extra			12	27	55	97

SILIMA® LINGERIE - SEGMENTS



SILIMA® Daily



SILIMA® Support



SILIMA® Lace



SILIMA® Medical

DAILY BRA LUCILLE



Bra Lucille

- Sizes: 75A-100A; 75-105 B, C; 80-100D
- colors: white and nude
- Art.-Nr. T57 580

- Seamless molded cup (molded),
- Wide padded straps,
- Good coverage back and front area provide high support
- → The comfortable "Allrounder"



DAILY BRA DIANA



Bra Diana

Classic design coupled with good support and coverage



- Skin friendly material out of microfibre with two-section cups
- Decorative lace band across the neckline add elegant flair
- Adjustable straps with shoulder padding
- Colors: skin, white, black, raspberry
- Sizes: 70 – 120 A – B, 75 – 110 C - D
- Ref. 57 255

Bra Flora

Exquisite bra with loving attention to detail



- Soft and smooth stretch lace with floral design, fully lined
- Soft shoulder padding, adjustable straps
- Separate lace inserts available for extra security, camisole effect

- Colors: Champaign, white, black
- Sizes: 70 – 95 AA – C, 70 – 90 D
- Ref. 57 208
- Separate lace inserts available ref. 64 207

**Thuasne
solution
Serena**

SERENA S

Thuasne solution - Serena®

Serena® comprises a variety of normal weight silicone breast forms and attractive lingerie.

Soft and skin friendly material have been used to fulfil the women's needs.

Serena® breast forms are particularly natural and have a silky soft surface, resistant to sea water and chlorine, allows the women after mastectomy to recover their natural comfort with full confidence.

Completed by attractive special purpose bras, designed to be worn with Serena® breast forms, holding them securely in place.



Breast forms Serena®

Serena® triangular - 9201: symmetric / anatomical form

symmetric full volume compensation, use left and right, Sizes: Cup 1-10



Serena® Heart - 9341: heart shape,

symmetric full volume compensation, slight lateral extensions, Sizes: Cup 1-15



Serena® oval - in two versions: 9301 round cavity or 9311 oval cavity

Symmetric with underarm extension to compensate tissue deficits,

versatile in use, Sizes: Cup 1-13



Lingerie Serena®

Victoria - 67100

Excelent support and maximum comfort with delicate lace

Colors: white, beige, black, **Sizes:** 75-105 A-E



Natalia - 67200

Delicate embroidery, three section cup provide firm support, even in larger sizes

Colors: white and beige, **Sizes:** 75-100 B-E



**CLEANING
CARE
STORAGE**



CLEANING AND CARE

Please...

- Wash lingerie in the care programme (with reduced spin cycle)
- Do not bleach do not tumble dry
- Do not iron
- Do not dry clean



Please do not use fabric softener
A laundry net is helpful



CLEANING AND CARE

Breast forms

- Wash with **mild soap and rise in warm water**
- Do use a mild cleansing lotion and rinse with water
- Dab gently with a **soft towel to dry**
- Do not place the form on a radiator, nor use a hair dryer
- Keep sharp or pointed away from the breastform
- Be careful with pet claws – they might damage the gossamer-thin PU film
- Use the **original box to store** or move the breast forms



SHELF LIFE AND STORAGE CONDITIONS

Storage conditions

- Should be stored in a clean, dry environment
- Temperature between 10 °C and 30 °C

Thank you!

Credits

Sylvain Leurent

Shutterstock

Studio Caterin

Wyrwa