

www.balton.p







Урология



The cystostomy kit is used for the performance of suprapubic cystostomy by means of puncture. Catheter made of superior quality polyurethane. This kind of material ensures optimally long time of maintenance in place.

Reference:	Catheter size	Catheter length [cm]	Split cannula size	Peel away needle length [mm]
ZCYS8F	8F	45	9F	80
ZCYS9F	9F		9F	80
ZCYS10F	10F		12F	120
ZCYS12F	12F		12F	120
ZCYS14F	14F		14F	120

When ordering use Ref. Number, e.g.: ZCYS12F.



The content of this document, including provided data, information, trademarks and logos are the sole property of BALTON. Any reproduction, whether in whole or in part, is strictly prohibited and will be considered a violation of BALTON's copyrights and other intellectual proprietary rights. ©2022 Balton Sp. z o.o. All rights reserved. This document and any images contained therein is intended solely for the use by healthcare professionals and BALTON's distributors.

