


Product Specifications

Tissue-Tek[®] Coverslips



Tissue-Tek[®] SF Glas[™] Coverslips

Product description

| Item Code | Product Description | Sales quantity |
|-----------|--|----------------|
| 9580 | SF Glas™ Coverslips 24 x 40 mm | 10 x 100/pcs |
| 9581 | SF Glas [™] Coverslips 24 x 50 mm | 10 x 100/pcs |
| 9579 | SF Glas™ Coverslips 24 x 55 mm | 10 x 100/pcs |
| 9582 | SF Glas [™] Coverslips 24 x 60 mm | 10 x 100/pcs |

Product specifications

| 1. Application | |
|----------------|--|
| | Glas coverslipping of (stained) slides |
| | |

| 2. Related Instrument | |
|-----------------------|---|
| Item Code | Product Description |
| 6502 | Tissue-Tek [®] Glas [™] g2 Coverslipper |

| 3. Regulatory information | on | | |
|---------------------------|--|----|--|
| Intended Use | Coverslipping tissue or the tissue will not be d | | las coverslips to ensure I factors. |
| CE status | 98/79/EC - CE IVD | | |
| MSDS available | Yes 🗆 | No | \boxtimes |
| | | | |
| | | | |

| 4. Packaging | | |
|-----------------------------|---|--|
| Package description | Packed per 100 pieces in plastic boxes. Packing unit 10x100 pieces, packed in plastic to extra avoid moisture. | |
| | Coverslips normally come with separating agent. This prevents the coverslips from adhering and therefore has a positive effect on the use in slipping machines. The separating agent has no influence on the assay results. | |
| Sales quantity | 10 boxes of 100/pcs | |
| Dimensions (LxWxH) | 75 x 35 x 30 mm | |
| Package dimensions (LxWxH) | 200 x 150 x 30 mm | |
| Shipping dimensions (LxWxH) | 200 x 150 x 30 mm | |
| Shipping weight (g) | 0.75 | |

| 5. Labeling | | |
|------------------|---------------------------|--|
| HIBC code | +H905 <i>Item code</i> 40 | |
| EAN code | | |
| Statistical code | 70109091 | |
| | | |



Label item 9580

| | mperature during storage and usage. cause condensation and lead to between the glasses. |
|-----------------|--|
| Dry environment | |
| Yes 🛛 | No 🗆 |
| Yes 🗆 | No 🛛 |
| | Avoid large variations in technologie Cooling of the product can condensed water forming technologie Dry environment Yes |

| 7. Material specification | IS | | |
|---------------------------|---|--|--|
| Components | D 263 M by Schott AG, colourless borosilicate glass. | | |
| | Chemical composition: $SiO_2 64.1\%$ $B_2O_3 8.4\%$ $Al_2O_3 4.2\%$ $Na_2O 6.4\%$ $K_2O 6.9\%$ ZnO 5.9% $TiO_2 4.0\%$ $Sb_2O_3 0.1\%$ | | |
| | Chemical properties of D 263 M: Hydrolytic resistance (DIN ISO 719) Hydrolytic class: HGB 1 Equivalent of alkali (Na2O) per gram of glass grains in µg/g: 20 | | |
| | Acid resistance (DIN 12 116): Acid class: S 2 Half surface weight loss after 6 hours mg/dm ² : 1.4 | | |
| | Alkali resistance (DIN ISO 695): Class: A 2 Surface weight loss after 3 hours mg/dm ² : 88 | | |
| Size | Item code: L x W x Thickness: | | |
| | 9580 L24 x W40 mm Thickness 0.13 - 0.16 mm 9581 L24 x W50 mm Thickness 0.13 - 0.16 mm 9579 L24 x W55 mm Thickness 0.13 - 0.16 mm 9582 L24 x W40 mm Thickness 0.13 - 0.16 mm | | |
| Refractive index | ηe (λ = 546.1 nm): 1.5255 ± 0.0015 ηD (λ = 589.3 nm): 1.5230 Abbe value ve: 55 | | |
| Light transmission | The light transmittance for a thickness of 0.15 mm is tVD65 in % $(d=0.15 \text{ mm}) = 91.7 \pm 0.3\%$. In the visible range of the spectrum D 263 M is without absorption. The excellent UV absorption properties make D 263 M an ideal material for use in fluorescence microscopy. | | |
| Cleaning | Single use only | | |
| Latex free | Yes* No Image: Constraint of the second sec | | |

| 6 months after receiving by buyer | |
|-----------------------------------|--|
| F.O.B. Warehouse, The Netherlands | |
| | |
| | |
| | |
| | |

9. Additional information