

ZENFLOW-14

Usable Length 100cm	Usable Length 130cm	Usable Length 140cm	Usable Length 150cm	Usable Length 180cm	Balloon Diameter [mm]	Balloon Length [mm]	NP [mm]	RBP [mm]
ZW100015002014B	—	—	—	ZW15001501514B	1.25	15	6	14
	—	—	—	ZW15001501514B	1.50	15	6	14
	—	—	—	ZW150015002014B	1.50	20	6	16
	ZW130015002014B	—	—	ZW180015002014B	1.50	—	—	—
	ZW130020002014B	—	—	—	2.00	—	—	—
	ZW130020002014B	—	—	—	2.50	—	—	—
	ZW130030002014B	—	—	—	3.00	—	—	—
	ZW130030002014B	—	—	—	3.50	—	—	—
	ZW130040002014B	—	—	—	4.00	—	—	—
	ZW130050002014B	—	—	—	5.00	—	—	—
ZW100015004014B	—	—	—	ZW150015004014B	1.50	—	—	—
	—	—	—	ZW150020004014B	2.00	—	—	—
	—	—	—	ZW150020004014B	2.50	—	—	—
	—	—	—	ZW150030004014B	3.00	—	—	—
	—	—	—	ZW150030004014B	3.50	—	—	—
	—	—	—	ZW150040004014B	4.00	—	—	—
	—	—	—	ZW150050004014B	5.00	—	—	—
	—	—	—	—	1.50	—	—	—
	—	—	—	—	2.00	—	—	—
	—	—	—	—	2.50	—	—	—
ZW100020008014B	—	—	—	—	3.00	60	6	16
	—	—	—	—	3.50	—	—	—
	—	—	—	—	4.00	—	—	—
	—	—	—	—	5.00	—	—	—
	—	—	—	—	1.50	—	—	—
	—	—	—	—	2.00	—	—	—
	—	—	—	—	2.50	—	—	—
	—	—	—	—	3.00	—	—	—
	—	—	—	—	3.50	—	—	—
	—	—	—	—	4.00	—	—	—
ZW100025008014B	—	—	—	—	5.00	—	—	—
	—	—	—	—	1.50	—	—	—
	—	—	—	—	2.00	—	—	—
	—	—	—	—	2.50	—	—	—
	—	—	—	—	3.00	—	—	—
	—	—	—	—	3.50	—	—	—
	—	—	—	—	4.00	—	—	—
	—	—	—	—	5.00	—	—	—
	—	—	—	—	1.50	—	—	—
	—	—	—	—	2.00	—	—	—
ZW100030008014B	—	—	—	—	2.50	—	—	—
	—	—	—	—	3.00	—	—	—
	—	—	—	—	3.50	—	—	—
	—	—	—	—	4.00	—	—	—
	—	—	—	—	5.00	—	—	—
	—	—	—	—	1.50	—	—	—
	—	—	—	—	2.00	—	—	—
	—	—	—	—	2.50	—	—	—
	—	—	—	—	3.00	—	—	—
	—	—	—	—	3.50	—	—	—
ZW100035008014B	—	—	—	—	4.00	—	—	—
	—	—	—	—	5.00	—	—	—
	—	—	—	—	1.50	—	—	—
	—	—	—	—	2.00	—	—	—
	—	—	—	—	2.50	—	—	—
	—	—	—	—	3.00	—	—	—
	—	—	—	—	3.50	—	—	—
	—	—	—	—	4.00	—	—	—
	—	—	—	—	5.00	—	—	—
	—	—	—	—	1.50	—	—	—
ZW100040008014B	—	—	—	—	2.00	—	—	—
	—	—	—	—	2.50	—	—	—
	—	—	—	—	3.00	—	—	—
	—	—	—	—	3.50	—	—	—
	—	—	—	—	4.00	—	—	—

Zhejiang Zylox Medical Device Co., Ltd
BLDG 1, No.18 Technology AVE., Wenyi West Road, Hangzhou, P. R. China, 311121
TEL: (+86) 571 8861 0082 FAX: (+86) 571 8861 0608 info@zyloxmedical.com

Zhejiang Zyxlo Medical Device Co., Ltd
BLDG 1, No.18 Technology AVE., Wenyi West Road, Hangzhou, P. R. China, 311121
TEL: (+86) 571 8861 0082 FAX: (+86) 571 8861 0608 info@zyxlomedical.com



Easy Access to The Most Challenging Anatomies



ZENFLOW™
PTA Balloon Catheter



zy03028-EN Rev02

Care Vision

Company Information

Zylox Medical Device Co., Ltd was founded in the beautiful city of Hangzhou in 2012 by a group of highly skillful and experienced returnees from the US and Europe. They have on average 5-20 years of working experiences in various functional areas such as research and development, management, and manufacturing in leading device makers including Johnson & Johnson, Medtronic, and Abbott before returning to China.

Since its inception, Zylox has developed an increasingly rich product pipeline of interventional and implantable peripheral vascular devices such as PTA balloons, PTA drug-coated Balloons (PTA DCB), superficial femoral artery (SFA) stent, drug-eluting SFA stent, and other vascular accessories. Zylox is also in the process of developing other sophisticated vascular devices such as covered stents and vena cava filters.

Zylox is certified by TUV according to ISO 13485 standard. Many of its products including PTA balloons and SFA stent are also certified with CE Marks. Zylox strives to continuously innovate and provide patients with high quality medical devices.

Innovate Mission Quality

ZENFLOW-35

Excellent 6 wings Folding and Re- folding



guarantee low crossing profile and
with durable hydrophilic coating
optimize balloon crossability


Dual-lumen Shaft Design



minimizes shaft OD while maintain-
ing a large inner lumen for quick
inflation and deflation times

Score-Flex Bal- loon Material

reduces balloon wall thickness
without compromising RBP and
enhances kink resistance



Product Specification

Catheter Type	Over-The-Wire (OTW)
Guidewire	0.035"
Catheter Tip	Flexible, short and cone-shaped
Balloon Material	Nylon
Folding	Three/six wings
Coating	Hydrophilic coating
Marker	2
Balloon Size	Diameter: 2.0~12.0mm ; length: 20~220mm
Shaft	4~8F
Usable length	75, 100, 110 and 130cm

Care
Vision
Innovate
Mission
Quality

Company Information

Zylox Medical Device Co., Ltd was founded in the beautiful city of Hangzhou in 2012 by a group of highly skillful and experienced returnees from the US and Europe. They have on average 5-20 years of working experiences in various functional areas such as research and development, management, and manufacturing in leading device makers including Johnson & Johnson, Medtronic, and Abbott before returning to China.

Since its inception, Zylox has developed an increasingly rich product pipeline of interventional and implantable peripheral vascular devices such as PTA balloons, PTA drug-coated Balloons (PTA DCB), superficial femoral artery (SFA) stent, drug-eluting SFA stent, and other vascular accessories. Zylox is also in the process of developing other sophisticated vascular devices such as covered stents and vena cava filters.

Zylox is certified by TUV according to ISO 13485 standard. Many of its products including PTA balloons and SFA stent are also certified with CE Marks. Zylox strives to continuously innovate and provide patients with high quality medical devices.

ZENFLOW-35

Excellent 6 wings
Folding and Re-
folding



guarantee low crossing profile and
with durable hydrophilic coating
optimize balloon crossability


Dual-lumen Shaft
Design



minimizes shaft OD while maintain-
ing a large inner lumen for quick
inflation and deflation times

Score-Flex Bal-
loon Material

reduces balloon wall thickness
without compromising RBP and
enhances kink resistance



Product Specification

Catheter Type	Over-The-Wire (OTW)
Guidewire	0.035"
Catheter Tip	Flexible, short and cone-shaped
Balloon Material	Nylon
Folding	Three/six wings
Coating	Hydrophilic coating
Marker	2
Balloon Size	Diameter: 2.0~12.0mm ; length: 20~220mm
Shaft	4~8F
Usable length	75, 100, 110 and 130cm

ZENFLOW-18

Low Entry Profile and Excellent Flexibility


conform the product to various anatomies and natural vessel tortuosity

Low Profile Introducer Sheath

results in lower possibility of puncture trauma


Specially-processed Delivery System

reduces kinking, achieves a greater pushability and allows you to treat your patients more effectively



Product Specification

Catheter Type	Over-The-Wire (OTW)
Guidewire	0.018"
Catheter Tip	Flexible, short and cone-shaped
Balloon Material	Nylon
Folding	Three/six wings
Coating	Hydrophilic coating
Marker	2
Balloon Size	Diameter: 2.0~9.0mm ; Length: 20~220mm
Shaft	4~6F
Usable Length	45, 75, 90, 130, 150 and 180cm



ZENFLOW-14

The Broadest Range of Sizes

provides the most comprehensive solution for below-the-knee(BTK) lesions

Advanced Crossability and Pushability

contribute to outstanding balloon performances


Long lengths up to 220mm

The design of long tapered tip enables you to achieve satisfactory outcomes for the most difficult distal lesions



Product Specification

Catheter Type	Over-The-Wire (OTW)
Guidewire	0.014"
Catheter Tip	Flexible, short and cone-shaped
Balloon Material	Nylon
Folding	Three wings
Coating	Hydrophilic coating
Marker	2
Balloon Size	Diameter: 1.25-5.0mm; Length: 1.5-220mm
Shaft	4F
Usable Length	100, 130, 140, 150 and 180cm



Ordering Information

Usable Length 45cm	Usable Length 75cm	Usable Length 90cm	Usable Length 130cm	Usable Length 150cm	Usable Length 180cm	Baloon Diameter [mm]	Baloon Length [mm]	NP [mm]	RBR [mm]
ZM4550000018B ZM455000004018B ZM455000006018B ZM455000008018B ZM455000010018B ZM455000012018B ZM455000014018B ZM455000016018B ZM455000018018B ZM455000020018B ZM455000022018B ZM455000024018B ZM455000026018B ZM455000028018B ZM455000030018B ZM455000032018B	ZM755020002018B	--	ZM130020002018B	--	--	2.00	--	--	16
	ZM755020002018B	--	ZM130020002018B	--	--	2.50	--	--	16
	ZM755030002018B	--	ZM130030002018B	--	--	3.00	--	--	16
	ZM755030002018B	--	ZM130030002018B	--	--	3.50	--	6	16
	ZM755040002018B	--	ZM130040002018B	--	--	4.00	--	--	16
	ZM755050002018B	--	ZM130050002018B	--	--	5.00	--	--	14
	ZM755060002018B	--	ZM130060002018B	--	--	6.00	--	--	14
	ZM755020004018B	ZM090020004018B	ZM130020004018B	--	ZM180020004018B	2.00	--	--	16
	ZM755020004018B	ZM090025004018B	ZM130025004018B	--	ZM180025004018B	2.50	--	--	16
	ZM755030004018B	ZM090030004018B	ZM130030004018B	--	ZM180030004018B	3.00	--	--	16
	ZM755030004018B	--	ZM130030004018B	--	--	3.50	--	--	16
	ZM755040004018B	--	ZM130040004018B	--	--	4.00	--	6	16
	ZM755050004018B	--	ZM130050004018B	--	--	5.00	--	40	14
	ZM755060004018B	--	ZM130060004018B	--	--	6.00	--	--	14
	ZM755070004018B	--	ZM130070004018B	--	--	7.00	--	--	12
	ZM755080004018B	--	ZM130080004018B	--	--	8.00	--	--	12
ZM455000006018B ZM455000008018B ZM455000010018B ZM455000012018B ZM455000014018B ZM455000016018B ZM455000018018B ZM455000020018B ZM455000022018B ZM455000024018B ZM455000026018B ZM455000028018B ZM455000030018B ZM455000032018B ZM455000034018B ZM455000036018B	ZM755020006018B	ZM090020006018B	ZM130020006018B	--	--	9.00	--	--	11
	ZM755020006018B	ZM090025006018B	ZM130025006018B	--	--	9.00	--	--	11
	ZM755030006018B	ZM090030006018B	ZM130030006018B	--	--	9.00	--	--	11
	ZM755030006018B	--	ZM130030006018B	--	--	9.00	--	--	11
	ZM755040006018B	--	ZM130040006018B	--	--	9.00	--	--	11
	ZM755050006018B	--	ZM130050006018B	--	--	9.00	--	--	11
	ZM755060006018B	--	ZM130060006018B	--	--	9.00	--	--	11
	ZM755070006018B	--	ZM130070006018B	--	--	9.00	--	--	11
	ZM755080006018B	--	ZM130080006018B	--	--	9.00	--	--	11
	ZM755020008018B	ZM090020008018B	ZM130020008018B	--	ZM180020008018B	2.00	--	--	15
	ZM755020008018B	ZM090025008018B	ZM130025008018B	--	ZM180025008018B	2.50	--	--	15
	ZM755030008018B	ZM090030008018B	ZM130030008018B	--	ZM180030008018B	3.00	--	--	15
	ZM755030008018B	--	ZM130030008018B	--	--	3.50	--	6	15
	ZM755040008018B	--	ZM130040008018B	--	--	4.00	--	80	15
	ZM755050008018B	--	ZM130050008018B	--	--	5.00	--	--	13
	ZM755060008018B	--	ZM130060008018B	--	--	6.00	--	--	12
ZM455000010018B ZM455000012018B ZM455000014018B ZM455000016018B ZM455000018018B ZM455000020018B ZM455000022018B ZM455000024018B ZM455000026018B ZM455000028018B ZM455000030018B ZM455000032018B ZM455000034018B ZM455000036018B ZM455000038018B ZM455000040018B	ZM755020010018B	--	ZM130020010018B	--	--	7.00	--	--	11
	ZM755020010018B	--	ZM130020010018B	--	--	7.00	--	--	11
	ZM755030010018B	--	ZM130030010018B	--	--	7.00	--	--	11
	ZM755030010018B	--	ZM130030010018B	--	--	7.00	--	--	11
	ZM755040010018B	--	ZM130040010018B	--	--	7.00	--	--	11
	ZM755050010018B	--	ZM130050010018B	--	--	7.00	--	--	11
	ZM755060010018B	--	ZM130060010018B	--	--	7.00	--	--	11
	ZM755070010018B	--	ZM130070010018B	--	--	7.00	--	--	11
	ZM755080010018B	--	ZM130080010018B	--	--	7.00	--	--	11
	ZM755020012018B	ZM090020012018B	ZM130020012018B	--	ZM180020012018B	2.00	--	--	15
	ZM755020012018B	ZM090025012018B	ZM130025012018B	--	ZM180025012018B	2.50	--	--	15
	ZM755030012018B	ZM090030012018B	ZM130030012018B	--	ZM180030012018B	3.00	--	--	15
	ZM755030012018B	--	ZM130030012018B	--	--	3.50	--	6	15
	ZM755040012018B	--	ZM130040012018B	--	--	4.00	--	100	15
	ZM755050012018B	--	ZM130050012018B	--	--	5.00	--	--	15
	ZM755060012018B	--	ZM130060012018B	--	--	6.00	--	--	13
ZM455000012018B ZM455000014018B ZM455000016018B ZM455000018018B ZM455000020018B ZM455000022018B ZM455000024018B ZM455000026018B ZM455000028018B ZM455000030018B ZM455000032018B ZM455000034018B ZM455000036018B ZM455000038018B ZM455000040018B ZM455000042018B	ZM755020012018B	--	ZM130020012018B	--	--	7.00	--	--	12
	ZM755020012018B	--	ZM130020012018B	--	--	7.00	--	--	12
	ZM755030012018B	--	ZM130030012018B	--	--	7.00	--	--	12
	ZM755030012018B	--	ZM130030012018B	--	--	7.00	--	--	12
	ZM755040012018B	--	ZM130040012018B	--	--	7.00	--	--	12
	ZM755050012018B	--	ZM130050012018B	--	--	7.00	--	--	12
	ZM755060012018B	--	ZM130060012018B	--	--	7.00	--	--	12
	ZM755070012018B	--	ZM130070012018B	--	--	7.00	--	--	12
	ZM755080012018B	--	ZM130080012018B	--	--	7.00	--	--	12
	ZM755020014018B	ZM090020014018B	ZM130020014018B	--	ZM180020014018B	2.00	--	--	15
	ZM755020014018B	ZM090025014018B	ZM130025014018B	--	ZM180025014018B	2.50	--	--	15
	ZM755030014018B	ZM090030014018B	ZM130030014018B	--	ZM180030014018B	3.00	--	--	15
	ZM755030014018B	--	ZM130030014018B	--	--	3.50	--	6	15
	ZM755040014018B	--	ZM130040014018B	--	--	4.00	--	--	15
	ZM755050014018B	--	ZM130050014018B	--	--	5.00	--	--	13
	ZM755060014018B	--	ZM130060014018B	--	--	6.00	--	--	13
ZM455000020018B ZM455000022018B ZM455000024018B ZM455000026018B ZM455000028018B ZM455000030018B ZM455000032018B ZM455000034018B ZM455000036018B ZM455000038018B ZM455000040018B ZM455000042018B ZM455000044018B ZM455000046018B ZM455000048018B ZM455000050018B	ZM755020020018B	ZM090020020018B	ZM130020020018B	--	--	6.00	--	--	12
	ZM755020020018B	ZM090025020018B	ZM130025020018B	--	ZM180025020018B	2.00	--	--	15
	ZM755030020018B	ZM090030020018B	ZM130030020018B	--	ZM180030020018B	2.50	--	--	15
	ZM755030020018B	--	ZM130030020018B	--	--	3.00	--	6	15
	ZM755040020018B	--	ZM130040020018B	--	--	3.50	--	--	15
	ZM755050020018B	--	ZM130050020018B	--	--	4.00	--	--	15
	ZM755060020018B	--	ZM130060020018B	--	--	5.00	--	--	13
	ZM755070020018B	--	ZM130070020018B	--	--	6.00	--	--	13
	ZM755080020018B	--	ZM130080020018B	--	--	6.00	--	--	12
	ZM755020022018B	ZM090020022018B	ZM130020022018B	--	ZM180020022018B	2.00	--	--	15
	ZM755020022018B	ZM090025022018B	ZM130025022018B	--	ZM180025022018B	2.50	--	--	15
	ZM755030022018B	ZM090030022018B	ZM130030022018B	--	ZM180030022018B	3.00	--	--	15
	ZM755030022018B	--	ZM130030022018B	--	--	3.50	--	6	15
	ZM755040022018B	--	ZM130040022018B	--	--	4.00	--	--	15
	ZM755050022018B	--	ZM130050022018B	--	--	5.00	--	--	13
	ZM755060022018B	--	ZM130060022018B	--	--	6.00	--	--	12