



**HAPPYcover<sup>®</sup>**  


**PAD**

DATA SHEET PAGE 1:2


<b>TRADING NAME</b>	<b>HAPPY COVER<sup>®</sup> PAD</b>
<b>DESCRIPTION</b>	Polyethylene BEDPAN COVER, for an easy and safe seal of its opened upper side, after use. Inside the bag there is a superabsorbent molecule powder named SGAIEL <sup>®</sup> that transforms liquids into solids. This powder is enclosed in a water-soluble film which acts in contact with liquids. Ready for use.
<b>ITEM CODE</b>	AM001N
<b>RDM (Italy)</b>	1570461
<b>CND (Italy)</b>	T0399
<b>PACKAGING</b>	pcs 20
<b>DIMENSIONS</b>	mm. 500X600
<b>EXPIRING TIME</b>	5 years
<b>BELONGING CLASS OF MD</b>	Class I medical device - Not sterile Pursuant to 93/42/EEC
<b>CE MARK</b>	Pursuant to 93/42/EEC and 2007/47 EEC of 09/02/2007
<b>INTENDED USE</b>	Disposable BEDPAN COVER for use in public and private hospitals, first aid, Emergency.
<b>RAW MATERIALS</b>	Polyethylene HDPE. White cross-linked sodium polyacrylate enclosed in a water-soluble film.



**HAPPYcover<sup>®</sup>**  


**PAD**

DATA SHEET PAGE 2:2

<p><b>INSTRUCTIONS FOR USE</b></p>	<p>Happy Cover<sup>®</sup> PAD is used to coat BEDPAN. SGAIEL<sup>®</sup> powder turns dejection into a gelatinous substance in about one minute. Remove bag from bedpan carefully in order to prevent spills. After use close the bag by adhesive liner or make a knot. Dispose as per current legislation.</p>
<p><b>CAUTION</b></p>	<p>Do not use sharp objects to open the package Do not disperse in the environment after use. Do not use if package is open or damaged.</p>
<p><b>CONTRAINDICATIONS</b></p>	<p>No contraindications by using these medical devices.</p>
<p><b>STORAGE</b></p>	<p>Store in a cool, dry place away from heat sources and sunlight.</p>
<p><b>TRANSPORTATION</b></p>	<p>Not regulated by EEC and Italia Transport Department as dangerous product. Terrestrial transport (RID/ADR): not dangerous. Maritime transport (IMG): not dangerous. Airplane transport (ICAD/IATA): not dangerous.</p>
<p><b>APPLICABLE REGULATION</b></p>	<p>93/42/EEC of 06/14/1993 2007/47 EEC of 09/05/2007</p> <p style="text-align: right;"><b>CE</b></p>
<p><b>GMDN</b></p>	<p>Code GMDN: 3923.21 Description: Polyethylene bag (BEDPAN COVER) with a folding and internal, double-sided adhesive pocket that allows an easy and safe sealing of its open upper side after use. Its interior is covered with a molecule powder named "SGAIEL<sup>®</sup>" with an absorbing function that turns liquids into solids. This powder is enclosed in a water-soluble film which acts upon contact with liquids. Ready for use.</p>
<p></p>	<p>Angimed S.r.l. Corso Umberto I, 590 ed.66 - 65015 Montesilvano (PE)</p>



**HAPPYcover<sup>®</sup>**  



**RIG**

DATA SHEET PAGE 1:2

<b>TRADING NAME</b>	<b>HAPPY COVER<sup>®</sup> RIG</b>
<b>DESCRIPTION</b>	Polyethylene vomit bag, for an easy and safe seal of its opened upper side, after use. Inside the bag there is a superabsorbent molecule-powder named SGAIEL <sup>®</sup> that transforms liquids into solids. This powder is enclosed in a water-soluble film which acts in contact with liquids. Ready for use.
<b>TEM CODE</b>	AM003
<b>RDM (Italy)</b>	1570463
<b>CND (Italy)</b>	A060399
<b>PACKAGING</b>	pcs 20
<b>DIMENSIONS</b>	mm. 250x250
<b>EXPIRING TIME</b>	5 years
<b>BELONGING CLASS OF MD</b>	Classe I medical device - Not sterile Pursuant to 93/42/EEC
<b>CE MARK</b>	Pursuant to 93/42/EEC and 2007/47 EEC of 09/02/2007
<b>INTENDED USE</b>	Disposable vomit system for use in public and private hospitals, first aid, Emergency and transport companies
<b>RAW MATERIALS</b>	Polyethylene HDPE. White cross-linked sodium polyacrylate enclosed in a water-soluble film.



DATA SHEET PAGE 2:2

<p><b>INSTRUCTIONS FOR USE</b></p>	<p>Happy Cover® RIG is not used in combination with other device. SGAIEL® powder turns vomit into a gelatinous substance in about one minute. After use close the bag by adhesive liner or make a knot. Dispose as per current legislation.</p>
<p><b>CAUTION</b></p>	<p>Do not use sharp objects to open the package Do not disperse in the environment after use. Do not use if package is open or damaged.</p>
<p><b>CONTRAINDICATIONS</b></p>	<p>No contraindications by using these medical devices.</p>
<p><b>STORAGE</b></p>	<p>Store in a cool, dry place away from heat sources and sunlight.</p>
<p><b>TRANSPORTATION</b></p>	<p>Not regulated by EEC and Italia Transport Department as dangerous product. Terrestrial transport (RID/ADR): not dangerous. Maritime transport (IMG): not dangerous. Airplane transport (ICAD/IATA): not dangerous.</p>
<p><b>APPLICABLE REGULATION</b></p>	<p>93/42/EEC of 06/14/1993 2007/47 EEC of 09/05/2007</p> <p style="text-align: right;"><b>CE</b></p>
<p><b>GMDN</b></p>	<p>Code GMDN: 3923.21 Description: Polyethylene sick bag with a folding and internal double-sided adhesive pocket that allows an easy and safe sealing of its open upper side after use. Its interior is covered with a molecule-powder named "SGAIEL®" with an absorbing function that turns liquids into solids. This powder is enclosed in a water-soluble film which acts upon contact with liquid. Ready for use.</p>
	<p>Angimed S.r.l. Corso Umberto I, 590 ed.66 - 65015 Montesilvano (PE)</p>

## DECLARATION OF CONFORMITY CE APPLICATIONS WITH DIRECTIVES 93/42 / EEC OF 14 JUNE 1993 - 2007/47 / EC OF SEPT EMBER 5, 2007

*Alessandro Di Teodoro, as the legal representative of Angimed S.r.l. ,  
declare*

*that medical devices class I below identified, meet the Essential Requirements of the Directive 93/42/EC of 14 June 1993 and 2007/47/EC of September 5, 2007 and are used for their intended use. This statement is based on the technical documentation referenced in the Procedure of Conformity Annex VII.*

ANGIMED SRL - MANUFACTURING PLANT:	
Corso Umberto I, 590 / ED.66 - 65015 Montesilvano (PE) Italy	
CODE	PRODUCT DESCRIPTION
AM001N	BEDPAN COVER
AM001	BEDPAN COVER
AM002	P.P. -URINE SYSTEM
AM002HM	P.P. -URINE SYSTEM WHIT LANYARD
AM003	RIG -VOMIT SYSTEM
AM003HM	RIG -VOMIT SYSTEM WHITH LANYARD
AM004	SGAIEL- POWDER SUPER ABSORBENT 400GR.
AM005	SGAIEL- POWDER SUPER ABSORBENT 200GR.
AM006	SGAIEL- POWDER SUPER ABSORBENT 1000GR.
AM007	SGAIEL- POWDER SUPER ABSORBENT WITH IDRO-SOLUBLE FILM - 05 GR.
AM008	SGAIEL- POWDER SUPER ABSORBENT WITH IDRO-SOLUBLE FILM - 10 GR.
AM009	SGAIEL- POWDER SUPER ABSORBENT WITH IDRO-SOLUBLE FILM - 15 GR.
AM013	SGAIEL- POWDER SUPER ABSORBENT WITH IDRO-SOLUBLE FILM - 13 GR.
AM020	SGAIEL- POWDER SUPER ABSORBENT WITH IDRO-SOLUBLE FILM - 20 GR.
AM025	SGAIEL- POWDER SUPER ABSORBENT WITH IDRO-SOLUBLE FILM - 25 GR.
AM050	SGAIEL- POWDER SUPER ABSORBENT WITH IDRO-SOLUBLE FILM - 50 GR.

*Alessandro Di Teodoro  
(C.E.O.)*