

BEL Engineering introduces the new series of semiMicro HPB balances with a resolution of 0.01 mg. Equipped with the new highly integrated electromagnetic cell, with automatic internal calibration on all models, they are characterized by a high level of precision and repeatability. The -ION models are equipped with an integrated ioniser on the back of the balance with the purpose to neutralize static electricity on charged samples like plastic parts, containers or films for acquiring the most accurate results. As in the tradition of BEL products, new semiMicro HPB balances deliver an excellent performance / price ratio, thanks to the continuous research and study of Bel Engineering in the high resolution and precision weighing field

Features

- Electromagnetic force restoration weighing system
- Highly integrated weighing cell
- Automatic internal calibration
- Efficient filtering to deliver quick and accurate results
- LCD with backlight with adjustable contrast
- Large glass draught shield with 3 sliding doors for easy access to the items being weighed
- Plastic protection cover for keypad

Accessories

- Portable ioniser Ion-A15 (BL0371)
- BEL Tablet System for balances. USB version (BL0329) or BLUETOOTH version (BL0330).
- Bluetooth printer DPP-250-BT for BEL Tablet system (AC013)
- Serial Printer TLP-50, with date/time (C054)
- Serial Printer DPP-250 (AC007)
- Density KIT for solids and liquids (AC002)
- Serial Printer STAT, with Statistics functions (AC022)
- Communication KIT for balances, only for 0,1mg res. (BL0476)
- Alphanumeric external keyboard (AC005)
- Serial to USB converter (E1002)
- Serial cable for serial output to printer or PC (E743)
- Factory Calibration certificate (BL0333)

HPB - Resolution 0,01mg-0,1mg



Model	Capacity (g)	Resolution (mg)	Pan (mm)	Repeatability (mg)	Linearity (mg)	Response time (sec.)
Automatic internal calibration						
HPB-425i	42	0,01	Ø 80	0,05	± 0,05	≤ 6
HPB-625i	62	0,01	Ø 80	0,05	± 0,08	≤ 6
HPB-105i	102	0,01	Ø 80	0,05	± 0,08	≤ 6
HPB-1245Di	42/120	0,01/0,1	Ø 80	0,02 / 0,1	± 0,05 / ± 0,2	≤ 6
HPB-1265Di	62/120	0,01/0,1	Ø 80	0,03 / 0,1	± 0,08 / ± 0,2	≤ 6
HPB-2245Di	42/220	0,01/0,1	Ø 80	0,02 / 0,1	± 0,1 / ± 0,2	≤ 6
HPB-2265Di	62/220	0,01/0,1	Ø 80	0,03 / 0,1	± 0,1/ ± 0,2	≤ 6
HPB-2285Di	82/220	0,01/0,1	Ø 80	0,05 / 0,1	± 0,1 / ± 0,2	≤ 6
HPB-22105Di	102/220	0,01/0,1	Ø 80	0,05 / 0,1	± 0,1 / ± 0,2	≤ 6
HPB 414Ai	410	0,1	Ø 80	0,1mg if weight ≤ 220g 0,3mg from 220g to 310g 0,35mg from 310g to 410g	± 0,0005	≤ 4
HPB 514Ai	510	0,1	Ø 80	0,1mg if weight ≤ 220g 0,3mg from 220g to 310g 0,35mg from 310g to 510g	± 0,0005	≤ 4
HPB 614Ai	610	0,1	Ø 80	0,1mg if weight ≤ 220g 0,3mg from 220g to 310g 0,35mg from 310g to 610g	± 0,0005	≤ 4

⁽¹⁾ If connected to optional STAT printer (AC022)



A-ION models equipped with integrated ioniser on the back of the balance.

The ioniser has two way of operating: 2 minutes or continuously for max 8 hours.

Balance technical data

Common features

- Span drift (+ 10...+ 30 °C): +/- 2ppm/°C
- Dimensions LxWxH (mm): 345x215x345
- Weighing chamber dimensions LxWxH (mm): 162x171x225
- Dimensions of external electronics box LxWxH: 175x105x50 (mm), weight 500g
- Power supply 110-230Vac, 50/60Hz; output 24V 500mA 13VA
- Net Weight: 7Kg
- Net Weight -ION models: 7,4 Kg

Ioniser technical data

- Operation modes: 2 minutes or continuously for max 8 hours
- Ozone concentration: 0-0.05 ppm (2cm from ion source)
- Distance sample-ion source: about 5-40 cm
- Discharge time: 9s/5cm, 13s/10cm
- Ambient conditions: 0-50°C, 20-80% air humidity (non-condensing)
- Ioniser Power supply: IN AC 100-240V 50/60Hz, OUT DC 12V, 500mA



Model	Capacity (g)	Resolution (mg)	Pan (mm)	Repeatability (mg)	Linearity (mg)	Response time (sec.)
-------	--------------	-----------------	----------	--------------------	----------------	----------------------

Automatic internal calibration

HPB-425i-ION	42	0,01	Ø 80	0,05	± 0,05	≤ 6
HPB-625i-ION	62	0,01	Ø 80	0,05	± 0,08	≤ 6
HPB-105i-ION	102	0,01	Ø 80	0,05	± 0,08	≤ 6
HPB-1245Di-ION	42/120	0,01/0,1	Ø 80	0,02 / 0,1	± 0,05 / ± 0,2	≤ 6
HPB-1265Di-ION	62/120	0,01/0,1	Ø 80	0,03 / 0,1	± 0,08 / ± 0,2	≤ 6
HPB-2245Di-ION	42/220	0,01/0,1	Ø 80	0,02 / 0,1	± 0,1 / ± 0,2	≤ 6
HPB-2265Di-ION	62/220	0,01/0,1	Ø 80	0,03 / 0,1	± 0,1 / ± 0,2	≤ 6
HPB-2285Di-ION	82/220	0,01/0,1	Ø 80	0,05 / 0,1	± 0,1 / ± 0,2	≤ 6
HPB-22105Di-ION	102/220	0,01/0,1	Ø 80	0,05 / 0,1	± 0,1 / ± 0,2	≤ 6
HPB 414Ai-ION	410	0,1	Ø 80	0,1mg if weight ≤ 220g 0,3mg from 220g to 310g 0,35mg from 310g to 410g	± 0,0005	≤ 4
HPB 514Ai-ION	510	0,1	Ø 80	0,1mg if weight ≤ 220g 0,3mg from 220g to 310g 0,35mg from 310g to 510g	± 0,0005	≤ 4
HPB 614Ai-ION	610	0,1	Ø 80	0,1mg if weight ≤ 220g 0,3mg from 220g to 310g 0,35mg from 310g to 610g	± 0,0005	≤ 4