

## Cleaning brush and valve cleaning brush

Reference	Designation	Working channel diameter (mm)	Catheter	Brush (mm)		
			Length and Colour (cm)	Single / Double	Length	Diameter
<b>Endoscope small diameter : bronchial and pediatric</b>						
<b>CJ-EDB-120-02</b>	Bronchial pediatric endoscope	1.5	120 White	●—●	12 / 12	2
<b>CJ-EDB-230-02</b>	For Air/Water channel	1.2 - 1.5	230 White	●—●	12 / 12	2
<b>CJ-EDB-120-03</b>	Bronchial endoscope	≥ 2.0	120 Yellow	●—●	20 / 20	3
<b>CJ-EDB-230-03</b>	Pediatric endoscope	≥ 2.0	230 Yellow	●—●	20 / 20	3
<b>CJ-EDB-300-03*</b>	Pediatric Enteroscope	≥ 2.2	300 Yellow	●—●	20 / 20	3
<b>Adult : Endoscope, enteroscope and cystoscope</b>						
<b>CJ-ESB-020-06</b>	For beginning of Air/Water channel	≥ 2.8 (2.8 to 5)	20 Light blue	—●	20	6
<b>CJ-ESB-230-06</b>	Endoscope adult	≥ 2.8 (2.8 to 5)	230 Blue	—●	20	6
<b>CJ-EDB-120-06</b>	Endoscope and cystoscope	≥ 2.8 (2.8 to 5)	230 Blue	●—●	20 / 20	6
<b>CJ-EDB-230-06</b>	Endoscope adult	≥ 2.8 (2.8 to 5)	230 Blue	●—●	20 / 20	6
<b>CJ-EDB-300-06*</b>	Enteroscope	≥ 2.8 (2.8 to 5)	300 Blue	●—●	20 / 20	6
<b>Valve cleaning brush</b>						
<b>CJ-GDB-11-05</b>	Valve cleaning brush	/	16 Sky blue	●=●	32 / 20	11 / 5

\* Special request