



CENGİZ MEDİKAL DIŞ TİC. LTD. ŞTİ.

İkitelli OSB Mutsan Sitesi M8 Blok No:17 Başakşehir / İstanbul

cengizmedikal34@gmail.com

Tel:+90 212 4854377-87

Date: 21.09.2022

Product Description	Size	Weight/ Roll	Product Code	Classification
Cengiz Boz GYPSUM (POP) BANDAGE FAST 3MINS	10cm x 270cm	100-110 gr	CB1027-922M3	Class I type1
Cengiz Boz GYPSUM (POP) BANDAGE FAST 3MINS	15cm x 270cm	155-165gr	CB1527-922M3	Class I type1
Cengiz Boz GYPSUM (POP) BANDAGE FAST 3MINS	20cm x 270cm	205-220gr	CB2027-922M3	Class I type1
Cengiz Boz GYPSUM (POP) BANDAGE SLOW MAX 12MINS	10cm x 270cm	100-110 gr	CB1027-922M12	Class I type1
Cengiz Boz GYPSUM (POP) BANDAGE SLOW MAX 12MINS	15cm x 270cm	155-165gr	CB1527-922M12	Class I type1
Cengiz Boz GYPSUM (POP) BANDAGE SLOW MAX 12MINS	20cm x 270cm	205-220gr	CB2027-922M12	Class I type1
FEATURES				

- Plaster of Paris bandage made of 100% cotton gauze.
- Gypsum covered in both sides with natural gypsum
- Pop bandages plastic tube rolls for easy usage.
- Application requires soaking in water 10 seconds in the water temperature not less than 20 degrees celcius.
- Initial setting time is in the fast freze products are 3 minutes, in the slow freze products are 8-12 minutes.
- Complete stiffening time takes place within 24 hours.
- Our pop bandages doesnt absorb Xray.
- Unsterile and single use, shelf life is 5 years.


Asim HOCAOGLU, Rh.D.
Head of Medical Devices
Registration and Coordination Department