



CT Scanner Accessories





1. Arm holders

► Arm up holder

Used to support the patient's arms when the arms are elevated.



CAAH-014A

Dimensions	Materials
W: 530 mm L: 378 mm H: 288 mm	Poly dichloropentadiene

The metal-free arm up holder is made of a material that minimizes artifacts. It is also light weight and washable. The metal-free arm up holder simplifies patient positioning.


► Handle

This handle is grasped by the patient to hold up his/her arms.




CTA-AHC

Dimensions	Materials
W: 350 mm L: 75 mm H: 165 mm	Iron Chrome-plated



► Arm down holder (one side)

These holders are mounted to the couch top to support the patient arm (one side only).




CAAH010A

For 470-mm-wide couch tops.

Dimensions	Materials
W: 335 mm L: 400 mm H: 48 mm	Polyvinyl chloride

► Arm down holder (both sides)

These holders are mounted to the couch top to support the patient's arms (both sides).




CAAH012A

For 470-mm-wide couch tops.

Dimensions	Materials
W: 700 mm L: 400 mm H: 53 mm	ABS resin

2. Baby bed / Child holder

► **Baby bed**



CTA-TMB

Placed on the couch top to restrain infant movement.

Dimensions	Materials
W: 230 mm L: 510 mm H: 80 mm	Acrylic resin Nylon Protective felt

► **Child holder**

This holder immobilizes infants for a wide range of body and head examinations.



CACH-001A

Dimensions	Materials
W: 370 mm L: 1000 mm H: 200 mm	Base: PET resin, Polyethylene Pillow: Polyethylene foam Mat, cushions: Urethane leather, Flexible urethane foam Bands, Body cover: Urethane leather, Polyester



3. Head rests

Used to immobilize the patient's head during scanning.

► **Head rests**



CAHR024A

The tilt angle can be set manually.

Dimensions	Materials
W: 200 mm L: 413 mm H: 148 mm	Carbon fiber Iron



CAHR-026A

For 3D subtraction.

Dimensions	Materials
W: 358 mm L: 629 mm H: 374 mm Forward/backward stroke: 160 mm Vertical stroke: 60 mm	Carbon fiber Polyvinyl chloride Iron



CAHR-028A

Dimensions	Materials
W: 190 mm L: 374 mm H: 208 mm	Carbon fiber reinforced plastic



CAHR-029A

Dimensions	Materials
W: 221 mm L: 478 mm H: 207 mm	Carbon fiber reinforced plastic



CAHR-030A

Dimensions	Materials
W: 220 mm L: 448 mm H: 207 mm	Carbon fiber reinforced plastic


3. Head rests

Used to immobilize the patient's head during scanning.

► Coronal head rests



CTA-HRJ (prone)	
Used when the patient is prone for coronal image acquisition.	
Dimensions	Materials
W: 230 mm	Carbon fiber Iron
L: 375 mm	
H: 120 mm	



CTA-HRO (supine)	
Used when the patient is supine for coronal image acquisition.	
Dimensions	Materials
W: 190 mm	Acrylic resin Nylon Iron
L: 360 mm	
H: 170 mm	

► Head rest for infant




CTA-HRN	
Used to immobilize the heads of infants.	
Dimensions	Materials
W: 150 mm	Acrylic resin Nylon Iron
L: 310 mm	
H: 137 mm	

4. Flat board

Placed on the couch top to make the couch-top upper surface flat.

► Flat board




CABM021A	
For 470-mm-wide couch tops.	
Dimensions	Materials
W: 440 mm	Acrylic foam Polyester Cotton
L: 1905 mm	
H: 40 mm	



5. Foot supports

Mounted to the end of the couch top to provide additional table length for feet support.

► Foot rests




CAEX-008A

400 mm extension
(according to foot mat sizes)

Dimensions	Materials
W: 470 mm	Carbon fiber
L: 450 mm	
H: 80 mm	

(To be exchanged with the head rest.)
(With foot-mat fastener)



CAEX-012A

310 mm extension
(according to foot mat sizes)

Dimensions	Materials
W: 465 mm	Carbon fiber reinforced plastic
L: 409 mm	
H: 74 mm	

CAEX-013A

400 mm extension
(according to foot mat sizes)

Dimensions	Materials
W: 446 mm	Carbon fiber reinforced plastic
L: 462 mm	
H: 56 mm	

6. Head mats

► Head mat (cushion)

Placed on the head rest.



CAHM-003A

Dimensions	Materials
W: 231 mm	Polyethylene foam
L: 323 mm	
H: 200 mm	



CAHM-004A

Dimensions	Materials
W: 390 mm	Polyethylene
L: 230 mm	
H: 20 mm	



CAHM-005A

Dimensions	Materials
W: 226 mm	Polyethylene
L: 281 mm	
H: 168 mm	

6. Head mats

► Wedge mats



CAWM-001A/1 10 degrees

Dimensions	Materials
W: 186 mm L: 180 mm H: 78 mm	Polyethylene foam

CAWM-002A/1 10 degrees

Dimensions	Materials
W: 189 mm L: 210 mm H: 83 mm	Polyethylene foam

CAWM-001A/2 20 degrees

Dimensions	Materials
W: 186 mm L: 180 mm H: 108 mm	Polyethylene foam

CAWM-002A/2 20 degrees

Dimensions	Materials
W: 190 mm L: 210 mm H: 121 mm	Polyethylene foam

CAWM-001A/3 30 degrees

Dimensions	Materials
W: 186 mm L: 180 mm H: 128 mm	Polyethylene foam

CAWM-002A/3 30 degrees

Dimensions	Materials
W: 186 mm L: 210 mm H: 149 mm	Polyethylene foam

► Table top pillow



CATP-001A

Dimensions	Materials	Color
W: 188 mm L: 274 mm H: 115 mm	Urethane foam Urethane resin	Terra Gray

CATP-002A

Dimensions	Materials	Color
W: 188 mm L: 274 mm H: 115 mm	Urethane foam Urethane resin	Angora Gray

► Head cushion

Placed between the patient's head and the head straps to immobilize the head.



CAHC-001A

Dimensions	Materials
W: 220 mm L: 80 mm H: 20 mm	Urethane foam Urethane sheet

► Side mats

Placed between the side of the patient's head and the head rest to immobilize the head.



CTA-HMC

Dimensions	Materials
W: 140 mm L: 80 mm H: 30 mm	Urethane Polyester Cotton

CASM-001A

Dimensions	Materials
W: 140 mm L: 80 mm H: 30 mm	Polyurethane Polyester Cotton

► Neck mat



CANM-001A (long) For CAHR-029A

Dimensions	Materials
W: 187 mm L: 470 mm H: 40 mm	Urethane foam Urethane sheet



CANM-001B (short)

Dimensions	Materials
W: 164 mm L: 470 mm H: 40 mm	Urethane foam Urethane sheet



7. Body mats

Placed on the couch top.

► Body mats



CABM-031A/1

For 2000-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2360 mm H: 20 mm	Polyurethane	Zenith Blue

CABM-031A/4

For 300 kg patient couch with a 1500-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 1860 mm H: 20 mm	Polyurethane	Zenith Blue



CABM-031A/5

For 2000-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2360 mm H: 20 mm	Polyurethane	Terra Gray

CABM-031A/6

For 1800-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2090 mm H: 20 mm	Polyurethane	Terra Gray

CABM-031A/7

For 1500-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 1790 mm H: 20 mm	Polyurethane	Terra Gray

CABM-031A/8

For 300 kg patient couch with a 1500-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 1860 mm H: 20 mm	Polyurethane	Terra Gray

CABM-032A/2

For 1700-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2070 mm H: 20 mm	Polyurethane	Terra Gray



CABM-041A

For 2000-mm-long couch top with mat covers on either side of the long side of the mat.

Dimensions	Materials	Color
W: 990 mm L: 2356 mm H: 20 mm	Urethane foam Urethane sheet	Angora Gray

► Body mat covers



CAMC-001A/1

Dimensions	Materials
W: 1020 mm L: 2900 mm	Polyvinyl chloride

CAMC-001A/2

Dimensions	Materials
W: 1020 mm L: 2600 mm	Polyvinyl chloride

7. Body mats

Placed on the couch top.

► Triangular mats

Placed under the patient's knees during scanning.



CAFM-002A/1

Dimensions	Materials	Color
W: 350 mm L: 570 mm H: 220 mm	Polyvinyl chloride Urethane	Win Blue



CAFM-002A/2

Dimensions	Materials	Color
W: 350 mm L: 570 mm H: 220 mm	Polyurethane	Terra Gray



CAFM-007A

Dimensions	Materials	Color
W: 350 mm L: 570 mm H: 220 mm	Polyurethane	Angora Gray

► Foot mats




CAFM-001B For CAEX-008A

Dimensions	Materials	Color
W: 470 mm L: 400 mm H: 20 mm	Polyvinyl chloride Urethane	Zenith Blue



CAFM-003A/2 For CAEX-008A

Dimensions	Materials	Color
W: 470 mm L: 400 mm H: 20 mm	Polyurethane	Terra Gray




CAFM-006A For CAEX-013A

Dimensions	Materials	Color
W: 400 mm L: 470 mm H: 20 mm	Urethane foam Urethane sheet	Angora Gray



8. Head straps

► Head straps (fastener)



CAHB-003A/1

Used to immobilize the head.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Blue



CAHB-003A/2

Used to immobilize the head.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Terra Gray

► Chin straps (fastener)



CAHB-002A/1

Used to immobilize the chin.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Blue



CAHB-002A/2


Used to immobilize the chin.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Terra Gray

9. Body straps

► Body straps (fastener)

Sliding type bands to be used to restrain patient movement.



CABB-008A/1, /2*

For 2000-mm-long couch top

CABB-008B/1, /2*

For 1800-mm-long couch top

CABB-008C/1, /2*

For 300 kg patient couch with a 1500-mm-long couch top

CABB-008D/1, /2*

For 1500-mm-long couch top

CABB-010A/1

For Aquilion Precision with a 2000-mm-long couch top

CABB-010A/2

For 1700-mm-long couch top

CABB-011A/1

For Aquilion ONE / INSIGHT Edition with a 2000-mm-long couch top



• Body band rails

For 2000-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2355 mm	Polyethylene

For 1800-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2085 mm	Polyethylene

For 300 kg patient couch with a 1500-mm-long couch top

Dimensions	Materials
W: 143 mm L: 1855 mm	Polyethylene

For 1500-mm-long couch top

Dimensions	Materials
W: 143 mm L: 1785 mm	Polyethylene

For Aquilion Precision with a 2000-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2305 mm	Polyethylene

For 1700-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2015 mm	Polyethylene

For Aquilion ONE / INSIGHT Edition with a 2000-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2405 mm	Polyethylene

• Body band 1

Dimensions	Materials
W: 300 mm L: 513 mm	Polyethylene

• Body band 2

Dimensions	Materials
W: 200 mm L: 600 mm	Polyethylene

• Body band 3

Dimensions	Materials
W: 200 mm L: 513 mm	Polyethylene

*: The color of the body straps for /1 differs from that for /2.
The color of the body straps for /2 is the same as that of body mats CABM-031A/5-8.



10-a. Compatibility chart

Item	Model	Aquilion Precision	Aquilion ONE /PRISM Edition	Aquilion ONE /GENESIS Edition	Aquilion Prime SP	Aquilion Lightning TSX-036A	Aquilion Lightning TSX-035A	Aquilion Exceed LB	Aquilion LB	Aquilion Start
1. Arm holders										
Arm up holder	CAAH-014A	◎	◎	◎	◎	○	○	◎	◎	○
Handle	CTA-AHC	—	—	—	—	—	○	—	○	—
Arm down holder (one side)	CAAH010A	—	—	—	—	—	○	—	○	—
Arm down holder (both sides)	CAAH012A	—	—	—	—	—	○	—	○	—
2. Baby bed / Child holder										
Baby bed	CTA-TMB	—	—	—	—	—	○	—	○	○
Child holder	CACH-001A	—	○	○	○	○	—	○	—	—
3. Head rests										
Head rest	CAHR024A	—	—	—	—	—	—	—	○	—
Head rest	CAHR-026A	—	—	—	—	—	○	—	○	—
Head rest	CAHR-028A	◎	—	—	—	—	—	—	—	—
Head rest	CAHR-029A	—	◎ ^{*1}	—	○ ^{*3}	○ ^{*3}	—	◎	—	—
Coronal head rest (prone)	CTA-HRJ	—	—	—	—	—	○	—	○	○
Coronal head rest (supine)	CTA-HRO	—	—	—	—	—	—	—	○	—
Head rest for infant	CTA-HRN	—	—	—	—	—	○	—	○	—
4. Flat board										
Flat board	CABM021A	—	○	○	○	○	○	○	○	○
5. Foot supports										
Foot rest	CAEX-008A	—	◎	◎	◎	○	○	◎	○	○
Foot rest	CAEX-012A	◎	—	—	—	—	—	—	—	—
6. Head mats										
Head mat (cushion)	CAHM-003A	—	◎ ^{*2}	◎	◎	◎	◎	—	◎	◎
Head mat (cushion)	CAHM-004A	◎	—	—	—	—	—	—	—	—
Head mat (cushion)	CAHM-005A	—	◎ ^{*1}	—	○ ^{*3}	○ ^{*3}	—	◎	—	—
Wedge mat (10 degrees)	CAWM-001A/1	○	—	—	—	—	—	—	—	—
Wedge mat (20 degrees)	CAWM-001A/2	◎	—	—	—	—	—	—	—	—
Wedge mat (30 degrees)	CAWM-001A/3	○	—	—	—	—	—	—	—	—
Wedge mat (10 degrees)	CAWM-002A/1	—	○	○	○	○	—	◎	—	—
Wedge mat (20 degrees)	CAWM-002A/2	—	◎	○	○	○	—	◎	—	—
Wedge mat (30 degrees)	CAWM-002A/3	—	◎	○	○	○	—	◎	—	—
Table top pillow	CATP-001A	◎	◎	—	—	—	—	○	—	—
Side mat	CTA-HMC	—	—	—	—	—	○	—	—	—
Side mat	CASM-001A	◎	◎	◎	◎	○	—	◎	◎	○
Neck mat (long)	CANM-001A/1	—	○ ^{*1}	—	○ ^{*3}	○ ^{*3}	—	○	—	—
Neck mat (short)	CANM-001B/1	—	—	○	○	○	—	—	—	—

Note: Please check with your sales representative as some of these systems may not be available in your country/region.

*1: Production on or after July 2020.

*2: Production on or before June 2020.

*3: When the Headrest (CAHR-029A) is purchased as an additional option, the compatible head mat is CAHM-005A and the neck mat is CANM-001A

◎: Standard
○: Option
—: Not available

10-a. Compatibility chart

Item	Model	Aquilion Precision	Aquilion ONE /PRISM Edition	Aquilion ONE /GENESIS Edition	Aquilion Prime SP	Aquilion Lightning TSX-036A	Aquilion Lightning TSX-035A	Aquilion Exceed LB	Aquilion LB	Aquilion Start
7. Body mats										
Body mat	CABM-031A/1	—	—	—	—	—	⊙ *1	—	⊙ *1	—
Body mat	CABM-031A/4	—	—	—	—	—	⊙ *4	—	⊙ *4	—
Body mat	CABM-031A/5	⊙ *1	⊙ *1	⊙ *1	⊙ *1	⊙ *1	—	⊙ *1	—	—
Body mat	CABM-031A/6	—	⊙ *2	⊙ *2	⊙ *2	⊙ *2	—	—	—	⊙ *2
Body mat	CABM-031A/7	—	⊙ *3	⊙ *3	⊙ *3	⊙ *3	—	—	—	⊙ *3
Body mat	CABM-031A/8	—	⊙ *4	⊙ *4	⊙ *4	⊙ *4	—	⊙ *4	—	—
Body mat	CABM-032A/2	⊙ *5	—	—	—	—	—	—	—	—
Body mat cover	CAMC-001A/1	⊙ *1	—	—	—	—	—	—	—	—
Body mat cover	CAMC-001A/2	⊙ *5	—	—	—	—	—	—	—	—
Triangular mat	CAFM-002A/1	—	—	—	—	—	○	—	⊙	—
Triangular mat	CAFM-002A/2	⊙	⊙	⊙	⊙	○	—	⊙	—	○
Foot mat for CAEX-008A	CAFM-001B	—	—	—	—	○	○	—	○	—
Foot mat for CAEX-008A	CAFM-003A/2	—	⊙	⊙	⊙	○	—	⊙	—	○
Foot mat for CAEX-012A	CAFM-005A	⊙	—	—	—	—	—	—	—	—
8. Head straps										
Head strap (fastener)	CAHB-003A/1	—	—	—	—	—	⊙	—	⊙	—
Head strap (fastener)	CAHB-003A/2	⊙	⊙	⊙	⊙	⊙	—	⊙	—	⊙
Chin strap (fastener)	CAHB-002A/1	—	—	—	—	—	⊙	—	⊙	—
Chin strap (fastener)	CAHB-002A/2	⊙	⊙	⊙	⊙	⊙	—	⊙	—	⊙
9. Body straps										
Body strap (fastener)	CABB-008A/1	—	—	—	—	—	⊙ *1	—	⊙ *1	—
Body strap (fastener)	CABB-008A/2	—	⊙ *1	⊙ *1	⊙ *1	⊙ *1	—	⊙ *1	—	—
Body strap (fastener)	CABB-008B/1	—	—	—	—	—	⊙ *2	—	⊙ *2	—
Body strap (fastener)	CABB-008B/2	—	⊙ *2	⊙ *2	⊙ *2	⊙ *2	—	—	—	⊙ *2
Body strap (fastener)	CABB-008C/1	—	—	—	—	—	⊙ *4	—	⊙ *4	—
Body strap (fastener)	CABB-008C/2	—	⊙ *4	⊙ *4	⊙ *4	⊙ *4	—	⊙ *4	—	—
Body strap (fastener)	CABB-008D/1	—	—	—	—	—	⊙ *3	—	⊙ *3	—
Body strap (fastener)	CABB-008D/2	—	⊙ *3	⊙ *3	⊙ *3	⊙ *3	—	—	—	⊙ *3
Body strap (fastener)	CABB-010A/1	⊙ *1	—	—	—	—	—	—	—	—
Body strap (fastener)	CABB-010A/2	⊙ *5	—	—	—	—	—	—	—	—

*1: Dedicated for patient couches with a 2000-mm long couch top

*2: Dedicated for patient couches with a 1800-mm long couch top

*3: Dedicated for patient couches with a 1500-mm long couch top

*4: Dedicated for 300 kg patient couches with a 1500-mm long couch top

*5: Dedicated for patient couches with a 1700-mm long couch top

⊙: Standard

○: Option

—: Not available



10-b. Compatibility Chart

Item	Model	Aquilion ONE /INSIGHT Edition	Aquilion Serve SP	Aquilion Serve
1. Arm holders				
Arm up holder	CAAH-014A	◎	○	○
2. Baby bed/Child holder				
Child Holder	CACH-001A	○	○	○
3. Head rests				
Head rest	CAHR-030A	◎	—	—
Head rest	CAHR-029A	—	◎	◎
4. Flat board				
Flat board	CABM021A	○	○	○
5. Foot supports				
Foot rest	CAEX-008A	—	○	○
Foot rest	CAEX-013A	◎	—	—
6. Head mats				
Head mat (cushion)	CAHM-005A	◎	◎	◎
Wedge mat (10 degrees)	CAWM-002A/1	○	○	○
Wedge mat (20 degrees)	CAWM-002A/2	◎	◎	◎
Wedge mat (30 degrees)	CAWM-002A/3	◎	○	○
Table top pillow	CATP-001A	—	—	○
Table top pillow	CATP-002A	◎	○	—
Head cushion	CAHC-001A	◎	○	—
Side mat	CASM-001A	—	—	○
Neck mat (long)	CANM-001A/1	—	○	○
Neck mat (short)	CANM-001B/1	◎	—	—
7. Body mats				
Body mat	CABM-031A/5	—	◎ *1	◎ *1
Body mat	CABM-031A/6	—	◎ *2	◎ *2
Body mat	CABM-031A/7	—	◎ *3	◎ *3
Body mat	CABM-031A/8	—	◎ *4	◎ *4
Body mat	CABM-041A/1	◎ *1	—	—
Triangular mat	CAFM-002A/2	—	—	○
Triangular mat	CAFM-007A/1	◎	○	—
Foot mat for CAEX-008A	CAFM-003A/2	—	○	○
Foot mat for CAEX-013A	CAFM-006A	◎	—	—
8. Head straps				
Head strap (fastener)	CAHB-003A/2	◎	◎	◎
Chin strap (fastener)	CAHB-002A/2	◎	◎	◎
9. Body straps				
Body strap (fastener)	CABB-008A/2	—	◎ *1	◎ *1
Body strap (fastener)	CABB-008B/2	—	◎ *2	◎ *2
Body strap (fastener)	CABB-008C/2	—	◎ *4	◎ *4
Body strap (fastener)	CABB-008D/2	—	◎ *3	◎ *3
Body strap (fastener)	CABB-011A/1	◎ *1	—	—

Note: Please check with your sales representatives as some of these systems may not be available in your country/region.

*1 Dedicated for patient couches with a 2000-mm long couch top

*2 Dedicated for patient couches with a 1800-mm long couch top

*3 Dedicated for patient couches with a 1500-mm long couch top

*4 Dedicated for 300kg patient couches with a 1500-mm long couch top

*5 Dedicated for patient couches with a 1700-mm long couch top

◎: Standard

○: Option

—: Not available



CANON MEDICAL SYSTEMS CORPORATION

<https://global.medical.canon>

©Canon Medical Systems Corporation 2003-2024. All rights reserved.
Design and specifications are subject to change without notice.
MCACT0093EAK 2024-05 CMSC/Produced in Japan.

Canon Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485. Canon Medical Systems Corporation meets the Environmental Management System standard ISO 14001.

Made For life